

A rapid method for estimation of Chemical Oxygen Demand.

Field of the Invention

The present invention relates to a rapid method for estimation of Chemical Oxygen Demand (COD) of water and wastewater. It particularly relates to determination of COD of industrial waste or domestic waste water or for determining the degree of pollution or develop designs for effluent treatment plants or efficiency of treatment plants or quick, rapid and onsite estimation of COD of water and wastewater or in untreated municipal wastewater or activation and pre-clarification tank inlets or cooling water and storm water reservoirs.

Background and prior art references

Oxygen demand is a significant parameter for determining the effect of organic pollutants in water. As microorganisms in the environment ingest the organic material, oxygen is depleted. This in turn can be harmful to fish and plant life. (<http://www.spectronic.com/spectron/spctech2.htm>). Wastewater from food processing is non-toxic but organic. High concentration of nutrients can be harmful for the environment. Extra quantities of nitrogen, fat and phosphorus require more oxygen for bacteria to decompose. If Chemical Oxygen Demand (COD) / Biochemical oxygen demand (B.O.D) content of the water is excessive, the oxygen supply in the water may be depleted below the level required to sustain aquatic life. (Ref: http://info.rfisk.is/verkefni/1077/hfe_4f2.htm).

The Chemical Oxygen Demand (COD) determination is a measure of the oxygen equivalent of that portion of the organic matter in a sample that is susceptible to oxidation by a strong chemical oxidant under controlled conditions (American Public Health Association (APHA). (1998). In Standard methods for examination of water and wastewater. 20th Edn. American Public Health Association (APHA), American Water Works Association (AWWA), Water Pollution Control Federation (WPCF), Washington, DC).

The limitation of the test lies in its inability to differentiate between the biologically oxidizable and biologically inert material. COD determination has an advantage over BOD test in that the results can be obtained in less than five hours where as BOD

requires 3 to 5 days. Further the test is relatively easy and with not much interference. (American Public Health Association (APHA), 1976). In: Standard methods for examination of water and wastewater. 14th edn American Public Health Association (APHA), American Water Works Association (AWWA), Water Pollution Control Federation (WPCF), Washington, DC).

The main chemical compounds in wastewater are Chemical Oxygen Demand (COD), nitrogen, phosphorus, fats, oils and grease. COD and BOD₅ are important parameters for measurement of organic matter content and oxygen needed to decompose the organic compounds. During the decomposition of organic matter there is less oxygen available in the sea and no oxygen in some places. It is possible to calculate COD or BOD into standard personal units, 60 g of oxygen to decompose the organic compounds from one person per day or equaling of 135 g of oxygen to decompose Chemical Oxygen Demand (COD) in waste water (COD= 2.25 x BOD). (Ref: http://info.rfisk.is/verkefni/1077/hfe_4f2.htm).

The calorimetric dichromate reflux method is commonly used method for determining the Chemical Oxygen Demand (COD) content in a sample, and has been preferred over procedures using other oxidants because of superior oxidizing ability, applicability to a wide variety of samples and ease of manipulation (American Public Health Association (APHA). (1989). In: Standard methods for examination of water and wastewater. 17th edn. American Public Health Association (APHA), American Water Works Association (AWWA), Water Pollution Control Federation (WPCF), Washington, DC). The basic principle involved is oxidizing the most types of organic matters by a boiling mixture of chromic and sulfuric acids. A sample is refluxed in strongly acid solution with a known excess of potassium dichromate ($K_2Cr_2O_7$). After digestion, the remaining unreduced $K_2Cr_2O_7$ is titrated with ferrous ammonium sulfate to determine the amount of $K_2Cr_2O_7$ consumed and the oxidizable organic matter is calculated in terms of oxygen equivalent.

The determination of Chemical Oxygen Demand (COD) is widely used in municipal and industrial laboratories to measure the general level of organic contamination in waste water (Ref: <http://www.chemetrics.com/InstProd/COD.I.htm>).

Wide ranges of instruments are available for Chemical Oxygen Demand (COD) estimation. CHEMetrics' - employs EPA (Environmental Protection Agency, USA) approved Dichromate reactor digestion method. (<http://www.chemetrics.com/InstProd/COD.I.htm>). Spectronic Instruments provide SPECTRONIC^R Spectrophotometer and the Bioscience ACCU-Test^R system based on APHA method. (American Public Health Association (APHA). (1989). In: Standard methods for examination of water and wastewater. 17th edn American Public Health Association (APHA), American Water Works Association (AWWA), Water Pollution Control Federation (WPCF), Washington, DC). It uses a semi-micro conversion of the Standard Methods procedure for the determination of Chemical Oxygen Demand (COD). (Ref: <http://www.spectronic.com/spectron/spctech2.htm>).

In the North Dakota Department of Health (NDDH), Chemistry Division Chemical Oxygen Demand (COD) SOP (Ref: <http://www.health.state.nd.us/lab/METHODS/I-I-4.HTM>) COD determination is through comparison of absorbance for the standards and the samples in HACH COD reactor and Sequonic Turner model 390 Spectrophotometer.

ISCO-STP COD analyzers and the Phoenix analyzer uses ozone for rapid oxidation and provides continuous Chemical Oxygen Demand (COD) measurement corresponding directly to result of the laboratory standard dichromate method. (<http://www.isco/html/prdprCOD.html>).

Among the reference methods of Chemical Oxygen Demand (COD) determination, EPA method 410.1/ 410.2 and the Standard Method SM5220-D employ closed refluxing (Ref: [http:// splash.metrokc.gov/wir/envlab/LABGUIDE/ com_ref.htm](http://splash.metrokc.gov/wir/envlab/LABGUIDE/com_ref.htm)). The biodegradability test - 5-day and 28 day includes COD 5 and 15 day determinations in the reference OCED method 3101D -- SOQ Appendix B Methods used by Burlington Research, Inc. (Ref: <http://199.72.5.39/soqapb.htm>).

PASTEL UV^(R) system works exclusively in UV range of light spectrum with UV analyzer (Spectrophotometer). It requires 16 photodiodes and specially developed

pulsed deuterium lamp. This analyzer requires calibration for accuracy within a very specific sample category such as industrial wastewater. For each additional application a new base line spectra has to be customized for specific samples. In this method the co-relation is only for BOD and TSS (Total Suspended Solids) which is 90% when compared to the traditional test methods. For COD the co-relation is valid only upto 6000 mg/L COD only. It is a indirect estimation based on reference data stored in internal UV based software (Ref.: <http://www.azurenv.com/ispec.htm>). The major limitations in different cases are the economics of the methods used, instruments employed and the time consumed. CHEMetrics' uses a standard digester block for heating the sample for 2 h at 150°C. In addition, the method needs to employ a Chemical Oxygen Demand (COD) Photometer (A-1051C) or Spectrophotometer which accepts a 16 mm diameter cell. The major limitation is that the results are not reportable to USEPA. On the other hand, 13 mm cells can be read in CHEMetrics VVR and A1051 photometer. However, the results are good for monitoring purposes only.

BioScience's EPA Accepted ACCU-TEST^R method needs specific type of Chemical Oxygen Demand (COD) reagent vials, a heating block, data management software package, etc. are required. It makes each test quite costly. BioScience's EPA Accepted ACCU-TEST^R is available only in a medium sensitivity ranges for Spectrophotometry 5 to 4500 mg/L Chemical Oxygen Demand (COD). In order to use the BioScience COD vials in the SPECTRONIC^R 401, both the test tube holder and the Light shield are required where as the GENESYSTM 2 or GENESYSTM 5 require the cell holder Platform and its Test tube holder (Ref:<http://www.spectronic.com/spectron/spctech2.htm>).

In the North Dakota Department of Health (NDDH) Chemistry Division Chemical Oxygen Demand (COD.) SOP (Ref: <http://www.health.state.nd.us/lab/METHODS/I-I-4.HTM>) Chemical Oxygen Demand (COD) determination is in a narrow range of 0 to 150 and 150 to 1500 mg/L. The calorimetrically determined Chemical Oxygen Demand (COD) value is based on measuring the consumed oxidant in the chrome (Cr⁶⁺) and chromous (Cr³⁺) valence state. The method employs HACH COD reactor

and Sequonic Turner model 390 Spectrophotometer. The spiking solutions in the 0 to 150 mg/L range are 10, 25, 50 and 150 where as for 150 to 1500 mg/L ranges, 5 concentration ranges (100, 250, 500, 1000 and 1500) are employed. It needs a sample volume of 10 ml, in addition to heating for 2 hours at 150°C.

In the Manganese III method for Chemical Oxygen Demand (COD) analysis (U.S. Patent document 5,556,787 Sept., 1996) determination is done by using an analysis reagent comprised of a mixture of stabilized Manganese III ion and an inorganic non-oxidizing acid such as sulfuric or phosphoric acid. The method involves titration and calorimetric determination. The method involves digestion of the test sample by heating between 100 to 200°C for up to 2 hours.

A method based on redox cell involved treatment of organic carbon with an excess of an oxidizing agent. The remaining oxidizing agent is determined in a redox cell to ascertain how much of the oxidizing agent was consumed by the organic carbon. This was the COD. This is multi-step process involving a large number of reagents (U.S. Patent 3,930,798 Jan., 1976). Another redox based method involves ozone as an oxidizing agent (U.S. Patent 5,324,666 Jun., 1994). The accuracy of the process is limited only by the control range of the ozone generating pumps. It can work within a range of less than 0.5 ozone load i.e. the ratio of chemical oxygen demand to added ozone. By increasing the ozone input the measuring range may be extended upwards but with reduced accuracy.

All the methods employed so far have long refluxing or digestion periods and are followed by titration or spectrophotometric reading of the standards and the samples at different wavelengths.

In the present invention the major limitations involved in estimation of chemical oxygen demand of water and wastewater have been over come. The novelty of the present invention is in use of a rapid, simple and effective method for estimation of COD of a wide range of samples. Another novelty of the present invention is the reduction in time period taken for estimation of COD, from 2h to less than 5 minutes. In the present invention, parameters have been studied for preparation of a reaction

mixture, which is easy to handle. Another novelty of the present invention is in the use of very small quantities of reagents and the test sample. In the present invention the various reaction mixtures can be read instantly for a very wide range of COD loads by using minor equipment e.g. photometer. Another novelty of the present invention is the ability to conduct the estimation on site, avoiding any possible changes in the test sample quality. Yet another novelty of the present invention is the COD color chart is used for determination of COD. over a wide range. Another novelty of the present invention is the flexibility to use reagents for determination of COD. over a wide range. Another novelty of the present invention is the simplicity for preparing reference standards easily and rapidly. Another novelty of the present invention is the stability of color based reference standards for long period without preserving them at low temperatures. These can thus be used repeatedly.

OBJECTS OF THE INVENTION

The main object of present invention is to provide a rapid method for estimation of Chemical Oxygen Demand (COD) of water and wastewater, which obviates the drawbacks listed above.

Another object of the present invention is to provide a cheap and simple method for quick estimation of COD of water and wastewater.

Yet, another object of present invention is to provide a rapid and sensitive method which, capable of determining COD in small quantities of test sample.

Yet another object of the present invention is to provide a color chart over a wide range of COD values.

Yet, another object of the present invention is to provide a kit for COD. estimation.

Yet, another object of the present invention is to provide a method which does not require much technical skill and sophisticated equipment.

Another object of the present invention is to provide an effective process for COD estimation with very little loss of sample.

SUMMARY OF THE INVENTION

The present invention has solved the problem of longer duration and consumption of large quantities of chemicals in COD estimation. The method of the invention has removed the need for heating the reaction mixture for nearly 2 h at 148°C, which is conventionally employed by standard methods and other commercially available kits.

One mL of sample solution is mixed with three reagents, consisting of 0.02 g mercuric sulphate, 0.5 mL of 0.25 N potassium dichromate solution and 1.5 mL of sulphuric acid - silver sulphate in a sequential manner. Its COD is checked with the help of color chart. Sample(s) showing a COD value of more than 10000 mg/L are diluted. After establishing its approximate COD value, the sample is diluted further, if necessary, to achieve a COD value in the range of 300 to 500 mg/L. Reagents listed above are added to these diluted sample and read its OD at 585 and 635 nm. These OD values are used for calculating precise COD values by comparing it with standard glucose solution.

The main utility of the present invention is for determining the degree of pollution, to develop designs for effluent treatment plants and to determine efficiency of treatment plants. Monitoring of COD is important for design and operation of wastewater treatment equipment. Another utility of the present invention is to provide an efficient method for quick, rapid and onsite estimation of COD of water and wastewater. Other utilities include applications in untreated municipal wastewater, activation and pre-clarification tank inlets, and even cooling water and storm water reservoirs. It is also used as a standard parameter for characterization of wastewater loads or for proof of a required purification level.

DETAILED DESCRIPTION OF THE INVENTION

The chemical oxygen demand (COD) determines the amount of oxygen required for chemical oxidation of organic matter using a strong chemical oxidant such as potassium dichromate under reflux conditions. The test is widely used to determine: 1) The degree of pollution in water bodies and their self purification capacity, 2) Efficiency of treatment plants, 3) Pollution loads, 4) Provides rough idea of B.O.D, which can be used for B.O.D estimation. The conventional estimation is based on

the principle that most of the organic matter is destroyed when boiled with a mixture of potassium dichromate and H_2SO_4 producing CO_2 and H_2O . A sample is refluxed with a known amount of potassium dichromate. It is then titrated against ferrous ammonium sulphate. The amount is proportional to O_2 required to oxidize the organic matter. The method takes 2 to 3 hours. Using standard COD kits available in the market can also do COD estimations. The standard COD kit method involves the use of thermoreactor and a photometer can also take 2 to 3 hours to complete. Since the initial COD of the sample cannot be guessed, kits for different COD ranges have to be tried and even sample may still has to be diluted. Dilution of sample leads to further use of more kits, it thus makes the estimation more expensive. However, a technique has been devised here, which gives good estimate of COD within 5 to 10 min and consumes very small quantities of reagents compared to conventional standard methods. It can be carried out easily on site. The reaction mixture develops a color, which can be easily read and distinguished. In the first stage, a sample is added without any dilution. The color developed with the undiluted sample gives a direct estimate of COD or gives a clear indication of the extent to which a sample needs to be diluted before adding to the reaction mixture. For samples with COD value of more than 10000, the color of the sample and the different reagents will be C10000 (brown). Make 4 different dilutions of the sample i.e. 10, 100, 200 and 500. The various diluted samples can be tested in increasing order of their dilution. Once a diluted sample and the reagent mixture show a color between C100 (yellow) and C500 (sea green). This dilution can be 10 to 300 times, depending upon the initial COD of the sample.

Accordingly, the present invention provides a rapid method for semi-quantitative estimation of Chemical Oxygen Demand (COD), which comprises:

- i) preparing standard solutions by dissolving glucose in water,
- ii) mixing the standard solutions with different reagents, mercuric sulphate (HgSO_4), potassium dichromate ($\text{K}_2\text{Cr}_2\text{O}_7$) and sulphuric acid (H_2SO_4) silver sulphate reagent in a sequential manner.



- iii) mixing the sample with different reagents, mercuric sulphate (HgSO_4), potassium dichromate ($\text{K}_2\text{Cr}_2\text{O}_7$) and sulphuric acid (H_2SO_4) silver sulphate reagent in a sequential manner.
- iv) diluting the sample to a desired level and adding the reagents in a sequential manner.
- v) noting down the color of the reaction mixture visually and record the color code using a color chart within 1 minute of incubation and reading optical density (O.D.) of the sample.

Accordingly the present invention provides a method for the preparation of COD chart useful for the estimation of COD in a sample, said method comprising the steps of:

- a) preparing standard glucose solution by dissolving glucose in distilled water at a concentration ranging between 200 mg/L to 1,00,000 mg/L with a COD concentration ranging between 213 mg/L to 1,06,700 mg/L,
- b) mixing 1 mL of standard glucose solutions from step (a), individually, with three reagents, consisting of 0.02 g mercuric sulfate, 0.5 mL of 0.25 N potassium dichromate solution and 1.5 mL of sulfuric acid - silver sulfate in a sequential manner, and
- c) preparing a COD color chart of different glucose concentration as shown in figure 1 of accompanying drawings based on the COD concentration ranging between 213 mg/L to 1,06,700 mg/L , and which chart is useful for rough estimation of COD values by comparing the colors of the chart with the color of the samples.

The present invention also provides a rapid method for the estimation of COD of an effluent from domestic, industrial, municipal and other sources, said method comprising the steps of

- a) preparing a set of standard glucose solutions by dissolving glucose in distilled water at a concentration ranging between 300 mg/L to 500 mg/L, with a COD concentration of 320 mg/L to 535 mg/L,

- b) mixing 1 mL of standard glucose solutions from step (a), individually, with three reagents, consisting of 0.02 g mercuric sulfate, 0.5 mL of 0.25 N potassium dichromate solution and 1.5 mL of sulfuric acid - silver sulfate in a sequential manner,
- c) mixing 1 mL of a sample to be tested with 0.01 to 0.03 g mercuric sulfate per ml of a standard glucose solution, 0.5 mL of 0.25 N to 0.30 N potassium dichromate and 1.5 mL of sulfuric acid – silver sulfate reagents in a sequential manner and noting down the color and if the color of the reaction mixture turns (color code C10000) brown,
- d) diluting the sample according to table 1 (a) with distilled water till a particular color range of C100 to C10000 is achieved as shown in Figure 1 of accompanying drawing,
- e) further diluting the sample with distilled water to obtain a COD range in between 320 to 535 mg/L wherein the dilution is determined according to tables 1 (a) & (b),
- f) mixing 1 mL of diluted sample from step (e), individually, with three reagents, consisting of 0.02 g mercuric sulfate, 0.5 mL of 0.25 N to 0.30 N potassium dichromate solution and 1.5 mL of sulfuric acid - silver sulfate in the above sequential manner,
- g) matching the color of the reaction mixture of the above sample (step f) with the COD color chart provided in Figure 1 of the accompanying drawing,
- h) noting the optical density (OD) of the reaction mixture of standard glucose solutions (step b) and sample (step f) at 585 nm or 635 nm, against air,
- i) drawing a correlation between the OD values of above sample with OD values of standard glucose solutions,
- j) multiplying the OD value of sample (step h) with a correction factor of 1.2, comparing this value (A) with the OD value of standard glucose solution of COD value 320 to 535 mg/l,
- k) multiplying value (A) with the compared COD value of standard Glucose solution and divided by the respective OD value of the glucose standard (B), and

- l) then multiply value (B) with a dilution factor (in case of diluted sample) to obtain the COD value of the sample in terms of mg/l.

In yet another embodiment, the tables 1(a) and 1(b) are as follows:

Table 1(a)

Color of the reaction mixture	Required dilution of sample	Color code	Expected COD of unknown sample (mg/L)
C 10000	10	C1000	10000
C 10000	100	C100-C500	10000 – 50000
C 10000	200	C250- C500	50000 – 100000
C 10000	500	C200-C2000	100000 - 1000000

Table 1(b)

Color of the reaction mixture	Expected COD (mg/L)	Dilution required to bring COD value (300-500 mg/L)
C9000	9000-9500	18-31
C6000	6000-8500	12-28
C2000	2000-5500	4-18
C1500	1500-1900	3-6
C1200	1200-1400	2-5
C600	600-1100	2-4
C500	500-550	NIL
C400	400-450	NIL
C300	300-350	NIL
C200	200-250	NIL
C150	100-150	NIL
C100	00 –100	NIL

In yet another embodiment, in step (d) the dilution of samples, provide dark brown reaction mixture, falling in the range of 10 to 500 for approximate cod values of 10000 mg/l and above as given in Table1(c) and further dilution of nil to 31 is done to bring the COD in the range of 300 to 500 mg/l as given in Table 1(d).

In yet another embodiment, in step (e), the samples provide yellowish to blackish brown colored reaction mixture falling in the range of 100 to 9500 mg/l COD and the

samples are further diluted to 31 times to obtain a diluted sample in the COD range of 300 to 500 mg/l, according to the Table 1(d).

In yet another embodiment, in step (d) the dilution of samples and color groups range from (A) to (G) for approximate COD values lying in the range of 100 to 9500 mg/l as given in table 1(d)

In yet another embodiment, tables 1(c) and 1(d) are shown below.

TABLE 1(c)

Color of the reaction mixture	Color code of the reaction mixture as per the provided chart	Dilution required to bring the COD values in the range of 100 – 1000 mg/l	Color obtained on dilution as per the provided color chart	Expected COD of the reaction mixture (mg/l)	Expected COD of unknown sample (mg/l)
Dark brown	C 10000 to C 100000	i) 10	C 1000	1000	10000
		ii) 100	C 100 – C 500	100 to 500	10000 to 50000
		iii) 200	C 250- C 500	250 to 500	50000 to 100000
		iv) 500	C 200 to C 2000	200 to 2000	100000 to 1000000

Table 1(d)

Color group	Color of the reaction mixture	Color code of the reaction mixture as per the chart provided	Expected COD of the reaction mixture	Dilution required to bring the COD values in the range of 300 to 500 mg/l
A	Yellowish	C 100 to C 250	1 to 250	Nil
B	Yellowish	C 300 to C 550	300 to 550	Nil
C	Greenish blue	C 600 to C 1000	600 to 1000	2 to 4
D	Blue to grayish blue	C 1500 to C 2500	1500 to 2500	3 to 6
E	Dirty brown	C 3000 to C 5500	3000 to 5500	6 to 18
F	Cola brown	C 6000 to C 8500	6000 to 8500	12 to 28
G	Blackish brown	C 9000 to 9500	9000 to 9500	18-31

One more embodiment of the invention provides a kit for estimation of Chemical Oxygen Demand (COD) which comprises:

- a) reagents A) HgSO_4 (solid), B) Glucose standard solutions ranging from 300 to 500 mg/L, C). $\text{K}_2\text{Cr}_2\text{O}_7$ solution (0.25 to 0.3 N), and D) H_2SO_4 - silver sulfate reagent,
- b) a glass vial,
- c) a photometric cell,
- d) a COD color chart as shown in fig 1 of the accompanying drawing and a photometer suitable for wave lengths in the range of 585 to 635 nm.

One more embodiment the present invention a method for the preparation of reference tables 1(e) and 1(f) of COD values, useful for the estimation of COD in a sample, said method comprising the steps of:

- a) preparing standard glucose solution by dissolving glucose in distilled water at a concentration ranging between 200 mg/L to 1,00,000 mg/L with a COD concentration ranging between 213 mg/L to 1,06,700 mg/L,
- b) mixing 1 mL of standard glucose solutions from step (a), individually, with three reagents, consisting of 0.02 g mercuric sulfate, 0.5 mL of 0.25 N potassium dichromate solution and 1.5 mL of sulfuric acid - silver sulfate in a sequential manner, and
- c) preparing reference tables 1(e) and 1(f) for COD values of different glucose concentration based on the COD concentration ranging between 213 mg/L to 1,06,700 mg/L, and which tables are useful for rough estimation of COD values by comparing the colors indicated in the tables with the color of the samples.

In another embodiment the present invention provides the reference tables 1(e) and 1(f) are as shown below.

TABLE 1(e)

Color of the reaction mixture	Dilution required to bring the COD values in the range of 100 – 1000 mg/l	Expected COD of the reaction mixture (mg/l)	Expected COD of unknown sample (mg/l)
Dark brown	v) 10	1000	10000
	vi) 100	100 to 500	10000 to 50000
	vii) 200	250 to 500	50000 to 100000
	viii) 500	200 to 2000	100000 to 1000000

Table 1(f)

Color group	Color of the reaction mixture	Expected COD of the reaction mixture	Dilution required to bring the COD values in the range of 300 to 500 mg/l
A	Yellowish	1 to 250	Nil
B	Yellowish	300 to 550	Nil
C	Greenish blue	600 to 1000	2 to 4
D	Blue to grayish blue	1500 to 2500	3 to 6
E	Dirty brown	3000 to 5500	6 to 18
F	Cola brown	6000 to 8500	12 to 28
G	Blackish brown	9000 to 9500	18-31

In another embodiment, the present invention provides a rapid method for the estimation of COD of an effluent from domestic, industrial, municipal and other sources using reference tables as given above, said method comprising the steps of

- preparing a set of standard glucose solutions by dissolving glucose in distilled water at a concentration ranging between 300 mg/L to 500 mg/L, with a COD concentration of 320 mg/L to 535 mg/L,
- mixing 1 mL of standard glucose solutions from step (a), individually, with three reagents, consisting of 0.02 g mercuric sulfate, 0.5 mL of 0.25 N

potassium dichromate solution and 1.5 mL of sulfuric acid - silver sulfate in a sequential manner,

- c) mixing 1 mL of a sample to be tested with 0.02 g mercuric sulfate, 0.5 mL of 0.25 N potassium dichromate and 1.5 mL of sulfuric acid – silver sulfate reagents in a sequential manner and noting down the color till the color of the reaction mixture turns brown,
- d) diluting the sample as per tables 1(e) and 1(f) with distilled water till a particular color range shown in the tables is achieved,
- e) further diluting the sample with distilled water to obtain a COD range in between 320 to 535 mg/L wherein the dilution is determined according to tables 1 (e) and 1(f),
- f) mixing 1 mL of diluted sample from step (e), individually, with three reagents, consisting of 0.02 g mercuric sulfate, 0.5 mL of 0.25 N to 0.30 N potassium dichromate solution and 1.5 mL of sulfuric acid - silver sulfate in the above sequential manner,
- g) matching the color of the reaction mixture of the above sample (step f) with the COD color provided in tables 1(e) and 1(f),
- h) noting the optical density (OD) of the reaction mixture of standard glucose solutions (step b) and sample (step f) at 585 nm or 635 nm, against air,
- i) drawing a correlation between the OD values of above sample with OD values of standard glucose solutions,
- j) multiplying the OD value of sample (step h) with a correction factor of 1.2, comparing this value (A) with the OD value of standard glucose solution of COD value 320 to 535 mg/l,
- k) multiplying value (A) with the compared COD value of standard Glucose solution and divided by the respective OD value of the glucose standard (B) and
- l) then multiply value (B) with a dilution factor (in case of diluted sample) to obtain the COD value of the sample in terms of mg/l.

In yet another embodiment, in step d, the dilution of samples provide dark brown reaction mixture falling in the range of 10 to 500 for approximate cod values of

10000 mg/l and above as given in Table 1(e) and further dilution of nil to 31 is done to bring the COD in the range of 300 to 500 mg/l as given in Table 1(f).

In yet another embodiment, in step (e), the samples provide yellowish to blackish brown coloured reaction mixture falling in the range of 100 to 9500 mg/l COD and the samples are further diluted to 31 times to obtain a diluted sample in the COD range of 300 to 500 mg/l, according to the Tables 1(e) and 1(f).

In yet another embodiment, in step (d) the dilution of samples and color groups range from (A) to (G) for approximate COD values lying in the range of 100 to 9500 mg/l as given in table 1(d).

In one embodiment of the present invention, the Dilution Factor is the total volume used divided by volume of the sample.

In another embodiment, the effluent used is selected from domestic, industrial, municipal, agricultural and other waste material sources.

In still another embodiment of the invention, the waste material is in the form of solid, liquid, semi-solid or viscous forms.

In yet another embodiment, the reagents used are HgSO_4 , glucose standard solutions, $\text{K}_2\text{Cr}_2\text{O}_7$ solution and H_2SO_4 –silver sulfate reagent.

In yet another embodiment, the volume of the test sample or standard solution used ranges from 20 to 100 ml.

In yet another embodiment, the incubation is carried out for a period in the range between 15 seconds to 1 minute.

In yet another embodiment, the estimated COD is in the range between 80 to 106700 mg/L.

In yet another embodiment, the glucose concentration tested is ranging between 50 to 100000 mg/L.

In yet another embodiment, the amount of HgSO_4 used in the test samples is from 10 mg to 30 mg.

In yet another embodiment, the concentration of $\text{K}_2\text{Cr}_2\text{O}_7$ used in the test samples is from 0.25 Normal to 0.30 Normal.

[illegible]

In yet another embodiment, the optical density readings are measured for the COD values ranging from 320 to 535 mg/L.

In yet another embodiment, the sample used for COD estimation comprises of agricultural waste, municipal market waste, fruit and food industry waste, beverages, chemicals, microbes and animal waste etc.

In yet another embodiment of the present invention, the kit for estimation of Chemical Oxygen Demand (COD) which may comprise of:

- 17

The file of this patent contains at least one drawing executed in color. Copies of this patent with color drawing(s) will be provided by the Patent & Trademark Office upon request and payment of the necessary fee.

BRIEF DESCRIPTION OF THE ACCOMPANYING DRAWING

Fig. 1 shows the color code chart of COD values.

Having now generally described this invention, a further understanding can be obtained by reference to certain specific examples, which are provided herein for purposes of illustration only and are not intended to be limiting unless otherwise specified.

EXAMPLE 1

Dissolved 0.5 and 5 g glucose in 50 mL distilled water, separately. Aliquots from the stock solutions were taken and diluted with water to achieve different concentrations in the range of 200 to 5000 (COD: 213 to 5335 mg/L) and 10000 to 100000 mg/L (COD: 10670 to 106700 mg/L), respectively. Weighed 0.02 g HgSO_4 in separate 15 mL and 50 mL test tubes. 1.0 mL of glucose sample was added to each test tube and shake well. 0.5 ml of $\text{K}_2\text{Cr}_2\text{O}_7$ was added to each sample. 1.5 mL H_2SO_4 silver sulphate reagent was added slowly with continuous shaking. The color of the reaction mixture was noted down visually within 1 min of incubation. The color varied with concentration. In 200 to 5000 mg/L glucose solutions, the color varied from C200 to C5000. At 5000 mg/L glucose and higher concentrations, the final color was always turbid. However, at concentration above 10000 mg/L glucose concentrations, the final color was always brown (C10000) and turbid. Hence for a preliminary estimate of COD value, color of the standard reaction mixture based on a range of glucose concentrations can be used as a reference. For samples, with initial COD concentration of 10000 or more, the sample needs to be diluted and its COD, re-estimated by taking the dilution factor and color of the reaction mixture into consideration.

Table 1: Correlation between COD concentration and color of the reaction mixture.

Glucose Concentration (mg/L)	COD of Glucose the solution (mg/L)	Color of the reaction mixture
00 (Blank)	0.0	C00
200	213.4	C200
250	266.7	C250
300	320.1	C300
350	373.4	C350
400	426.8	C400
450	480.1	C450
500	533.5	C500
600	640.2	C600
700	746.9	C700
800	853.6	C800
900	960.3	C900
1000	1067.0	C100
1500	1600.5	C1500
2000	2134.0	C2000
2500	2667.5	C2500
3000	3201.0	C3000
3500	3734.5	C3500
4000	4268.0	C4000
4500	4801.5	C4500
5000	5335.0	C5000
10000	10670.0	C10000
50000	53350.0	C50000
100000	106700.0	C100000

EXAMPLE 2

Dissolved 0.5 and 5 g glucose in 50 mL and 100 mL distilled water, respectively. Aliquots from the stock solutions were taken and diluted with water to achieve different concentrations in the range of 50 to 5000 (COD: 80 to 5335 mg/L) and 5500 to 100000 mg/L (COD: 5868 to 106700 mg/L), respectively. Weighed 0.02 g HgSO_4 in separate 15 mL and 50 mL test tubes. 1.0 mL of glucose sample was added to each test tube and shake well. 0.5 mL of $\text{K}_2\text{Cr}_2\text{O}_7$ was added to each sample. 1.5 mL H_2SO_4 silver sulphate reagent was added slowly with continuous shaking. The color of the reaction mixture was noted down visually within 1 min of incubation. The color varied with concentration. In 50 to 9500 mg/L glucose solutions, the color varied from C50 to C9500. At 5500 mg/L glucose and higher concentrations, the final color was always turbid. However, at concentration above 10000 mg/L glucose

concentrations, the final color was always C10000. Hence, for a preliminary estimate of COD value, color of the standard reaction mixture based on a range of glucose concentrations can be used as a reference. For samples, with initial COD concentration of 10000 or more, the sample needs to be diluted and its COD, re-estimated by taking the dilution factor and color of the reaction mixture into consideration.

Table 2: Correlation between COD concentration and color of the reaction mixture.

Glucose Concen tration (mg/L)	COD of Glucose solution (mg/L)	Color of the reaction mixture
00 (Blank)	0.0	C00
50	80.3	C50
100	106.7	C100
150	241.0	C150
200	213.4	C200
250	266.7	C250
300	320.1	C300
350	373.4	C350
400	426.8	C400
450	480.1	C450
500	533.5	C500
550	583.8	C550
600	640.2	C600
650	693.5	C650
700	746.9	C700
750	800.2	C750
800	853.6	C800
850	906.9	C850
900	960.3	C900
950	1013.6	C950
1000	1067.0	C1000
1100	1173.7	C1100
1200	1280.4	C1200
1300	1387.1	C1300
1400	1493.8	C1400
1500	1600.5	C1500
1600	1707.2	C1600
1700	1813.9	C1700
1800	1920.6	C1800
1900	2027.3	C1900
2000	2134.0	C2000
2500	2667.5	C2500
3000	3201.0	C3000
3500	3734.5	C3500
4000	4268.0	C4000

4500	4801.5	C4500
5000	5335.0	C5000
5500	5868.5	C5500
6000	6402.0	C6000
6500	6935.5	C6500
7000	7469.0	C7000
7500	8002.5	C7500
8000	8536.0	C8000
8500	9069.5	C8500
9000	9603.0	C9000
9500	10136.5	C9500
10000	10670.0	C10000
15000	16005.0	C15000
20000	21340.0	C20000
30000	32010.0	C30000
40000	42680.0	C40000
50000	53350.0	C50000
100000	106700.0	C100000

EXAMPLE 3

Dissolved 0.5 g glucose in 50ml distilled water. Aliquots from the stock solution were taken and diluted with water to achieve different concentrations in the range of 100 to 2000 mg/L (COD: 106.7 to 2134 mg/L). Weighed 0.02 g HgSO_4 in separate 15 mL and 25mL test tubes. 1.0 ml of glucose sample was added to each test tube and shake well. 0.5 ml of $\text{K}_2\text{Cr}_2\text{O}_7$ was added to each sample. 1.5 ml H_2SO_4 silver sulphate reagent was added slowly with continuous shaking. The color of the reaction mixture was noted down visually within 1 min of incubation and also recorded spectrophotometrically at 585 and 635 nm against air. The color varied with concentration. In 100 to 2000 mg/L glucose solutions, the color varied from C100 to C2000. At 100 to 800 mg/L glucose concentrations, OD 585 varied from 0.104 to 0.307 and OD 635 varied from 0.092 to 0.276. In the COD range of 106.7 to 640.2 mg/L the OD 585 and OD 635 there was an incremental difference of 0.030 and 0.027 for each 106.7 mg/L COD increase, respectively. At different wavelengths, OD became stable at 1000 mg/L glucose and higher concentrations. Hence for a preliminary estimate of COD value, color of the standard reaction mixture based on a range of glucose concentrations can be used as a reference. For samples, with initial OD concentration of 100 to 640 mg/L, OD 585 or OD 635 can be recorded for obtaining more values that are precise.

Table 3: Correlation between COD concentration, color of the reaction mixture and their OD at different wave lengths.

Glucose Concen tration (mg/L)	COD of Color Glucose solution (mg/L)	OD 585 nm of the reaction mixture	OD 635 nm of the reaction mixture
100	106.7	C100	0.104
200	213.4	C200	0.133
300	320.1	C300	0.164
400	426.8	C400	0.196
500	533.5	C500	0.222
600	640.2	C600	0.252
800	853.6	C800	0.307
1000	1067.0	C1000	0.362
1200	1280.4	C1200	0.347
1400	1493.8	C1400	0.351
1600	1707.2	C1600	0.354
1800	1920.6	C1800	0.344
2000	2134.0	C2000	0.358

EXAMPLE 4:

Dissolved 1 g glucose in 100mL distilled water. Aliquots from the stock solution were taken and diluted with water to achieve different concentrations in the range of 100 to 2000 mg/L (COD: 106.7 to 2134 mg/L). Weighed 0.02 g HgSO_4 in separate 15 mL and 25mL test tubes. 1.0 mL of glucose sample was added to each test tube and shake well. 0.5 ml of $\text{K}_2\text{Cr}_2\text{O}_7$ was added to each sample. 1.5 ml H_2SO_4 silver sulphate reagent was added slowly with continuous shaking. The color of the reaction mixture was noted down visually within 30 sec of incubation and also recorded spectrophotometrically at 585 and 635 nm against air. The color varied with concentration. In 100 to 2000 mg/L glucose solutions, the color varied from C100 to C2000. At 100 to 1000 mg/L glucose concentrations, OD 585 varied from 0.109 to 0.347 and OD 635 varied from 0.096 to 0.317. In the COD range of 106.7 to 640.2 mg/L the OD 585 and OD 635 there was an incremental difference of 0.027 and 0.026 for each 106.7 mg/L COD increase, respectively. At different wavelengths, OD became stable at 1000 mg/L glucose and higher concentrations. Hence, for a preliminary estimate of COD value, color of the standard reaction mixture based on a range of glucose concentrations can be used as a reference. For samples, with

initial COD concentration of 100 to 640 mg/L, OD 585 or OD 635 can be recorded for obtaining more precise values.

Table 4: Correlation between COD concentration, color of the reaction mixture and their OD at different wave lengths.

Glucose Concen- tration (mg/L)	COD of Color Glucose solution (mg/L)	of the reaction mixture	OD 585 nm of the reaction mixture	OD 635 nm of the reaction mixture
100	106.7	C100	0.109	0.096
200	213.4	C200	0.132	0.116
300	320.1	C300	0.161	0.142
400	426.8	C400	0.193	0.164
500	533.5	C500	0.222	0.195
600	640.2	C600	0.246	0.226
800	853.6	C800	0.300	0.274
1000	1067.0	C1000	0.347	0.317
1200	1280.4	C1200	0.350	0.321
1400	1493.8	C1400	0.341	0.326
1600	1707.2	C1600	0.352	0.332
1800	1920.6	C1800	0.354	0.331
2000	2134.0	C2000	0.354	0.349

EXAMPLE 5:

Dissolved 0.25 g glucose in 25 mL distilled water. Aliquots from the stock solution were taken and diluted with water to achieve different concentrations in the range of 100 to 500 mg/L (COD: 106.7 to 533.5 mg/L). Weighed 0.02 g HgSO_4 in separate 15 mL test tubes. 1.0 mL of glucose sample was added to each test tube and shake well. 0.5 mL of $\text{K}_2\text{Cr}_2\text{O}_7$ was added to each sample. 1.5 mL H_2SO_4 silver sulphate reagent was added slowly with continuous shaking. The color of the reaction mixture was noted down visually within 30 sec of incubation and also recorded spectrophotometrically at 585 and 635 nm against air. The color varied with concentration. In 100 to 500 mg/L glucose solutions, the color varied from C100 to C500. At 100 to 500 mg/L glucose concentrations, OD 585 varied from 0.106 to 0.223 and OD 635 varied from 0.089 to 0.199. In the COD range of 106.7 to 533.5 mg/L the OD 585 and OD 635 there was an incremental difference of 0.029 and

0.027 for each 106.7 mg/L COD increase, respectively. Hence, for a preliminary estimate of COD value, color of the standard reaction mixture based on a range of glucose concentrations can be used as a reference. For samples, with initial COD concentration of 106.7 to 533.5 mg/L, OD 585 or OD 635 can be recorded for obtaining more precise values.

Table 5: Correlation between COD concentration, color of the reaction mixture and their OD at different wave lengths.

Glucose Concen- tration (mg/L)	COD of Glucose solution (mg/L)	Color of the reaction mixture	OD 585 nm of the reaction mixture	OD 635 nm of the reaction mixture
100	106.7	C100	0.106	0.089
200	213.4	C200	0.133	0.116
300	320.1	C300	0.162	0.144
400	426.8	C400	0.192	0.166
500	533.5	C500	0.223	0.199

EXAMPLE 6:

Dissolved 0.2 g glucose in 50 mL distilled water. Aliquots from the stock solution were taken and diluted with water to achieve different concentrations in the range of 100 to 500 mg/L (COD: 106.7 to 533.5 mg/L). Weighed 0.02 g HgSO_4 in separate 15 mL test tubes. 1.0 ml of glucose sample was added to each test tube and shake well. 0.5 ml of $\text{K}_2\text{Cr}_2\text{O}_7$ was added to each sample. 1.5 mL H_2SO_4 silver sulphate reagent was added slowly with continuous shaking. The color of the reaction mixture was noted down visually within 15 sec of incubation and also recorded spectrophotometrically at 585 and 635 nm against air. The color varied with concentration. In 100 to 500 mg/L glucose solutions, the color varied from C100 to C500. At 100 to 500 mg/L glucose concentrations, OD 585 varied from 0.104 to 0.225 and OD 635 varied from 0.089 to 0.196. In the COD range of 106.7 to 533.5 mg/L the OD 585 and OD 635 there was an incremental difference of 0.030 and 0.028 for each 106.7 mg/L COD increase, respectively. Hence, for a preliminary estimate of COD value, color of the standard reaction mixture based on a range of glucose concentrations can be used as a reference. For samples, with initial COD

concentration of 106.7 to 533.5 mg/L, OD 585 or OD 635 can be recorded for obtaining more accurate values.

Table 6: Correlation between COD concentration, color of the reaction mixture and their OD at different wave lengths.

Glucose Concen tration (mg/L)	COD of Glucose solution (mg/L)	Color of the reaction mixture	OD 585 nm of the reaction mixture	OD 635 nm of the reaction mixture
100	106.7	C100	0.104	0.089
200	213.4	C200	0.132	0.115
300	320.1	C300	0.164	0.145
400	426.8	C400	0.196	0.165
500	533.5	C500	0.225	0.196

EXAMPLE 7:

Dissolved 0.3 g glucose in 50ml distilled water. Aliquots from the stock solution were taken and diluted with water to achieve different concentrations in the range of 300 to 500 mg/L (COD: 320.1 to 533.5 mg/L). Weighed 0.02 g HgSO_4 in separate 15 mL test tubes. 1.0 ml of glucose sample was added to each test tube and shake well. 0.5 mL of $\text{K}_2\text{Cr}_2\text{O}_7$ was added to each sample. 1.5 mL H_2SO_4 silver sulphate reagent was added slowly with continuous shaking. The color of the reaction mixture was noted down visually within 15 sec of incubation and also recorded spectrophotometrically at 585 against air. The color varied with concentration. In 300 to 500 mg/L glucose solutions, the color varied from C300 to C500. At 300 to 500 mg/L glucose concentrations, OD 585 varied on an average from 0.158 to 0.210. In the COD range of 320.1 to 533.5 mg/L the OD 585 there was an incremental difference of 0.026 for each 106.7 mg/L COD increase. Hence, for a preliminary estimate of COD value, color of the standard reaction mixture based on a range of glucose concentrations can be used as a reference. For samples, with initial COD concentration of 320.1 to 533.5 mg/L, OD 585 can be recorded for obtaining values that are more precise.

Table 7: Correlation between COD concentration, color of the reaction mixture and their OD at 585 nm.

S.	Glucose	COD of	Color	OD 585 nm
----	---------	--------	-------	-----------

No.	Concentration (mg/L)	Glucose solution (mg/L)	of the reaction mixture	of the reaction mixture
1.	300	320.1	C300	0.154
2.	300	320.1	C300	0.168
3.	300	320.1	C300	0.162
4.	300	320.1	C300	0.155
5.	300	320.1	C300	0.153
6.	300	320.1	C300	0.158
7.	300	320.1	C300	0.157
1.	400	426.8	C400	0.178
2.	400	426.8	C400	0.207
3.	400	426.8	C400	0.189
4.	400	426.8	C400	0.183
5.	400	426.8	C400	0.181
6.	400	426.8	C400	0.185
7.	400	426.8	C400	0.186
1.	500	533.5	C500	0.193
2.	500	533.5	C500	0.234
3.	500	533.5	C500	0.207
4.	500	533.5	C500	0.215
5.	500	533.5	C500	0.207
6.	500	533.5	C500	0.208
7.	500	533.5	C500	0.210

EXAMPLE 8:

Dissolved 0.5 g glucose in 100 mL distilled water. Aliquots from the stock solution were taken and diluted with water to achieve different COD concentrations in the range of 320 to 533 mg/L. A 20 mL sample of anaerobically digested damaged wheat grain slurry initially inoculated with *Aspergillus niger* for 10 days was taken. COD of the sample was estimated to be 52300 mg/L. (By Merck method and read on photometer.) Based on the estimated COD value, diluted samples were prepared in water to achieve a final COD in the range of 415 to 520 mg/L. Weighed

0.02 g HgSO_4 in separate 15 mL test tubes. 1.0 mL of glucose solution and sample were added to each test tube separately and mixed well. 0.5 mL of $\text{K}_2\text{Cr}_2\text{O}_7$ was added to each sample. 1.5 mL H_2SO_4 silver sulphate reagent was added slowly with continuous shaking. The color of the reaction mixture was noted down visually within 15 sec of incubation and also recorded spectrophotometrically at 585 and 635 nm against air. The color varied with concentration. In 320 to 533 mg/L glucose solutions, the color varied from C300 to C550. At 320 to 533 mg/L glucose concentrations, average OD 585 varied from 0.152 to 0.192 and OD 635 varied from 0.141 to 0.183. In the COD range of 320.1 to 533.5 mg/L the OD 585 and OD 635 there was an incremental difference of 0.026 and 0.024 for each 106.7 mg/L COD increase, respectively. On the other hand, color of the final sample reaction mixture varied from C400 to C500.

At 415.5 to 519.3 mg/L COD concentrations of the sample, average OD 585 varied from 0.145 to 0.169 and OD 635 varied from 0.129 to 0.148. In this COD range, there was an incremental difference of 0.024 and 0.020 for each 103.8 mg/L COD increase at OD 585 and OD 635, respectively.

Hence, for a preliminary estimate of COD value, color of the standard reaction mixture based on a range of glucose concentrations can be used as a reference. For samples, with COD concentration of 415.5 to 519.3 mg/L, OD 585 or OD 635 can be recorded for obtaining values that is more precise.

Table 8: Correlation between COD. concentration, color of the reaction mixture and their OD at different wave lengths.

COD (mg/L)	Color of the reaction mixture	OD of the reaction mixture	
		585 nm	635 nm
Glucose:			
320	C300	0.154	0.142
	C300	0.150	0.139
	C300	0.153	0.142
426	C400	0.178	0.166
	C400	0.175	0.162
	C400	0.181	0.169

533	C550	0.192	0.183
	C550	0.190	0.180
	C550	0.194	0.185

Sample: Damaged wheat grain slurry

415.5	C400	0.142	0.128
	C400	0.152	0.132
	C400	0.141	0.128
467.9	C450	0.146	0.128
	C450	0.149	0.135
	C450	0.147	0.128
519.3	C500	0.180	0.161
	C500	0.166	0.140
	C500	0.162	0.142

EXAMPLE 9:

Dissolved 0.15 g glucose in 25 mL distilled water. Aliquots from the stock solution were taken and diluted with water to achieve different COD concentrations in the range of 320 to 533 mg/L. A 30 mL sample of anaerobically digested damaged wheat grain slurry initially inoculated with *Bacillus licheniformis* was taken. COD of the sample was estimated to be 54700 mg/L. (By Merck method and read on photometer.) Based on the estimated COD value, diluted samples were prepared in water to achieve a final COD in the range of 300 to 500 mg/L. Weighed 0.02 g HgSO₄ in separate 15 mL test tubes. 1.0 mL of glucose solution and sample were added to each test tube separately and mixed well. 0.5 mL of K₂Cr₂O₇ was added to each sample. 1.5 mL H₂SO₄ silver sulphate reagent was added slowly with continuous shaking. The color of the reaction mixture was noted down visually within 15 sec of incubation and also recorded spectrophotometrically at 585 and 635 nm against air. The color varied with concentration. In 320 to 533 mg/L glucose solutions, the color varied from C300 to C550. At 320 to 533 mg/L glucose concentrations, average OD 585 varied from 0.168 to 0.234 and OD 635 varied from 0.151 to 0.213. In the COD range of 320.1 to 533.5 mg/L the OD 585 and OD 635

there was an incremental difference of 0.033 for each 106.7 mg/L COD increase. On the other hand, color of the final sample reaction mixture varied from C300 to C550. At 300 to 500 mg/L COD concentrations of the sample, average OD 585 varied from 0.136 to 0.175 and OD 635 varied from 0.123 to 0.156. In this COD range, there was an incremental difference of 0.020 and 0.017 for each 100 mg/L COD increase at OD 585 and OD 635, respectively.

Hence, for a preliminary estimate of COD value, color of the standard reaction mixture based on a range of glucose concentrations can be used as a reference. For samples, with COD concentration of 300 to 500 mg/L, OD 585 or OD 635 can be recorded for obtaining values that is more precise.

Table 9: Correlation between COD concentration, color of the reaction mixture and their OD at different wave lengths.

COD (mg/L)	Color of the reaction mixture	OD of the reaction mixture	
		585 nm	635 nm
Glucose:			
320	C300	0.164	0.146
	C300	0.167	0.152
	C300	0.172	0.154
426	C400	0.198	0.182
	C400	0.214	0.194
	C400	0.209	0.193
533	C550	0.235	0.212
	C550	0.236	0.217
	C550	0.231	0.211
Sample: Damaged wheat grain slurry			
300.7	C300	0.140	0.129
	C300	0.135	0.121
	C300	0.134	0.121
411.6	C400	0.158	0.146
	C400	0.156	0.140
	C400	0.157	0.141
500.6	C500	0.174	0.155



C500
C500



0.175 0.160
0.175 0.155

EXAMPLE 10:

Dissolved 0.2-g glucose in 50 mL distilled water. Aliquots from the stock solution were taken and diluted with water to achieve different COD concentrations in the range of 320 to 533 mg/L. A 25-mL sample of anaerobically digested damaged wheat grain slurry was taken. COD of the sample was estimated to be 61300 mg/L. (By Merck method and read on photometer.) Based on the estimated COD value, diluted samples were prepared in water to achieve a final COD in the range of 335 to 505 mg/L. Weighed 0.02 g HgSO_4 in separate 15-mL test tubes. 1.0 mL of glucose solution and sample were added to each test tube separately and mixed well. 0.5 mL of $\text{K}_2\text{Cr}_2\text{O}_7$ was added to each sample. 1.5 mL H_2SO_4 silver sulphate reagent was added slowly with continuous shaking. The color of the reaction mixture was noted down visually within 15 sec of incubation and recorded spectrophotometrically at 585 and 635 nm against air. The color varied with concentration. In 320 to 533 mg/L glucose solutions, the color varied from C300 to C550. At 320 to 533 mg/L glucose concentrations, average OD 585 varied from 0.162 to 0.207 and OD 635 varied from 0.145 to 0.186. In the COD range of 320.1 to 533.5 mg/L, the OD 585 and OD 635 there was an incremental difference of 0.023 and 0.025 for each 106.7 mg/L COD increase, respectively. On the other hand, color of the final sample reaction mixture varied from C350 to C500.

At 335 to 505 mg/L COD concentrations of the sample, average OD 585 varied from 0.132 to 0.169 and OD 635 varied from 0.116 to 0.151. In this COD range, there was an incremental difference of 0.019 and 0.022 for each 85 mg/L COD increase at OD 585 and OD 635, respectively.

Hence, for a preliminary estimate of COD value, color of the standard reaction mixture based on a range of glucose concentrations can be used as a reference. For samples, with COD concentration of 335 to 505 mg/L, OD 585 or OD 635 can be recorded for obtaining values that is more precise.

Table 10: Correlation between COD concentration, color of the reaction mixture and their OD at different wavelengths.

COD (mg/L)	Color of the reaction mixture	OD of the reaction mixture	
		585 nm	635 nm
Glucose:			
320	C300	0.160	0.143
	C300	0.164	0.147
	C300	0.162	0.145
426	C400	0.190	0.184
	C400	0.191	0.176
	C400	0.187	0.168
533	C550	0.207	0.190
	C550	0.206	0.184
	C550	0.207	0.185
Sample: Damaged wheat grain slurry			
337.1	C350	0.128	0.112
	C350	0.134	0.119
	C350	0.135	0.117
421.3	C400	0.153	0.135
	C400	0.151	0.134
	C400	0.155	0.137
504.6	C500	0.168	0.149
	C500	0.169	0.149
	C500	0.169	0.154

EXAMPLE 11:

Dissolved 0.1 g glucose in 50 mL distilled water. Aliquots from the stock solution were taken and diluted with water to achieve different COD concentrations in the range of 320 to 533 mg/L. A 40-mL sample of anaerobically digested damaged wheat grain slurry initially treated with *Aspergillus niger* for 10 days and inoculated with *Bacillus subtilis* was taken. COD of the sample was estimated to be 6600 mg/L. (By Merck method and read on photometer.) Based on the estimated COD value, diluted samples were prepared in water to achieve a final COD in the range of 330 to

530 mg/L. Weighed 0.02 g HgSO_4 in separate 15-mL test tubes. 1.0 mL of glucose solution and sample were added to each test tube separately and mixed well. 0.5 mL of $\text{K}_2\text{Cr}_2\text{O}_7$ was added to each sample. 1.5 mL H_2SO_4 silver sulfate reagent was added slowly with continuous shaking. The color of the reaction mixture was noted down visually within 15 sec of incubation and recorded spectrophotometrically at 585 and 635 nm against air. The color varied with concentration. In 320 to 533 mg/L glucose solutions, the color varied from C300 to C550. At 320 to 533 mg/L glucose concentrations, average OD 585 varied from 0.155 to 0.215 and OD 635 varied from 0.144 to 0.194. In the COD range of 320.1 to 533.5 mg/L, the OD 585 and OD 635 there was an incremental difference of 0.030 and 0.025 for each 106.7-mg/L COD increase, respectively. On the other hand, color of the final sample reaction mixture varied from C350 to C550.

At 330 to 530 mg/L COD concentrations of the sample, average OD 585 varied from 0.137 to 0.174 and OD 635 varied from 0.121 to 0.153. In this COD range, there was an incremental difference of 0.019 and 0.016 for each 100-mg/L COD increase at OD 585 and OD 635, respectively.

Hence, for a preliminary estimate of COD value, color of the standard reaction mixture based on a range of glucose concentrations can be used as a reference. For samples, with COD concentration of 330 to 530 mg/L, OD 585 or OD 635 can be recorded for obtaining values that is more precise.

Table 11: Correlation between COD concentration, color of the reaction mixture and their OD at different wavelengths.

COD (mg/L)	Color of the reaction mixture	OD of the reaction mixture	
		585 nm	635 nm
Glucose:			
320	C300	0.154	0.143
	C300	0.155	0.142
	C300	0.157	0.146
426	C400	0.182	0.170
	C400	0.182	0.168
	C400	0.184	0.170

533	C550	0.206	0.187
	C550	0.208	0.190
	C550	0.230	0.205

Sample: Damaged wheat grain slurry

331	C350	0.134	0.118
	C350	0.138	0.126
	C350	0.138	0.120
429	C400	0.157	0.139
	C400	0.153	0.135
	C400	0.154	0.134
529	C500	0.172	0.152
	C500	0.171	0.151
	C500	0.178	0.156

EXAMPLE 12:

Dissolved 0.2 g glucose in 50 mL distilled water. Aliquots from the stock solution were taken and diluted with water to achieve different COD concentrations in the range of 320 to 533 mg/L. A 30 mL sample of partially digested vegetable waste slurry was taken. COD of the sample was estimated to be 9640 mg/L (By Merck method and read on photometer). Based on the estimated COD value, diluted samples were prepared in water to achieve a final COD in the range of 290 to 480 mg/L. Weighed 0.02 g HgSO_4 in separate 15-mL test tubes. 1.0 mL of glucose solution and sample were added to each test tube separately and mixed well. 0.5 mL of $\text{K}_2\text{Cr}_2\text{O}_7$ was added to each sample. 1.5 mL H_2SO_4 silver sulfate reagent was added slowly with continuous shaking. The color of the reaction mixture was noted down visually within 15 sec of incubation and also recorded spectrophotometrically at 585 and 635 nm against air. The color varied with concentration. In 320 to 533 mg/L glucose solutions, the color varied from C300 to C550. At 320 to 533 mg/L glucose concentrations, average OD 585 varied from 0.153 to 0.207 and OD 635 varied from 0.142 to 0.188. In the COD range of 320.1 to 533.5 mg/L the OD 585 and OD 635 there was an incremental difference of 0.027 and 0.023 for each 106.7 mg/L COD increase, respectively. On the other hand, color of the final sample reaction mixture varied from C300 to C500.

At 290 to 480 mg/L COD concentrations of the sample, average OD 585 varied from 0.138 to 0.157 and OD 635 varied from 0.123 to 0.141. In this COD range, there was an incremental difference of 0.009 for each 95 mg/L COD increase at OD 585 and OD 635.

Hence, for a preliminary estimate of COD value, color of the standard reaction mixture based on a range of glucose concentrations can be used as a reference. For samples, with COD concentration of 290 to 480 mg/L, OD 585 or OD 635 can be recorded for obtaining values that is more accurate.

Table 12: Correlation between COD concentration, color of the reaction mixture and their OD at different wave lengths.

COD (mg/L)	Color of the reaction mixture	OD of the reaction mixture	
		585 nm	635 nm
Glucose:			
320	C300	0.152	0.142
	C300	0.156	0.145
	C300	0.151	0.139
426	C400	0.181	0.169
	C400	0.183	0.170
	C400	0.180	0.168
533	C550	0.205	0.186
	C550	0.212	0.194
	C550	0.204	0.184
Sample: Vegetable waste slurry			
290	C300	0.140	0.123
	C300	0.132	0.120
	C300	0.141	0.127
386	C400	0.158	0.143
	C400	0.151	0.140
	C400	0.144	0.128
482	C500	0.162	0.146
	C500	0.156	0.141
	C500	0.153	0.137

EXAMPLE 13:

Dissolved 0.25 g glucose in 50 mL distilled water. Aliquots from the stock solution were taken and diluted with water to achieve different COD concentrations in the range of 320 to 533 mg/L. A 50 mL sample of partially digested vegetable waste slurry was mixed with anaerobically digested damaged wheat grain slurry. COD of the sample was estimated to be 9400 mg/L (By Merck method and read on photometer). Based on the estimated COD value, diluted samples were prepared in water to achieve a final COD in the range of 297 to 533 mg/L. Weighed 0.02 g HgSO_4 in separate 15 mL test tubes. 1.0 mL of glucose solution and sample were added to each test tube separately and mixed well. 0.5 mL of $\text{K}_2\text{Cr}_2\text{O}_7$ was added to each sample. 1.5 mL H_2SO_4 silver sulphate reagent was added slowly with continuous shaking. The color of the reaction mixture was noted down visually within 15 sec of incubation and also recorded spectrophotometrically at 585 and 635 nm against air. The color varied with concentration. In 320 to 533 mg/L glucose solutions, the color varied from C300 to C500. At 320 to 533 mg/L glucose concentrations, average OD 585 varied from 0.158 to 0.208 and OD 635 varied from 0.138 to 0.187. In the COD range of 320.1 to 533.5 mg/L the OD 585 and OD 635 there was an incremental difference of 0.025 and 0.025 for each 106.7 mg/L COD increase, respectively. On the other hand, color of the final sample reaction mixture varied from C300 to C500.

At 297 to 533 mg/L COD concentrations of the sample, average OD 585 varied from 0.132 to 0.179 and OD 635 varied from 0.118 to 0.162. In this COD range, there was an incremental difference of 0.023 and 0.027 for each 118 mg/L COD increase at OD 585 and OD 635, respectively.

Hence, for a preliminary estimate of COD value, color of the standard reaction mixture based on a range of glucose concentrations can be used as a reference. For samples, with COD concentration of 297 to 533 mg/L, OD 585 or OD 635 can be recorded for obtaining more precise values.

Table 13: Correlation between COD concentration, color of the reaction mixture and their OD at different wave lengths.

COD (mg/L)	Color of the reaction mixture	OD of the reaction mixture	
		585 nm	635 nm

 Glucose:

320	C300	0.167	0.148
	C300	0.154	0.135
	C300	0.152	0.131
426	C400	0.188	0.167
	C400	0.185	0.162
	C400	0.181	0.159
533	C550	0.215	0.195
	C550	0.206	0.184
	C550	0.202	0.181

Sample: Vegetable waste slurry + Damaged wheat grain slurry

297	C300	0.126	0.111
	C300	0.138	0.124
	C300	0.133	0.118
445	C450	0.168	0.152
	C450	0.164	0.148
	C450	0.168	0.152
533	C550	0.188	0.169
	C550	0.180	0.163
	C550	0.169	0.153

EXAMPLE 14:

Dissolved 0.15 g glucose in 50 mL distilled water. Aliquots from the stock solution were taken and diluted with water to achieve different COD concentrations in the range of 320 to 533 mg/L. A 30 mL sample of apple pomace soaked for 1 day was taken. COD of the sample was estimated to be 51700 mg/L (By Merck method and read on photometer). Based on the estimated COD value, diluted samples were prepared in water to achieve a final COD in the range of 310 to 517 mg/L. Weighed 0.02 g HgSO_4 in separate 15 mL test tubes. 1.0 mL of glucose solution and sample were added to each test tube separately and mixed well. 0.5 ml of $\text{K}_2\text{Cr}_2\text{O}_7$ was added to each sample. 1.5 mL H_2SO_4 silver sulphate reagent was added slowly with continuous shaking. The color of the reaction mixture was noted down visually within 15 sec of incubation and also recorded spectrophotometrically at 585 and 635 nm

against air. The color varied with concentration. In 320 to 533 mg/L glucose solutions, the color varied from C300 to C550. At 320 to 533 mg/L glucose concentrations, average OD 585 varied from 0.157 to 0.210 and OD 635 varied from 0.142 to 0.189. In the COD range of 320.1 to 533.5 mg/L the OD 585 and OD 635 there was an incremental difference of 0.026 and 0.024 for each 106.7 mg/L COD increase, respectively. On the other hand, color of the final sample reaction mixture varied from C300 to C500.

At 310 to 517 mg/L COD concentrations of the sample, average OD 585 varied from 0.141 to 0.186 and OD 635 varied from 0.126 to 0.167. In this COD range, there was an incremental difference of 0.022 and 0.020 for each 104 mg/L COD increase at OD 585 and OD 635, respectively.

Hence for a preliminary estimate of COD value, color of the standard reaction mixture based on a range of glucose concentrations can be used as a reference. For samples, with COD concentration of 310 to 517 mg/L, OD 585 or OD 635 can be recorded for obtaining more precise values.

Table 14: Correlation between COD concentration, color of the reaction mixture and their OD at different wave lengths.

COD (mg/L)	Color of the reaction mixture	OD of the reaction mixture	
		585 nm	635 nm
Glucose:			
320	C300	0.157	0.142
	C300	0.154	0.139
	C300	0.159	0.146
426	C400	0.186	0.168
	C400	0.185	0.162
	C400	0.187	0.170
533	C550	0.216	0.196
	C550	0.209	0.188
	C550	0.207	0.183
Sample: Apple pomace			
310	C300	0.140	0.126
	C300	0.142	0.126
	C300	0.140	0.126
413	C400	0.171	0.152
	C400	0.165	0.148
	C400	0.165	0.147
517	C500	0.186	0.166
	C500	0.187	0.168
	C500	0.186	0.167

EXAMPLE 15:

Dissolved 0.1 g glucose in 50 mL distilled water. Aliquots from the stock solution were taken and diluted with water to achieve different concentrations in the range of 320 to 533 mg/L. A 50 mL sample was prepared by mixing apple pomace soaked for 1 day with digested damaged wheat grains and partially digested vegetable market waste. COD of the sample was estimated to be 40300 mg/L (By Merck method and read on photometer). Based on the estimated COD value, diluted samples were prepared in water to achieve a final COD in the range of 310 to 497 mg/L. Weighed 0.02 g HgSO_4 in separate 15 mL test tubes. 1.0 mL of glucose solution and sample

were added to each test tube separately and mixed well. 0.5 ml of $K_2Cr_2O_7$ was added to each sample. 1.5 mL H_2SO_4 silver sulphate reagent was added slowly with continuous shaking. The color of the reaction mixture was noted down visually within 15 sec of incubation and also recorded spectrophotometrically at 585 and 635 nm against air. The color varied with concentration. In 320 to 533 mg/L glucose solutions, the color varied from C300 to C550. At 320 to 533 mg/L glucose concentrations, average OD 585 varied from 0.159 to 0.212 and OD 635 varied from 0.142 to 0.190. In the COD range of 320.1 to 533.5 mg/L the OD 585 and OD 635 there was an incremental difference of 0.026 and 0.024 for each 106.7 mg/L COD increase, respectively. On the other hand, color of the final sample reaction mixture varied from C300 to C500.

At 310 to 497 mg/L COD concentrations of the sample, average OD 585 varied from 0.125 to 0.169 and OD 635 varied from 0.112 to 0.155. In this COD range, there was an incremental difference of 0.022 and 0.021 for each 93 mg/L COD increase at OD 585 and OD 635, respectively.

Hence for a preliminary estimate of COD value, color of the standard reaction mixture based on a range of glucose concentrations can be used as a reference. For samples, with COD concentration of 310 to 497 mg/L, OD 585 or OD 635 can be recorded for obtaining more precise values.

Table 15: Correlation between COD concentration, color of the reaction mixture and their OD at different wave lengths.

COD (mg/L)	Color of the reaction mixture	OD of the reaction mixture	
		585 nm	635 nm
Glucose:			
320	C300	0.160	0.144
	C300	0.157	0.141
	C300	0.159	0.142
426	C400	0.187	0.170
	C400	0.185	0.167
	C400	0.189	0.172
535	C550	0.220	0.199
	C550	0.209	0.187
	C550	0.207	0.184

Sample : Apple pomace+ Vegetable waste slurry + Damaged wheat grain slurry

310	C300	0.126	0.114
	C300	0.126	0.112
	C300	0.124	0.111
404	C400	0.148	0.136
	C400	0.154	0.146
	C400	0.152	0.141
497	C500	0.171	0.159
	C500	0.168	0.154
	C500	0.167	0.152

EXAMPLE 16:

Dissolved 0.15 g glucose in 25 mL distilled water. Aliquots from the stock solution were taken and diluted with water to achieve different COD concentrations in the range of 320 to 535 mg/L. A 30 mL sample of partially digested vegetable waste slurry mixed with potassium hydrogen phosphate salt was taken. Weighed 0.02 g HgSO_4 in separate 15 mL test tubes. 1.0 mL of glucose solution and sample were added to each test tube separately and mixed well. 0.5 mL of $\text{K}_2\text{Cr}_2\text{O}_7$ was added to each sample. 1.5 mL H_2SO_4 silver sulphate reagent was added slowly with continuous shaking. The color of the reaction mixture was noted down visually within 15 sec of incubation and also recorded spectrophotometrically at 585 and 635 nm against air. An approximate COD value of the sample was estimated from the color of standard reaction mixture prepared with glucose as reference material in the COD range of 200 to 100000 mg/L. OD 585 and OD 635 of the reaction mixture was read spectrophotometrically OD 585 and 635 nm in the range of 0.130 to 0.175 and 0.120 to 0.160 respectively were considered and calculated with reference to glucose standard. (Glucose COD in the range of 320 to 535 mg/L corresponds to OD 585 and 635 nm in the range of 0.157 to 0.207 and 0.142 to 0.189, respectively). COD of the sample was estimated to be 777 mg/L. (By Merck method and read on photometer).

Based on the color of the reaction mixture, COD of the sample was estimated to be around 750 mg/L. On the basis of the dilution method and the spectrophotometric

reading at OD 585 and OD 635 of the sample, COD was estimated to be 824 and 852 mg/L, respectively. A comparison with the COD value calculated through Merck method reveals the following relation i.e. Our COD value (at OD 585) x 0.94 = COD value (Merck Method) or our COD value (at OD 635) x 0.91 = COD value (Merck Method).

Hence for a preliminary estimate of COD value, color of the standard reaction mixture based on a range of glucose concentrations can be used as a reference. For samples, with COD concentration of 300 to 500 mg/L, OD 585 or OD 635 can be recorded for obtaining more precise values.

Table 16: Correlation between COD concentration, color of the reaction mixture and their OD at different wave lengths.

COD (mg/L)	Sample dilution (Times)	Color of the reaction mixture	OD	Estimated COD (mg/L)	Estimated final COD (mg/L) (Col. 2x5)
1	2	3	4	5	6
<u>Values at OD 585 nm</u>					
Glucose:					
320	NA	C300	0.162	NA	NA
426	NA	C400	0.186	NA	NA
535	NA	C550	0.206	NA	NA
Sample: Vegetable waste (filtrate) + Potassium hydrogen phosphate					
	Nil	C1500	0.240	≅1500#	≅1500
	5	C150	0.101	150#	750
	2	C400	0.150	412	824
<u>Values at OD 635 nm</u>					
Glucose:					
320	NA	C300	0.150	NA	NA
426	NA	C400	0.167	NA	NA
535	NA	C550	0.187	NA	NA
Sample: Vegetable waste (filtrate) + Potassium hydrogen phosphate					
	Nil	C1500	0.40	≅1500#	≅1500
	5	C150	0.101	150#	750
	2	C400	0.150	426	852

#: COD as deducible from Table 2 given in Example 2. D: Dark. T: Turbid.

*: OD out of range. NA: Not applicable. Samples are diluted to get an OD 585 in the range of 0.130 to 0.175 and OD 635 in the range of 0.120 to 0.160, to obtain a COD value in the range of 320 to 535 mg/L. (Based on the Tables 3 to 15 given in Examples 3 to 15).

EXAMPLE 17:

Dissolved 0.15 g glucose in 25 mL distilled water. Aliquots from the stock solution were taken and diluted with water to achieve different COD concentrations in the range of 320 to 535 mg/L. A 30 mL sample of partially digested palm oil mill effluent (filtrate) was taken. Weighed 0.02 g HgSO_4 in separate 15 mL test tubes. 1.0 mL of glucose solution and sample were added to each test tube separately and mixed well. 0.5 mL of $\text{K}_2\text{Cr}_2\text{O}_7$ was added to each sample. 1.5 mL H_2SO_4 silver sulphate reagent was added slowly with continuous shaking. The color of the reaction mixture was noted down visually within 15 sec of incubation and also recorded spectrophotometrically at 585 and 635 nm against air. An approximate COD value of the sample was estimated from the color of standard reaction mixture prepared with glucose as reference material in the COD range of 200 to 100000 mg/L. OD 585 and OD 635 of the reaction mixture was read spectrophotometrically. OD 585 and 635 nm in the range of 0.130 to 0.175 and 0.120 to 0.160 respectively were considered and calculated with reference to glucose standard. (Glucose COD in the range of 320 to 535 mg/L corresponds to OD 585 and 635 nm in the range of 0.159 to 0.205 and 0.149 to 0.183, respectively).

COD of the sample was estimated to be 687 mg/L. (By Merck method and read on photometer).

Based on the color of the reaction mixture, COD of the sample was estimated to be around 500 mg/L. On the basis of the dilution method and the spectrophotometric reading at OD 585 and OD 635 of the sample, COD was estimated to be 874 and 882 mg/L, respectively. A comparison with the COD value calculated through Merck method reveals the following relation i.e. Our COD value (at OD 585) $\times 0.79 =$ COD value (Merck Method) or our COD value (at OD 635) $\times 0.79 =$ COD value (Merck Method).

Hence for a preliminary estimate of COD value, color of the standard reaction mixture based on a range of glucose concentrations can be used as a reference. For samples, with COD concentration of 300 to 500 mg/L, OD 585 or OD 635 can be recorded for obtaining more precise values.

Table 17: Correlation between COD concentration, color of the reaction mixture and their OD at different wave lengths.

COD (mg/L)	Sample dilution (Times)	Color of the reaction mixture	OD	Estimated COD (mg/L)	Estimated final COD (mg/L) (Col. 2x5)
1	2	3	4	5	6
<u>Values at OD 585 nm</u>					
Glucose:					
320	NA	C300	0.159	NA	NA
426	NA	C400	0.186	NA	NA
535	NA	C550	0.205	NA	NA
Sample: Palm oil mill effluent (filtrate)					
	Nil	C500	0.236	≈500#	≈500
	2	C450	0.159	437	874
<u>Values at OD 635 nm</u>					
Glucose:					
320	NA	C300	0.149	NA	NA
426	NA	C400	0.166	NA	NA
535	NA	C550	0.183	NA	NA
Sample: Palm oil mill effluent (filtrate)					
	Nil	C500	0.213	≈500#	≈500
	2	C450	0.143	441	882

#: COD as deducible from Table 2 given in Example 2. D: Dark. T: Turbid.

*: OD out of range. NA: Not applicable. Samples are diluted to get an OD 585 in the range of 0.130 to 0.175 and OD 635 in the range of 0.120 to 0.160, to obtain a COD value in the range of 320 to 535 mg/L. (Based on the Tables 3 to 15 given in Examples 3 to 15).

EXAMPLE 18:

Dissolved 0.15 g glucose in 25 mL distilled water. Aliquots from the stock solution were taken and diluted with water to achieve different COD concentrations in the range of 320 to 535 mg/L. A 30 mL sample of partially digested vegetable waste slurry mixed was taken. Weighed 0.02 g HgSO_4 in separate 15 mL test tubes. 1.0 mL of glucose solution and sample were added to each test tube separately and mixed well. 0.5 mL of $\text{K}_2\text{Cr}_2\text{O}_7$ was added to each sample. 1.5 mL H_2SO_4 silver sulphate reagent was added slowly with continuous shaking. The color of the reaction mixture was noted down visually within 15 sec of incubation and also recorded spectrophotometrically at 585 and 635 nm against air. An approximate COD value of the sample was estimated from the color of standard reaction mixture prepared with glucose as reference material in the COD range of 200 to 100000 mg/L. OD 585 and OD 635 of the reaction mixture was read spectrophotometrically. OD 585 and 635 nm in the range of 0.130 to 0.175 and 0.120 to 0.160 respectively were considered and calculated with reference to glucose standard. (Glucose COD in the range of 320 to 535 mg/L corresponds to OD 585 and 635 nm in the range of 0.157 to 0.207 and 0.142 to 0.189, respectively).

COD of the sample was estimated to be 1930 mg/L. (By Merck method and read on photometer).

Based on the color of the reaction mixture, COD of the sample was estimated to be around 1400 mg/L. On the basis of the dilution method and the spectrophotometric reading at OD 585 and OD 635 of the sample, COD was estimated to be 1280 and 1288 mg/L, respectively. A comparison with the COD value calculated through Merck method reveals the following relation i.e. Our COD value (at OD 585) $\times 1.51 =$ COD value (Merck Method) or our COD value (at OD 635) $\times 1.50 =$ COD value (Merck Method).

Hence for a preliminary estimate of COD value, color of the standard reaction mixture based on a range of glucose concentrations can be used as a reference. For samples, with COD concentration of 300 to 500 mg/L, OD 585 or OD 635 can be recorded for obtaining more precise values.

Table 18: Correlation between COD concentration, color of the reaction mixture and their OD at different wave lengths.

COD (mg/L)	Sample dilution (Times)	Color of the reaction mixture	OD	Estimated COD (mg/L)	Estimated final COD (mg/L) (Col. 2x5)
1	2	3	4	5	6

Values at OD 585 nm

Glucose:

320	NA	C300	0.161	NA	NA
426	NA	C400	0.187	NA	NA
535	NA	C550	0.202	NA	NA

Sample: Vegetable waste slurry

Nil	C1400	0.294	≅ 1400#	≅ 1400
4	C320	0.134	320	1280

Values at OD 635 nm

Glucose:

320	NA	C300	0.144	NA	NA
426	NA	C400	0.167	NA	NA
535	NA	C550	0.184	NA	NA

Sample: Vegetable waste slurry

Nil	C1400	0.268	≅ 1400#	≅ 1400
4	C320	0.121	322	1288

#: COD as deducible from Table 2 given in Example 2. D: Dark. T: Turbid.

*: OD out of range. NA: Not applicable. Samples are diluted to get an OD 585 in the range of 0.130 to 0.175 and OD 635 in the range of 0.120 to 0.160, to obtain a COD value in the range of 320 to 535 mg/L. (Based on the Tables 4 to 15 given in Examples 4 to 15.

EXAMPLE 19:

Dissolved 0.15 g glucose in 25 mL distilled water. Aliquots from the stock solution were taken and diluted with water to achieve different COD concentrations in the range of 320 to 535 mg/L. A 30 mL sample of partially digested Vegetable waste

slurry (filtrate) mixed was taken. Weighed 0.02 g HgSO_4 in separate 15 mL test tubes. 1.0 mL of glucose solution and sample were added to each test tube separately and mixed well. 0.5 mL of $\text{K}_2\text{Cr}_2\text{O}_7$ was added to each sample. 1.5 mL H_2SO_4 silver sulphate reagent was added slowly with continuous shaking. The color of the reaction mixture was noted down visually within 15 sec of incubation and also recorded spectrophotometrically at 585 and 635 nm against air. An approximate COD value of the sample was estimated from the color of standard reaction mixture prepared with glucose as reference material in the COD range of 200 to 100000 mg/L. OD 585 and OD 635 of the reaction mixture was read spectrophotometrically. OD 585 and 635 nm in the range of 0.130 to 0.175 and 0.120 to 0.160 respectively were considered and calculated with reference to glucose standard. (Glucose COD in the range of 320 to 535 mg/L corresponds to OD 585 and 635 nm in the range of 0.157 to 0.207 and 0.142 to 0.189, respectively).

COD of the sample was estimated to be 1033 mg/L. (By Merck method and read on photometer).

Based on the color of the reaction mixture, COD of the sample was estimated to be around 1500 mg/L. On the basis of the dilution method and the spectrophotometric reading at OD 585 and OD 635 of the sample, COD was estimated to be 1632 and 1652 mg/L, respectively. A comparison with the COD value calculated through Merck method reveals the following relation i.e. Our COD value (at OD 585) $\times 0.63 =$ COD value (Merck Method) or our COD value (at OD 635) $\times 0.62 =$ COD value (Merck Method).

Hence for a preliminary estimate of COD value, color of the standard reaction mixture based on a range of glucose concentrations can be used as a reference. For samples, with COD concentration of 300 to 500 mg/L, OD 585 or OD 635 can be recorded for obtaining more precise values.

Table 19: Correlation between COD concentration, color of the reaction mixture and their OD at different wave lengths.

COD (mg/L)	Sample dilution (Times)	Color of the reaction mixture	OD	Estimated COD (mg/L)	Estimated final COD (mg/L) (Col. 2x5)
1	2	3	4	5	6
<u>Values at OD 585 nm</u>					
Glucose:					
320	NA	C300	0.159	NA	NA
426	NA	C400	0.186	NA	NA
535	NA	C550	0.205	NA	NA
Sample: Vegetable waste slurry (filtrate)					
	Nil	C1500	0.257	≅1500#	≅1500
	4	C400	0.148	408	1632
<u>Values at OD 635 nm</u>					
Glucose:					
320	NA	C300	0.145	NA	NA
426	NA	C400	0.166	NA	NA
535	NA	C550	0.187	NA	NA
Sample: Vegetable waste slurry (filtrate)					
	Nil	C1500	0.236	≅1500#	≅1500
	4	C400	0.134	413	1652

#: COD as deducible from Table 2 given in Example 2. D: Dark. T: Turbid.

*: OD out of range. NA: Not applicable. Samples are diluted to get an OD 585 in the range of 0.130 to 0.175 and OD 635 in the range of 0.120 to 0.160, to obtain a COD value in the range of 320 to 535 mg/L. (Based on the Tables 3 to 15 given in Examples 3 to 15).

EXAMPLE 20:

Dissolved 0.15 g glucose in 25 mL distilled water. Aliquots from the stock solution were taken and diluted with water to achieve different COD concentrations in the range of 320 to 535 mg/L. A 30 mL sample of Damaged wheat grain (filtrate) mixed with vegetable waste slurry (filtrate) was taken. Weighed 0.02 g HgSO₄ in separate

15 mL test tubes. 1.0 mL of glucose solution and sample were added to each test tube separately and mixed well. 0.5 mL of $K_2Cr_2O_7$ was added to each sample. 1.5 mL H_2SO_4 silver sulphate reagent was added slowly with continuous shaking. The color of the reaction mixture was noted down visually within 15 sec of incubation and also recorded spectrophotometrically at 585 and 635 nm against air. An approximate COD value of the sample was estimated from the color of standard reaction mixture prepared with glucose as reference material in the COD range of 200 to 100000 mg/L. OD 585 and OD 635 of the reaction mixture was read spectrophotometrically. OD 585 and 635 nm in the range of 0.130 to 0.175 and 0.120 to 0.160 respectively were considered and calculated with reference to glucose standard. (Glucose COD in the range of 320 to 535 mg/L corresponds to OD 585 and 635 nm in the range of 0.157 to 0.207 and 0.142 to 0.189, respectively).

COD of the sample was estimated to be 4067 mg/L. (By Merck method and read on photometer).

Based on the color of the reaction mixture, COD of the sample was estimated to be around 3000 mg/L. On the basis of the dilution method and the spectrophotometric reading at OD 585 and OD 635 of the sample, COD was estimated to be 4170 and 4180 mg/L, respectively. A comparison with the COD value calculated through Merck method reveals the following relation i.e. Our COD value (at OD 585) $\times 0.97 =$ COD value (Merck Method) or our COD value (at OD 635) $\times 0.97 =$ COD value (Merck Method).

Hence for a preliminary estimate of COD value, color of the standard reaction mixture based on a range of glucose concentrations can be used as a reference. For samples, with COD concentration of 300 to 500 mg/L, OD 585 or OD 635 can be recorded for obtaining more precise values.

Table 20: Correlation between COD concentration, color of the reaction mixture and their OD at different wave lengths.

COD (mg/L)	Sample dilution (Times)	Color of the reaction mixture	OD	Estimated COD (mg/L)	Estimated final COD (mg/L) (Col. 2x5)
1	2	3	4	5	6

Values at OD 585 nm

Glucose:

320	NA	C300	0.161	NA	NA
426	NA	C400	0.186	NA	NA
535	NA	C550	0.206	NA	NA

Sample: Damaged wheat grain (filtrate) + Vegetable waste slurry (filtrate)

Nil	C1400	1.031	≅1400#	≅ 1400
5	C600	0.234	≅ 600#	≅ 3000
10	C400	0.152	417	4170

Values at OD 635 nm

Glucose:

320	NA	C300	0.146	NA	NA
426	NA	C400	0.167	NA	NA
535	NA	C550	0.186	NA	NA

Sample: Damaged wheat grain (filtrate) + Vegetable waste slurry (filtrate)

Nil	C1400	0.784	≅1400#	≅ 1400
5	C600	0.213	≅ 600#	≅ 3000
10	C400	0.137	418	4180

#: COD as deducible from Table 2 given in Example 2. D: Dark. T: Turbid.

*: OD out of range. NA: Not applicable. Samples are diluted to get an OD 585 in the range of 0.130 to 0.175 and OD 635 in the range of 0.120 to 0.160, to obtain a COD value in the range of 320 to 535 mg/L. (Based on the Tables 3 to 15 given in Examples 3 to 15).

EXAMPLE 21:

Dissolved 0.15 g glucose in 25 mL distilled water. Aliquots from the stock solution were taken and diluted with water to achieve different COD concentrations in the

range of 320 to 535 mg/L. A 30 mL sample of Damaged wheat grain (filtrate) mixed with vegetable waste slurry was taken. Weighed 0.02 g HgSO_4 in separate 15 mL test tubes. 1.0 mL of glucose solution and sample were added to each test tube separately and mixed well. 0.5 ml of $\text{K}_2\text{Cr}_2\text{O}_7$ was added to each sample. 1.5 mL H_2SO_4 silver sulphate reagent was added slowly with continuous shaking. The color of the reaction mixture was noted down visually within 15 sec of incubation and also recorded spectrophotometrically at 585 and 635 nm against air. An approximate COD value of the sample was estimated from the color of standard reaction mixture prepared with glucose as reference material in the COD range of 200 to 100000 mg/L. OD 585 and OD 635 of the reaction mixture was read spectrophotometrically. OD 585 and 635 nm in the range of 0.130 to 0.175 and 0.120 to 0.160 respectively were considered and calculated with reference to glucose standard. (Glucose COD in the range of 320 to 535 mg/l corresponds to OD 585 and 635 nm in the range of 0.157 to 0.207 and 0.142 to 0.189, respectively).

COD of the sample was estimated to be 4059 mg/L. (By Merck method and read on photometer).

Based on the color of the reaction mixture, COD of the sample was estimated to be around 2900 mg/L. On the basis of the dilution method and the spectrophotometric reading at OD 585 and OD 635 of the sample, COD was estimated to be 4440 and 4490 mg/L, respectively. A comparison with the COD value calculated through Merck method reveals the following relation i.e. Our COD value (at OD 585) $\times 0.91 =$ COD value (Merck Method) or our COD value (at OD 635) $\times 0.90 =$ COD value (Merck Method).

Hence for a preliminary estimate of COD value, color of the standard reaction mixture based on a range of glucose concentrations can be used as a reference. For samples, with COD concentration of 300 to 500 mg/L, OD 585 or OD 635 can be recorded for obtaining more precise values.

Table 21: Correlation between COD concentration, color of the reaction mixture and their OD at different wave lengths.

COD (mg/L)	Sample dilution (Times)	Color of the reaction mixture	OD	Estimated COD (mg/L)	Estimated final COD (mg/L) (Col. 2x5)
1	2	3	4	5	6

Values at OD 585 nm

Glucose:

320	NA	C300	0.159	NA	NA
426	NA	C400	0.186	NA	NA
535	NA	C300	0.205	NA	NA

Sample: Damaged wheat grain (filtrate) + Vegetable waste

Nil	C2000	0.438	≈2000#	≈2000
5	C600	0.229	≈ 580#	≈2900
10	C450	0.162	444	4440

Values at OD 635 nm

Glucose:

320	NA	C300	0.145	NA	NA
426	NA	C400	0.166	NA	NA
535	NA	C550	0.187	NA	NA

Sample: Damaged wheat grain (filtrate) + Vegetable waste

Nil	C2000	0.378	≈2000#	≈2000
5	C600	0.206	≈ 580#	≈2900
10	C450	0.146	449	4490

#: COD as deducible from Table 2 given in Example 2. D: Dark. T: Turbid.

*: OD out of range. NA: Not applicable. Samples are diluted to get an OD 585 in the range of 0.130 to 0.175 and OD 635 in the range of 0.120 to 0.160, to obtain a COD value in the range of 320 to 535 mg/L. (Based on the Tables 3 to 15 given in Examples 3 to 15).

EXAMPLE 22:

Dissolved 0.15 g glucose in 25 mL distilled water. Aliquots from the stock solution were taken and diluted with water to achieve different COD concentrations in the range of 320 to 535 mg/L. A 30 mL sample of Damaged wheat grain (filtrate) mixed with palm oil mill effluent (filtrate) was taken. Weighed 0.02 g HgSO₄ in separate 15

mL test tubes. 1.0 mL of glucose solution and sample were added to each test tube separately and mixed well. 0.5 mL of $K_2Cr_2O_7$ was added to each sample. 1.5 mL H_2SO_4 silver sulphate reagent was added slowly with continuous shaking. The color of the reaction mixture was noted down visually within 15 sec of incubation and also recorded spectrophotometrically at 585 and 635 nm against air. An approximate COD value of the sample was estimated from the color of standard reaction mixture prepared with glucose as reference material in the COD range of 200 to 100000 mg/L. OD 585 and OD 635 of the reaction mixture was read spectrophotometrically. OD 585 and 635 nm in the range of 0.130 to 0.175 and 0.120 to 0.160 respectively were considered and calculated with reference to glucose standard. (Glucose COD in the range of 320 to 535 mg/L corresponds to OD 585 and 635 nm in the range of 0.157 to 0.207 and 0.142 to 0.189, respectively).

COD of the sample was estimated to be 5391 mg/L. (By Merck method and read on photometer).

Based on the color of the reaction mixture, COD of the sample was estimated to be around 3700 mg/L. On the basis of the dilution method and the spectrophotometric reading at OD 585 and OD 635 of the sample, COD was estimated to be 5344 mg/L. A comparison with the COD value calculated through Merck method reveals the following relation i.e. Our COD value (at OD 585 or OD 635) $\times 1.01$ = COD value (Merck Method).

Hence for a preliminary estimate of COD value, color of the standard reaction mixture based on a range of glucose concentrations can be used as a reference. For samples, with COD concentration of 300 to 500 mg/L, OD 585 or OD 635 can be recorded for obtaining more precise values.

Table 22: Correlation between COD concentration, color of the reaction mixture and their OD at different wave lengths.

COD (mg/L)	Sample dilution (Times)	Color of the reaction mixture	OD	Estimated COD (mg/L)	Estimated final COD (mg/L) (Col. 2x5)
1	2	3	4	5	6
<u>Values at OD 585 nm</u>					
Glucose:					
320	NA	C300	0.161	NA	NA
426	NA	C400	0.186	NA	NA
535	NA	C550	0.206	NA	NA
Sample: Damaged wheat grain (filtrate) + Palm oil mill effluent (filtrate)					
	Nil	C2000	0.533	≈2000#	≈2000
	5	C600	0.258	≈ 580#	≈2900
	10	C400	0.183	≈ 373#	≈3730
	16.7	C300	0.133	320	5344
<u>Values at OD 635 nm</u>					
Glucose:					
320	NA	C300	0.146	NA	NA
426	NA	C400	0.167	NA	NA
535	NA	C550	0.186	NA	NA
Sample: Damaged wheat grain (filtrate) + Palm oil mill effluent (filtrate)					
	Nil	C2000	0.442	≈2000#	≈2000
	5	C600	0.234	≈ 580#	≈2900
	10	C400	0.165	≈ 373#	≈3730
	16.7	C300	0.120	320	5344

#: COD as deducible from Table 2 given in Example 2. D: Dark. T: Turbid.

*: OD out of range. NA: Not applicable. Samples are diluted to get an OD 585 in the range of 0.130 to 0.175 and OD 635 in the range of 0.120 to 0.160, to obtain a COD value in the range of 320 to 535 mg/L. (Based on the Tables 3 to 15 given in Examples 3 to 15).

EXAMPLE 23:

Dissolved 0.15 g glucose in 25 mL distilled water. Aliquots from the stock solution were taken and diluted with water to achieve different COD concentrations in the range of 320 to 535 mg/L. A 30 mL sample of Damaged wheat grain (filtrate) was taken. Weighed 0.02 g HgSO_4 in separate 15 mL test tubes. 1.0 mL of glucose solution and sample were added to each test tube separately and mixed well. 0.5 ml of $\text{K}_2\text{Cr}_2\text{O}_7$ was added to each sample. 1.5 mL H_2SO_4 silver sulphate reagent was added slowly with continuous shaking. The color of the reaction mixture was noted down visually within 15 sec of incubation and also recorded spectrophotometrically at 585 and 635 nm against air. An approximate COD value of the sample was estimated from the color of standard reaction mixture prepared with glucose as reference material in the COD range of 200 to 100000 mg/L. OD 585 and OD 635 of the reaction mixture was read spectrophotometrically. OD 585 and 635 nm in the range of 0.130 to 0.175 and 0.120 to 0.160 respectively were considered and calculated with reference to glucose standard. (Glucose COD in the range of 320 to 535 mg/L corresponds to OD 585 and 635 nm in the range of 0.157 to 0.207 and 0.142 to 0.189, respectively).

COD of the sample was estimated to be 6250 mg/L. (By Merck method and read on photometer).

Based on the color of the reaction mixture, COD of the sample was estimated to be around 5800 mg/L. On the basis of the dilution method and the spectrophotometric reading at OD 585 and OD 635 of the sample, COD was estimated to be 8020 and 6700 mg/L, respectively. A comparison with the COD value calculated through Merck method reveals the following relation i.e. Our COD value (at OD 585) $\times 0.78 =$ COD value (Merck Method) or our COD value (at OD 635) $\times 0.93 =$ COD value (Merck Method).

Hence for a preliminary estimate of COD value, color of the standard reaction mixture based on a range of glucose concentrations can be used as a reference. For samples, with COD concentration of 300 to 500 mg/L, OD 585 or OD 635 can be recorded for obtaining more precise values.

Table 23: Correlation between COD concentration, color of the reaction mixture and their OD at different wave lengths.

COD (mg/L)	Sample dilution (Times)	Color of the reaction mixture	OD	Estimated COD (mg/L)	Estimated final COD (mg/L) (Col. 2x5)
1	2	3	4	5	6
<u>Values at OD 585 nm</u>					
Glucose:					
320	NA	C300	0.159	NA	NA
426	NA	C400	0.186	NA	NA
535	NA	C550	0.205	NA	NA
Sample: Damaged wheat grain (filtrate)					
	Nil	C1900	1.325	≅1900#	≅1900
	10	C600	0.204	≅ 580#	≅5800
	20	C400	0.146	401	8020
<u>Values at OD 635 nm</u>					
Glucose:					
320	NA	C300	0.149	NA	NA
426	NA	C400	0.166	NA	NA
535	NA	C550	0.183	NA	NA
Sample: Damaged wheat grain (filtrate)					
	Nil	C1900	1.146	≅1900#	≅1900
	10	C600	0.184	≅ 580#	≅5800
	20	C350	0.130	335	6700

#: COD as deducible from Table 2 given in Example 2. D: Dark. T: Turbid.

*: OD out of range. NA: Not applicable. Samples are diluted to get an OD 585 in the range of 0.130 to 0.175 and OD 635 in the range of 0.120 to 0.160, to obtain a COD value in the range of 320 to 535 mg/L. (Based on the Tables 3 to 15 given in Examples 3 to 15).

EXAMPLE 24:

Dissolved 0.15 g glucose in 25 mL distilled water. Aliquots from the stock solution were taken and diluted with water to achieve different COD concentrations in the range of 320 to 535 mg/L. A 30 mL sample of Damaged wheat grain (filtrate) slurry was taken. Weighed 0.02 g HgSO_4 in separate 15 mL test tubes. 1.0 mL of glucose

solution and sample were added to each test tube separately and mixed well. 0.5 mL of $K_2Cr_2O_7$ was added to each sample. 1.5 mL H_2SO_4 silver sulphate reagent was added slowly with continuous shaking. The color of the reaction mixture was noted down visually within 15 sec of incubation and also recorded spectrophotometrically at 585 and 635 nm against air. An approximate COD value of the sample was estimated from the color of standard reaction mixture prepared with glucose as reference material in the COD range of 200 to 100000 mg/L. OD 585 and OD 635 of the reaction mixture was read spectrophotometrically. OD 585 and 635 nm in the range of 0.130 to 0.175 and 0.120 to 0.160 respectively were considered and calculated with reference to glucose standard. (Glucose COD in the range of 320 to 535 mg/L corresponds to OD 585 and 635 nm in the range of 0.157 to 0.207 and 0.142 to 0.189, respectively).

COD of the sample was estimated to be 9794 mg/L. (By Merck method and read on photometer).

Based on the color of the reaction mixture, COD of the sample was estimated to be around 10,000 mg/L. On the basis of the dilution method and the spectrophotometric reading at OD 585 and OD 635 of the sample, COD was estimated to be 6840 and 6700 mg/L, respectively. A comparison with the COD value calculated through Merck method reveals the following relation i.e. Our COD value (at OD 585) $\times 1.43 =$ COD value (Merck Method) or our COD value (at OD 635) $\times 1.46 =$ COD value (Merck Method).

Hence for a preliminary estimate of COD value, color of the standard reaction mixture based on a range of glucose concentrations can be used as a reference. For samples, with COD concentration of 300 to 500 mg/L, OD 585 or OD 635 can be recorded for obtaining more precise values.

Table 24: Correlation between COD concentration, color of the reaction mixture and their OD at different wave lengths.

COD (mg/L)	Sample dilution (Times)	Color of the reaction mixture	OD	Estimated COD (mg/L)	Estimated final COD (mg/L) (Col. 2x5)
1	2	3	4	5	6
<u>Values at OD 585 nm</u>					
Glucose:					
320	NA	C300	0.159	NA	NA
426	NA	C400	0.186	NA	NA
535	NA	C550	0.205	NA	NA
Sample: Damaged wheat grain (filtrate)					
	Nil	C10000	0.834	$\cong 10000\#$	$\cong 10000$
	20	C350	0.142	342	6840
<u>Values at OD 635 nm</u>					
Glucose:					
320	NA	C300	0.145	NA	NA
426	NA	C400	0.166	NA	NA
535	NA	C550	0.187	NA	NA
Sample: Damaged wheat grain (filtrate)					
	Nil	C10000	0.654	$\cong 10000\#$	$\cong 10000$
	20	C350	0.127	335	6700

#: COD as deducible from Table 2 given in Example 2. D: Dark. T: Turbid.

*: OD out of range. NA: Not applicable. Samples are diluted to get an OD 585 in the range of 0.130 to 0.175 and OD 635 in the range of 0.120 to 0.160, to obtain a COD value in the range of 320 to 535 mg/L (Based on the Tables 3 to 15 given in Examples 3 to 15).

EXAMPLE 25:

Dissolved 0.15 g glucose in 25 mL distilled water. Aliquots from the stock solution were taken and diluted with water to achieve different COD concentrations in the range of 320 to 535 mg/L. A 30 mL sample of Damaged wheat grain (filtrate) slurry was taken. Weighed 0.02 g HgSO_4 in separate 15 mL test tubes. 1.0 mL of glucose

solution and sample were added to each test tube separately and mixed well. 0.5 mL of $K_2Cr_2O_7$ was added to each sample. 1.5 mL H_2SO_4 silver sulphate reagent was added slowly with continuous shaking. The color of the reaction mixture was noted down visually within 15 sec of incubation and also recorded spectrophotometrically at 585 and 635 nm against air. An approximate COD value of the sample was estimated from the color of standard reaction mixture prepared with glucose as reference material in the COD range of 200 to 100000 mg/L. OD 585 and OD 635 of the reaction mixture was read spectrophotometrically. OD 585 and 635 nm in the range of 0.130 to 0.175 and 0.120 to 0.160 respectively were considered and calculated with reference to glucose standard. (Glucose COD in the range of 320 to 535 mg/L corresponds to OD 585 and 635 nm in the range of 0.157 to 0.207 and 0.142 to 0.189, respectively).

COD of the sample was estimated to be 15502 mg/L. (By Merck method and read on photometer).

Based on the color of the reaction mixture, COD of the sample was estimated to be around 10,000 mg/L. On the basis of the dilution method and the spectrophotometric reading at OD 585 and OD 635 of the sample, COD was estimated to be 8000 and 7725 mg/L, respectively. A comparison with the COD value calculated through Merck method reveals the following relation i.e. Our COD value (at OD 585) $\times 1.94 =$ COD value (Merck Method) or our COD value (at OD 635) $\times 2.01 =$ COD value (Merck Method).

Hence for a preliminary estimate of COD value, color of the standard reaction mixture based on a range of glucose concentrations can be used as a reference. For samples, with COD concentration of 300 to 500 mg/L, OD 585 or OD 635 can be recorded for obtaining more precise values.

Table 25: Correlation between COD concentration, color of the reaction mixture and their OD at different wave lengths.

COD (mg/L)	Sample dilution (Times)	Color of the reaction mixture	OD	Estimated COD (mg/L)	Estimated final COD (mg/L) (Col. 2x5)
1	2	3	4	5	6

Values at OD 585 nm

Glucose:

320	NA	C300	0.164	NA	NA
426	NA	C400	0.186	NA	NA
535	NA	C550	0.205	NA	NA

Sample: Damaged wheat grain (filtrate)

Nil	C10000	0.964	>10000#	> 10000
10	C600	0.267	> 600#	> 6000
25	C300	0.137	320	8000

Values at OD 635 nm

Glucose:

320	NA	C300	0.152	NA	NA
426	NA	C400	0.166	NA	NA
535	NA	C550	0.184	NA	NA

Sample: Damaged wheat grain (filtrate)

Nil	C10000	0.714	>10000#	> 10000
10	C600	0.252	> 600#	> 6000
25	C300	0.122	309	7725

#: COD as deducible from Table 2 given in Example 2. D: Dark. T: Turbid.

*: OD out of range. NA: Not applicable. Samples are diluted to get an OD 585 in the range of 0.130 to 0.175 and OD 635 in the range of 0.120 to 0.160, to obtain a COD value in the range of 320 to 535 mg/L. (Based on the Tables 3 to 15 given in Examples 3 to 15).

EXAMPLE 26:

Dissolved 0.15 g glucose in 25 mL distilled water. Aliquots from the stock solution were taken and diluted with water to achieve different COD concentrations in the

range of 320 to 535 mg/L. A 30 mL sample of Cattle dung slurry was taken. Weighed 0.02 g HgSO_4 in separate 15 mL test tubes. 1.0 mL of glucose solution and sample were added to each test tube separately and mixed well. 0.5 mL of $\text{K}_2\text{Cr}_2\text{O}_7$ was added to each sample. 1.5 mL H_2SO_4 silver sulphate reagent was added slowly with continuous shaking. The color of the reaction mixture was noted down visually within 15 sec of incubation and also recorded spectrophotometrically at 585 and 635 nm against air. An approximate COD value of the sample was estimated from the color of standard reaction mixture prepared with glucose as reference material in the COD range of 200 to 100000 mg/l. OD 585 and OD 635 of the reaction mixture was read spectrophotometrically. OD 585 and 635 nm in the range of 0.130 to 0.175 and 0.120 to 0.160 respectively were considered and calculated with reference to glucose standard. (Glucose COD in the range of 320 to 535 mg/L corresponds to OD 585 and 635 nm in the range of 0.157 to 0.207 and 0.142 to 0.189, respectively). COD of the sample was estimated to be 18515 mg/L. (By Merck method and read on photometer).

Based on the color of the reaction mixture, COD of the sample was estimated to be around 10,000 mg/L. On the basis of the dilution method and the spectrophotometric reading at OD 585 and OD 635 of the sample, COD was estimated to be 8660 mg/L. A comparison with the COD value calculated through Merck method reveals the following relation i.e. Our COD value (at OD 585 or OD 635) $\times 2.14 =$ COD value (Merck Method).

Hence for a preliminary estimate of COD value, color of the standard reaction mixture based on a range of glucose concentrations can be used as a reference. For samples, with COD concentration of 300 to 500 mg/L, OD 585 or OD 635 can be recorded for obtaining more precise values.

Table 26: Correlation between COD concentration, color of the reaction mixture and their OD at different wave lengths.

COD (mg/L)	Sample dilution (Times)	Color of the reaction mixture	OD	Estimated COD (mg/L)	Estimated final COD (mg/L) (Col. 2x5)
1	2	3	4	5	6

Values at OD 585 nm

Glucose:

320	NA	C300	0.162	NA	NA
426	NA	C400	0.185	NA	NA
535	NA	C550	0.205	NA	NA

Sample: Cattle dung slurry

Nil	C10000	1.878	≅10000#	≅ 10000
20	C450	0.157	433	8660

Values at OD 635 nm

Glucose:

320	NA	C300	0.148	NA	NA
426	NA	C400	0.171	NA	NA
535	NA	C550	0.189	NA	NA

Sample: Cattle dung slurry

Nil	C10000	1.637	≅10000#	≅ 10000
20	C450	0.145	433	8660

#: COD as deducible from Table 2 given in Example 2. D: Dark. T: Turbid.

*: OD out of range. NA: Not applicable. Samples are diluted to get an OD 585 in the range of 0.130 to 0.175 and OD 635 in the range of 0.120 to 0.160, to obtain a COD value in the range of 320 to 535 mg/L. (Based on the Tables 3 to 15 given in Examples 3 to 15).

EXAMPLE 27:

Dissolved 0.15 g glucose in 25 mL distilled water. Aliquots from the stock solution were taken and diluted with water to achieve different COD concentrations in the range of 320 to 535 mg/L. A 30 mL sample of Palm Oil Mill Effluent (Digested) with Damaged wheat grains (filtrate) slurry was taken. Weighed 0.02 g HgSO₄ in

separate 15 mL test tubes. 1.0 mL of glucose solution and sample were added to each test tube separately and mixed well. 0.5 mL of $K_2Cr_2O_7$ was added to each sample. 1.5 mL H_2SO_4 silver sulphate reagent was added slowly with continuous shaking. The color of the reaction mixture was noted down visually within 15 sec of incubation and also recorded spectrophotometrically at 585 and 635 nm against air. An approximate COD value of the sample was estimated from the color of standard reaction mixture prepared with glucose as reference material in the COD range of 200 to 100000 mg/l. OD 585 and OD 635 of the reaction mixture was read spectrophotometrically. OD 585 and 635 nm in the range of 0.130 to 0.175 and 0.120 to 0.160 respectively were considered and calculated with reference to glucose standard. (Glucose COD in the range of 320 to 535 mg/L corresponds to OD 585 and 635 nm in the range of 0.157 to 0.207 and 0.142 to 0.189, respectively). COD of the sample was estimated to be 19941 mg/L. (By Merck method and read on photometer).

Based on the color of the reaction mixture, COD of the sample was estimated to be around 12,000 mg/L. On the basis of the dilution method and the spectrophotometric reading at OD 585 and OD 635 of the sample, COD was estimated to be 10656 and 11355 mg/L, respectively. A comparison with the COD value calculated through Merck method reveals the following relation i.e. Our COD value (at OD 585) $\times 1.87 =$ COD value (Merck Method) or our COD value (at OD 635) $\times 1.76 =$ COD value (Merck Method).

Hence for a preliminary estimate of COD value, color of the standard reaction mixture based on a range of glucose concentrations can be used as a reference. For samples, with COD concentration of 300 to 500 mg/L, OD 585 or OD 635 can be recorded for obtaining more precise values.

Table 27: Correlation between COD concentration, color of the reaction mixture and their OD at different wave lengths.

COD (mg/L)	Sample dilution (Times)	Color of the reaction mixture	OD	Estimated COD (mg/L)	Estimated final COD (mg/L) (Col. 2x5)
1	2	3	4	5	6
<u>Values at OD 585 nm</u>					
Glucose:					
320	NA	C300	0.158	NA	NA
426	NA	C400	0.185	NA	NA
535	NA	C550	0.205	NA	NA
Sample: Palm Oil Mill Effluent (Digested) + Damaged wheat grains (filtrate)					
	Nil	C10000	6.000*	>10000#	> 10000
	20	C600	0.195	≅ 600#	≅ 12000
	33.3	C300	0.131	320	10656
<u>Values at OD 635 nm</u>					
Glucose:					
320	NA	C300	0.137	NA	NA
426	NA	C400	0.163	NA	NA
535	NA	C550	0.185	NA	NA
Sample: Palm Oil Mill Effluent (Digested) + Damaged wheat grains (filtrate)					
	Nil	C10000	6.000*	>10000#	> 10000
	20	C600	0.182	≅ 600#	≅ 12000
	33.3	C300	0.122	341	11355

#: COD as deducible from Table 2 given in Example 2. D: Dark. T: Turbid.

*: OD out of range. NA: Not applicable. Samples are diluted to get an OD 585 in the range of 0.130 to 0.175 and OD 635 in the range of 0.120 to 0.160, to obtain a COD value in the range of 320 to 535 mg/L (Based on the Tables 3 to 15 given in Examples 3 to 15).

EXAMPLE 28:

Dissolved 0.15 g glucose in 25 mL distilled water. Aliquots from the stock solution were taken and diluted with water to achieve different COD concentrations in the range of 320 to 535 mg/L. A 30 mL sample of Palm Oil Mill Effluent (Digested) was taken. Weighed 0.02 g HgSO_4 in separate 15 mL test tubes. 1.0 mL of glucose

solution and sample were added to each test tube separately and mixed well. 0.5 mL of $K_2Cr_2O_7$ was added to each sample. 1.5 mL H_2SO_4 silver sulphate reagent was added slowly with continuous shaking. The color of the reaction mixture was noted down visually within 15 sec of incubation and also recorded spectrophotometrically at 585 and 635 nm against air. An approximate COD value of the sample was estimated from the color of standard reaction mixture prepared with glucose as reference material in the COD range of 200 to 100000 mg/l. OD 585 and OD 635 of the reaction mixture was read spectrophotometrically. OD 585 and 635 nm in the range of 0.130 to 0.175 and 0.120 to 0.160 respectively were considered and calculated with reference to glucose standard. (Glucose COD in the range of 320 to 535 mg/l corresponds to OD 585 and 635 nm in the range of 0.157 to 0.207 and 0.142 to 0.189, respectively).

COD of the sample was estimated to be 22287 mg/L. (By Merck method and read on photometer).

Based on the color of the reaction mixture, COD of the sample was estimated to be around 20,000 mg/L. On the basis of the dilution method and the spectrophotometric reading at OD 585 and OD 635 of the sample, COD was estimated to be 15840 and 13560 mg/l, respectively. A comparison with the COD value calculated through Merck method reveals the following relation i.e. Our COD value (at OD 585) $\times 1.41 =$ COD value (Merck Method) or our COD value (at OD 635) $\times 1.64 =$ COD value (Merck Method).

Hence for a preliminary estimate of COD value, color of standard reaction mixture can be used as a reference. For samples, with COD concentration of 300 to 500 mg/L, OD 585 or OD 635 can be recorded for obtaining more precise values.

Table 28: Correlation between COD concentration, color of the reaction mixture and their OD at different wave lengths.

COD (mg/L)	Sample dilution (Times)	Color of the reaction mixture	OD	Estimated COD (mg/L)	Estimated final COD (mg/L) (Col. 2x5)
1	2	3	4	5	6
<u>Values at OD 585 nm</u>					
Glucose:					
320	NA	C300	0.162	NA	NA
426	NA	C400	0.185	NA	NA
535	NA	C550	0.205	NA	NA
Sample: Palm Oil Mill Effluent (Digested)					
	Nil	C10000	6.000*	>10000#	> 10000
	10	C2000	0.345	> 2000#	> 20000
	50	C250	0.119	≅ 250	≅ 12500
	40	C400	0.143	396	15840
<u>Values at OD 635 nm</u>					
Glucose:					
320	NA	C300	0.148	NA	NA
426	NA	C400	0.171	NA	NA
535	NA	C550	0.189	NA	NA
Sample: Palm Oil Mill Effluent (Digested)					
	Nil	C10000	6.000*	>10000#	> 10000
	10	C2000	0.330	> 2000#	> 20000
	50	C250	0.109	≅ 250	≅ 12500
	40	C400	0.131	339	13560

#: COD as deducible from Table 2 given in Example 2. D: Dark. T: Turbid.

*: OD out of range. NA: Not applicable. Samples are diluted to get an OD 585 in the range of 0.130 to 0.175 and OD 635 in the range of 0.120 to 0.160, to obtain a COD value in the range of 320 to 535 mg/L. (Based on the Tables 3 to 15 given in Examples 3 to 15).

EXAMPLE 29:

Dissolved 0.15 g glucose in 25 mL distilled water. Aliquots from the stock solution were taken and diluted with water to achieve different COD concentrations in the

range of 320 to 535 mg/L. A 30 mL sample of Palm Oil Mill Effluent (Digested) with Apple pomace was taken. Weighed 0.02 g HgSO_4 in separate 15 mL test tubes. 1.0 mL of glucose solution and sample were added to each test tube separately and mixed well. 0.5 mL of $\text{K}_2\text{Cr}_2\text{O}_7$ was added to each sample. 1.5 mL H_2SO_4 silver sulphate reagent was added slowly with continuous shaking. The color of the reaction mixture was noted down visually within 15 sec of incubation and also recorded spectrophotometrically at 585 and 635 nm against air. An approximate COD value of the sample was estimated from the color of standard reaction mixture prepared with glucose as reference material in the COD range of 200 to 100000 mg/L. OD 585 and OD 635 of the reaction mixture was read spectrophotometrically. OD 585 and 635 nm in the range of 0.130 to 0.175 and 0.120 to 0.160 respectively were considered and calculated with reference to glucose standard. (Glucose COD in the range of 320 to 535 mg/L corresponds to OD 585 and 635 nm in the range of 0.157 to 0.207 and 0.142 to 0.189, respectively).

COD of the sample was estimated to be 44620 mg/L. (By Merck method and read on photometer).

Based on the color of the reaction mixture, COD of the sample was estimated to be around 40,000 mg/L. On the basis of the dilution method and the spectrophotometric reading at OD 585 and OD 635 of the sample, COD was estimated to be 38640 and 40080 mg/L, respectively. A comparison with the COD value calculated through Merck method reveals the following relation i.e. Our COD value (at OD 585) $\times 1.15 =$ COD value (Merck Method) or our COD value (at OD 635) $\times 1.11 =$ COD value (Merck Method).

Hence for a preliminary estimate of COD value, color of the standard reaction mixture based on a range of glucose concentrations can be used as a reference. For samples, with COD concentration of 300 to 500 mg/L, OD 585 or OD 635 can be recorded for obtaining more precise values.

Table 29: Correlation between COD concentration, color of the reaction mixture and their OD at different wave lengths.

COD (mg/L)	Sample dilution (Times)	Color of the reaction mixture	OD	Estimated COD (mg/L)	Estimated final COD (mg/L) (Col. 2x5)
1	2	3	4	5	6

Values at OD 585 nm

Glucose:

320	NA	C300	0.155	NA	NA
426	NA	C400	0.185	NA	NA
535	NA	C550	0.205	NA	NA

Sample: Palm Oil Mill Effluent (Digested) + Apple pomace

Nil	C10000	6.000*	>10000#	> 10000
20	C2000	0.352	> 2000#	> 40000
100	C2000	0.220	≅ 250#	≅ 25000
120	C300	0.132	320	38640

Values at OD 635 nm

Glucose:

320	NA	C300	0.137	NA	NA
426	NA	C400	0.163	NA	NA
535	NA	C550	0.185	NA	NA

Sample: Palm Oil Mill Effluent (Digested) + Apple pomace

Nil	C10000	6.000*	>10000#	> 10000
20	C2000	0.322	> 2000#	> 40000
100	C2000	0.190	≅ 250#	≅ 25000
120	C300	0.119	334	40080

#: COD as deducible from Table 2 given in Example 2. D: Dark. T: Turbid.

*: OD out of range. NA: Not applicable. Samples are diluted to get an OD 585 in the range of 0.130 to 0.175 and OD 635 in the range of 0.120 to 0.160, to obtain a COD value in the range of 320 to 535 mg/L. (Based on the Tables 3 to 15 given in Examples 3 to 15).

EXAMPLE 30:

Dissolved 0.15 g glucose in 25 mL distilled water. Aliquots from the stock solution were taken and diluted with water to achieve different COD concentrations in the range of 320 to 535 mg/L. A 30 mL sample of Vegetable waste slurry with Apple

pomace was taken. Weighed 0.02 g HgSO_4 in separate 15 mL test tubes. 1.0 mL of glucose solution and sample were added to each test tube separately and mixed well. 0.5 mL of $\text{K}_2\text{Cr}_2\text{O}_7$ was added to each sample. 1.5 mL H_2SO_4 silver sulphate reagent was added slowly with continuous shaking. The color of the reaction mixture was noted down visually within 15 sec of incubation and also recorded spectrophotometrically at 585 and 635 nm against air. An approximate COD value of the sample was estimated from the color of standard reaction mixture prepared with glucose as reference material in the COD range of 200 to 100000 mg/l. OD 585 and OD 635 of the reaction mixture was read spectrophotometrically. OD 585 and 635 nm in the range of 0.130 to 0.175 and 0.120 to 0.160 respectively were considered and calculated with reference to glucose standard. (Glucose COD in the range of 320 to 535 mg/L corresponds to OD 585 and 635 nm in the range of 0.157 to 0.207 and 0.142 to 0.189, respectively).

COD of the sample was estimated to be 35484 mg/L. (By Merck method and read on photometer).

Based on the color of the reaction mixture, COD of the sample was estimated to be around 40,000 mg/L. On the basis of the dilution method and the spectrophotometric reading at OD 585 and OD 635 of the sample, COD was estimated to be 40000 and 34000 mg/L, respectively. A comparison with the COD value calculated through Merck method reveals the following relation i.e. Our COD value (at OD 585) $\times 0.89 =$ COD value (Merck Method) or our COD value (at OD 635) $\times 1.04 =$ COD value (Merck Method).

Hence for a preliminary estimate of COD value, color of the standard reaction mixture based on a range of glucose concentrations can be used as a reference. For samples, with COD concentration of 300 to 500 mg/L, OD 585 or OD 635 can be recorded for obtaining more precise values.

Table 30: Correlation between COD concentration, color of the reaction mixture and their OD at different wave lengths.

COD (mg/L)	Sample dilution (Times)	Color of the reaction mixture	OD	Estimated COD (mg/L)	Estimated final COD (mg/L) (Col. 2x5)
1	2	3	4	5	6

Values at OD 585 nm

Glucose:

320	NA	C300	0.163	NA	NA
426	NA	C400	0.185	NA	NA
535	NA	C550	0.202	NA	NA

Sample: Vegetable waste slurry+ Apple pomace

Nil	C10000	6.00*	>10000#	>10000
20	C2000	0.349	> 2000#	≅ 40000
100	C400	0.145	400	40000

Values at OD 635 nm

Glucose:

320	NA	C300	0.147	NA	NA
426	NA	C400	0.170	NA	NA
535	NA	C550	0.187	NA	NA

Sample: Vegetable waste slurry+ Apple pomace

Nil	C10000	6.00*	>10000#	>10000
20	C2000	0.318	> 2000#	≅ 40000
100	C400	0.130	340	34000

#: COD as deducible from Table 2 given in Example 2. D: Dark. T: Turbid.

*: OD out of range. NA: Not applicable. Samples are diluted to get an OD 585 in the range of 0.130 to 0.175 and OD 635 in the range of 0.120 to 0.160, to obtain a COD value in the range of 320 to 535 mg/L. (Based on the Tables 3 to 15 given in Examples 3 to 15).

EXAMPLE 31:

Dissolved 0.15 g glucose in 25 mL distilled water. Aliquots from the stock solution were taken and diluted with water to achieve different COD concentrations in the range of 320 to 535 mg/L. A 30 mL sample of Damaged wheat grain (filtrate) with Apple pomace was taken. Weighed 0.02 g HgSO₄ in separate 15 mL test tubes. 1.0

mL of glucose solution and sample were added to each test tube separately and mixed well. 0.5 mL of $K_2Cr_2O_7$ was added to each sample. 1.5 mL H_2SO_4 silver sulphate reagent was added slowly with continuous shaking. The color of the reaction mixture was noted down visually within 15 sec of incubation and also recorded spectrophotometrically at 585 and 635 nm against air. An approximate COD value of the sample was estimated from the color of standard reaction mixture prepared with glucose as reference material in the COD range of 200 to 100000 mg/L. OD 585 and OD 635 of the reaction mixture was read spectrophotometrically. OD 585 and 635 nm in the range of 0.130 to 0.175 and 0.120 to 0.160 respectively were considered and calculated with reference to glucose standard. (Glucose COD in the range of 320 to 535 mg/L corresponds to OD 585 and 635 nm in the range of 0.157 to 0.207 and 0.142 to 0.189, respectively).

COD of the sample was estimated to be 36018 mg/L. (By Merck method and read on photometer).

Based on the color of the reaction mixture, COD of the sample was estimated to be around 30,000 mg/L. On the basis of the dilution method and the spectrophotometric reading at OD 585 and OD 635 of the sample, COD was estimated to be 42600 mg/l. A comparison with the COD value calculated through Merck method reveals the following relation i.e. Our COD value (at OD 585 or OD 635) $\times 0.84$ = COD value (Merck Method).

Hence for a preliminary estimate of COD value color of the standard reaction mixture based on a range of glucose concentrations can be used as a reference. For samples, with COD concentration of 300 to 500 mg/L, OD 585 or OD 635 can be recorded for obtaining more precise values.

Table 31: Correlation between COD concentration, color of the reaction mixture and their OD at different wave lengths.

COD (mg/L)	Sample dilution (Times)	Color of the reaction mixture	OD	Estimated COD (mg/L)	Estimated final COD (mg/L) (Col. 2x5)
1	2	3	4	5	6

Values at OD 585 nm

Glucose:

320	NA	C300	0.164	NA	NA
426	NA	C400	0.186	NA	NA
535	NA	C550	0.205	NA	NA

Sample: Damaged wheat grain (filtrate)+ Apple pomace

Nil	C10000	6.000*	>10000#	> 10000
10	C1500	0.442	> 1500#	> 15000
50	C600	0.236	> 600#	> 30000
100	C400	0.155	426	42600

Values at OD 635 nm

Glucose:

320	NA	C300	0.152	NA	NA
426	NA	C400	0.166	NA	NA
535	NA	C550	0.184	NA	NA

Sample: Damaged wheat grain (filtrate)+ Apple pomace

Nil	C10000	6.000*	>10000#	> 10000
10	C1500	0.377	> 1500#	> 15000
50	C600	0.212	> 600#	> 30000
100	C400	0.138	426	42600

#: COD as deducible from Table 2 given in Example 2. D: Dark. T: Turbid.

*: OD out of range. NA: Not applicable. Samples are diluted to get an OD 585 in the range of 0.130 to 0.175 and OD 635 in the range of 0.120 to 0.160, to obtain a COD value in the range of 320 to 535 mg/L. (Based on the Tables 3 to 15 given in Examples 3 to 15).

EXAMPLE 32:

Dissolved 0.15 g glucose in 25 mL distilled water. Aliquots from the stock solution were taken and diluted with water to achieve different COD concentrations in the

range of 320 to 535 mg/L. A 30 mL sample of Damaged wheat grain slurry with Tea with Sodium sulphate was taken. Weighed 0.02 g HgSO_4 in separate 15 mL test tubes. 1.0 mL of glucose solution and sample were added to each test tube separately and mixed well. 0.5 mL of $\text{K}_2\text{Cr}_2\text{O}_7$ was added to each sample. 1.5 mL H_2SO_4 silver sulphate reagent was added slowly with continuous shaking. The color of the reaction mixture was noted down visually within 15 sec of incubation and also recorded spectrophotometrically at 585 and 635 nm against air. An approximate COD value of the sample was estimated from the color of standard reaction mixture prepared with glucose as reference material in the COD range of 200 to 100000 mg/L. OD 585 and OD 635 of the reaction mixture was read spectrophotometrically. OD 585 and 635 nm in the range of 0.130 to 0.175 and 0.120 to 0.160 respectively were considered and calculated with reference to glucose standard. (Glucose COD in the range of 320 to 535 mg/L corresponds to OD 585 and 635 nm in the range of 0.157 to 0.207 and 0.142 to 0.189, respectively).

COD of the sample was estimated to be 53867 mg/L. (By Merck method and read on photometer).

Based on the color of the reaction mixture, COD of the sample was estimated to be around 100,000 mg/L. On the basis of the dilution method and the spectrophotometric reading at OD 585 and OD 635 of the sample, COD was estimated to be 49896 and 51744 mg/L, respectively. A comparison with the COD value calculated through Merck method reveals the following relation i.e. Our COD value (at OD 585) $\times 1.08 = \text{COD value (Merck Method)}$ or our COD value (at OD 635) $\times 1.04 = \text{COD value (Merck Method)}$.

Hence for a preliminary estimate of COD value, color of the standard reaction mixture based on a range of glucose concentrations can be used as a reference. For samples, with COD concentration of 300 to 500 mg/L, OD 585 or OD 635 can be recorded for obtaining more precise values.

Table 32: Correlation between COD concentration, color of the reaction mixture and their OD at different wave lengths.

COD (mg/L)	Sample dilution (Times)	Color of the reaction mixture	OD	Estimated COD (mg/L)	Estimated final COD (mg/L) (Col. 2x5)
1	2	3	4	5	6

Values at OD 585 nm

Glucose:

320	NA	C300	0.163	NA	NA
426	NA	C400	0.183	NA	NA
535	NA	C550	0.205	NA	NA

Sample: Damaged wheat grain + Tea + Sodium sulphate

Nil	C10000	6.00*	>10000#	> 10000
10	C10000	0.801*	≅10000#	≅100000
154	C300	0.138	324	49896

Values at OD 635 nm

Glucose:

320	NA	C300	0.150	NA	NA
426	NA	C400	0.167	NA	NA
535	NA	C550	0.187	NA	NA

Sample: Damaged wheat grain + Tea + Sodium sulphate

Nil	C10000	6.00*	>10000#	> 10000
10	C10000	0.801*	≅10000#	≅100000
154	C300	0.124	336	51744

#: COD as deducible from Table 2 given in Example 2. D: Dark. T: Turbid.

*: OD out of range. NA: Not applicable. Samples are diluted to get an OD 585 in the range of 0.130 to 0.175 and OD 635 in the range of 0.120 to 0.160, to obtain a COD value in the range of 320 to 535 mg/L. (Based on the Tables 3 to 15 given in Examples 3 to 15).

EXAMPLE 33:

Dissolved 0.15 g glucose in 25 mL distilled water. Aliquots from the stock solution were taken and diluted with water to achieve different COD concentrations in the

range of 320 to 535 mg/L. A 30 mL sample of Apple pomace was taken. Weighed 0.02 g HgSO_4 in separate 15 mL test tubes. 1.0 mL of glucose solution and sample were added to each test tube separately and mixed well. 0.5 mL of $\text{K}_2\text{Cr}_2\text{O}_7$ was added to each sample. 1.5 mL H_2SO_4 silver sulphate reagent was added slowly with continuous shaking. The color of the reaction mixture was noted down visually within 15 sec of incubation and also recorded spectrophotometrically at 585 and 635 nm against air. An approximate COD value of the sample was estimated from the color of standard reaction mixture prepared with glucose as reference material in the COD range of 200 to 100000 mg/L. OD 585 and OD 635 of the reaction mixture was read spectrophotometrically. OD 585 and 635 nm in the range of 0.130 to 0.175 and 0.120 to 0.160 respectively were considered and calculated with reference to glucose standard. (Glucose COD in the range of 320 to 535 mg/L corresponds to OD 585 and 635 nm in the range of 0.157 to 0.207 and 0.142 to 0.189, respectively). COD of the sample was estimated to be 55890 mg/L. (By Merck method and read on photometer).

Based on the color of the reaction mixture, COD of the sample was estimated to be around 50,000 mg/L. On the basis of the dilution method and the spectrophotometric reading at OD 585 and OD 635 of the sample, COD was estimated to be 65875 and 66500 mg/L, respectively. A comparison with the COD value calculated through Merck method reveals the following relation i.e. Our COD value (at OD 585) $\times 0.85 = \text{COD value (Merck Method)}$ or our COD value (at OD 635) $\times 0.84 = \text{COD value (Merck Method)}$.

Hence for a preliminary estimate of COD value, color of the standard reaction mixture based on a range of glucose concentrations can be used as a reference. For samples, with COD concentration of 300 to 500 mg/L, OD 585 or OD 635 can be recorded for obtaining more precise values.

Table 33: Correlation between COD concentration, color of the reaction mixture and their OD at different wave lengths.

COD (mg/L)	Sample dilution (Times)	Color of the reaction mixture	OD	Estimated COD (mg/L)	Estimated final COD (mg/L) (Col. 2x5)
1	2	3	4	5	6
<u>Values at OD 585 nm</u>					
Glucose:					
320	NA	C300	0.161	NA	NA
426	NA	C400	0.186	NA	NA
535	NA	C550	0.206	NA	NA
Sample: Apple pomace					
	Nil	C10000	6.00*	>10000#	>10000
	10	C5000	0.525	≅ 5000#	≅ 50000
	125	C500	0.171	527	65875
<u>Values at OD 635 nm</u>					
Glucose:					
320	NA	C300	0.150	NA	NA
426	NA	C400	0.167	NA	NA
535	NA	C550	0.188	NA	NA
Sample: Apple pomace					
	Nil	C10000	6.00*	>10000#	>10000
	10	C5000	0.428	≅ 5000#	≅ 50000
	125	C500	0.157	532	66500

#: COD as deducible from Table 2 given in Example 2. D: Dark. T: Turbid.

*: OD out of range. NA: Not applicable. Samples are diluted to get an OD 585 in the range of 0.130 to 0.175 and OD 635 in the range of 0.120 to 0.160, to obtain a COD value in the range of 320 to 535 mg/L (Based on the Tables 3 to 15 given in Examples 3 to 15).

EXAMPLE 34:

Dissolved 0.15 g glucose in 25 mL distilled water. Aliquots from the stock solution were taken and diluted with water to achieve different COD concentrations in the range of 320 to 535 mg/L. A 30 mL sample of Damaged wheat grain with Tea was

taken. Weighed 0.02 g HgSO_4 in separate 15 mL test tubes. 1.0 mL of glucose solution and sample were added to each test tube separately and mixed well. 0.5 mL of $\text{K}_2\text{Cr}_2\text{O}_7$ was added to each sample. 1.5 mL H_2SO_4 silver sulphate reagent was added slowly with continuous shaking. The color of the reaction mixture was noted down visually within 15 sec of incubation and also recorded spectrophotometrically at 585 and 635 nm against air. An approximate COD value of the sample was estimated from the color of standard reaction mixture prepared with glucose as reference material in the COD range of 200 to 100000 mg/l. OD 585 and OD 635 of the reaction mixture was read spectrophotometrically. OD 585 and 635 nm in the range of 0.130 to 0.175 and 0.120 to 0.160 respectively were considered and calculated with reference to glucose standard. (Glucose COD in the range of 320 to 535 mg/L corresponds to OD 585 and 635 nm in the range of 0.157 to 0.207 and 0.142 to 0.189, respectively).

COD of the sample was estimated to be 58760 mg/L. (By Merck method and read on photometer).

Based on the color of the reaction mixture, COD of the sample was estimated to be around 64800 mg/L. On the basis of the dilution method and the spectrophotometric reading at OD 585 and OD 635 of the sample, COD was estimated to be 78848 and 84546 mg/L, respectively. A comparison with the COD value calculated through Merck method reveals the following relation i.e. Our COD value (at OD 585) $\times 0.74 =$ COD value (Merck Method) or our COD value (at OD 635) $\times 0.70 =$ COD value (Merck Method).

Hence for a preliminary estimate of COD value, color of the standard reaction mixture based on a range of glucose concentrations can be used as a reference. For samples, with COD concentration of 300 to 500 mg/L, OD 585 or OD 635 can be recorded for obtaining more precise values.

Table 34: Correlation between COD concentration, color of the reaction mixture and their OD at different wave lengths.

COD (mg/L)	Sample dilution (Times)	Color of the reaction mixture	OD	Estimated COD (mg/L)	Estimated final COD (mg/L) (Col. 2x5)
1	2	3	4	5	6

Values at OD 585 nm

Glucose:

320	NA	C300	0.162	NA	NA
426	NA	C400	0.186	NA	NA
535	NA	C550	0.206	NA	NA

Sample: Damaged wheat grain + Tea

Nil	C10000	6.00*	>10000#	> 10000
10	C10000	6.00*	>10000#	≅100000
100	C1500	0.272	≅1500#	≅150000
270	C250	0.124	≅ 240#	≅ 64800
154	C500	0.164	512	78848

Values at OD 635 nm

Glucose:

320	NA	C300	0.150	NA	NA
426	NA	C400	0.167	NA	NA
535	NA	C500	0.187	NA	NA

Sample: Damaged wheat grain + Tea

Nil	C10000	6.00*	>10000#	> 10000
10	C10000	6.00*	>10000#	≅100000
100	C1500	0.272	≅1500#	≅150000
270	C250	0.124	≅ 240#	≅ 64800
154	C500	0.160	549	84546

#: COD as deducible from Table 2 given in Example 2. D: Dark. T: Turbid.

*: OD out of range. NA: Not applicable. Samples are diluted to get an OD 585 in the range of 0.130 to 0.175 and OD 635 in the range of 0.120 to 0.160, to obtain a COD value in the range of 320 to 535 mg/L. (Based on the Tables 3 to 15 given in Examples 3 to 15).

EXAMPLE 35:

Dissolved 0.15 g glucose in 25 mL distilled water. Aliquots from the stock solution were taken and diluted with water to achieve different COD concentrations in the range of 320 to 535 mg/L. A 30 mL sample of Trizyme with Tea was taken. Weighed 0.02 g HgSO_4 in separate 15 mL test tubes. 1.0 mL of glucose solution and sample were added to each test tube separately and mixed well. 0.5 mL of $\text{K}_2\text{Cr}_2\text{O}_7$ was added to each sample. 1.5 mL H_2SO_4 silver sulphate reagent was added slowly with continuous shaking. The color of the reaction mixture was noted down visually within 15 sec of incubation and also recorded spectrophotometrically at 585 and 635 nm against air. An approximate COD value of the sample was estimated from the color of standard reaction mixture prepared with glucose as reference material in the COD range of 200 to 100000 mg/L. OD 585 and OD 635 of the reaction mixture was read spectrophotometrically. OD 585 and 635 nm in the range of 0.130 to 0.175 and 0.120 to 0.160 respectively were considered and calculated with reference to glucose standard. (Glucose COD in the range of 320 to 535 mg/L corresponds to OD 585 and 635 nm in the range of 0.157 to 0.207 and 0.142 to 0.189, respectively). COD of the sample was estimated to be 61088 mg/L. (By Merck method and read on photometer).

Based on the color of the reaction mixture, COD of the sample was estimated to be around 100,000 mg/L. On the basis of the dilution method and the spectrophotometric reading at OD 585 and OD 635 of the sample, COD was estimated to be 79156 and 67914 mg/L, respectively. A comparison with the COD value calculated through Merck method reveals the following relation i.e. Our COD value (at OD 585) $\times 0.77 = \text{COD value (Merck Method)}$ or our COD value (at OD 635) $\times 0.90 = \text{COD value (Merck Method)}$.

Hence for a preliminary estimate of COD value color of the standard reaction mixture based on a range of glucose concentrations can be used as a reference. For samples, with COD concentration of 300 to 500 mg/L, OD 585 or OD 635 can be recorded for obtaining more precise values.

Table 35: Correlation between COD concentration, color of the reaction mixture and their OD at different wave lengths.

COD (mg/L)	Sample dilution (Times)	Color of the reaction mixture	OD	Estimated COD (mg/L)	Estimated final COD (mg/L) (Col. 2x5)
1	2	3	4	5	6

Values at OD 585 nm

Glucose:

320	NA	C300	0.163	NA	NA
426	NA	C400	0.183	NA	NA
535	NA	C550	0.205	NA	NA

Sample: Trizyme + Tea

Nil	C10000	6.00*	>10000#	> 10000
50	C2000	0.338	≅ 2000#	≅100000
154	C500	0.164	514	79156

Values at OD 635 nm

Glucose:

320	NA	C300	0.152	NA	NA
426	NA	C400	0.170	NA	NA
535	NA	C550	0.184	NA	NA

Sample: Trizyme + Tea

Nil	C10000	6.00*	>10000#	> 10000
50	C2000	0.305	≅ 2000#	≅100000
154	C500	0.147	441	67914

#: COD as deducible from Table 2 given in Example 2. D: Dark. T: Turbid.

*: OD out of range. NA: Not applicable. Samples are diluted to get an OD 585 in the range of 0.130 to 0.175 and OD 635 in the range of 0.120 to 0.160, to obtain a COD value in the range of 320 to 535 mg/L. (Based on the Tables 3 to 15 given in Examples 3 to 15).

EXAMPLE 36:

Dissolved 0.15 g glucose in 25 mL distilled water. Aliquots from the stock solution were taken and diluted with water to achieve different COD concentrations in the range of 320 to 535 mg/L. A 30 mL sample of Tea was taken. Weighed 0.02 g

HgSO₄ in separate 15 mL test tubes. 1.0 mL of glucose solution and sample were added to each test tube separately and mixed well. 0.5 mL of K₂Cr₂O₇ was added to each sample. 1.5 mL H₂SO₄ silver sulphate reagent was added slowly with continuous shaking. The color of the reaction mixture was noted down visually within 15 sec of incubation and also recorded spectrophotometrically at 585 and 635 nm against air. An approximate COD value of the sample was estimated from the color of standard reaction mixture prepared with glucose as reference material in the COD range of 200 to 100000 mg/L. OD 585 and OD 635 of the reaction mixture was read spectrophotometrically. OD 585 and 635 nm in the range of 0.130 to 0.175 and 0.120 to 0.160 respectively were considered and calculated with reference to glucose standard. (Glucose COD in the range of 320 to 535 mg/L corresponds to OD 585 and 635 nm in the range of 0.157 to 0.207 and 0.142 to 0.189, respectively). COD of the sample was estimated to be 1,36,608 mg/L. (By Merck method and read on photometer).

Based on the color of the reaction mixture, COD of the sample was estimated to be around 200,000 to 220,000 mg/L. On the basis of the dilution method and the spectrophotometric reading at OD 585 and OD 635 of the sample, COD was estimated to be 173983 and 172649 mg/L, respectively. A comparison with the COD value calculated through Merck method reveals the following relation i.e. Our COD value (at OD 585) x 0.78 = COD value (Merck Method) or our COD value (at OD 635) x 0.79 = COD value (Merck Method).

Hence for a preliminary estimate of COD value, color of the standard reaction mixture based on a range of glucose concentrations can be used as a reference. For samples, with COD concentration of 300 to 500 mg/L, OD 585 or OD 635 can be recorded for obtaining more precise values.

Table 36: Correlation between COD concentration, color of the reaction mixture and their OD at different wave lengths.

COD (mg/L)	Sample dilution (Times)	Color of the reaction mixture	OD	Estimated COD (mg/L)	Estimated final COD (mg/L) (Col. 2x5)
1	2	3	4	5	6

Values at OD 585 nm

Glucose:

320	NA	C300	0.161	NA	NA
426	NA	C400	0.187	NA	NA
535	NA	C550	0.202	NA	NA

Sample: Tea

Nil	C10000	6.00*	>10000#	>10000
20	C10000	1.648	>10000#	≅ 200000
200	C1100	0.216	≅ 1100#	≅ 220000
333.3	C500	0.164	522	173983

Values at OD 635 nm

Glucose:

320	NA	C300	0.144	NA	NA
426	NA	C400	0.167	NA	NA
535	NA	C550	0.184	NA	NA

Sample: Tea

Nil	C10000	6.00*	>10000#	>10000
20	C10000	1.531	>10000#	≅ 200000
200	C1100	0.196	≅ 1100#	≅ 220000
333.3	C500	0.148	518	172649

#: COD as deducible from Table 2 given in Example 2. D: Dark. T: Turbid.

*: OD out of range. NA: Not applicable. Samples are diluted to get an OD 585 in the range of 0.130 to 0.175 and OD 635 in the range of 0.120 to 0.160, to obtain a COD value in the range of 320 to 535 mg/L. (Based on the Tables 3 to 15 given in Examples 3 to 15).

EXAMPLE 37:

Dissolved 0.15 g glucose in 25 mL distilled water. Aliquots from the stock solution were taken and diluted with water to achieve different COD concentrations in the range of 320 to 535 mg/L. A 30 mL sample of Damaged wheat grain (as such) was taken. Weighed 0.02 g HgSO_4 in separate 15 mL test tubes. 1.0 mL of glucose solution and sample were added to each test tube separately and mixed well. 0.5 mL of $\text{K}_2\text{Cr}_2\text{O}_7$ was added to each sample. 1.5 mL H_2SO_4 silver sulphate reagent was added slowly with continuous shaking. The color of the reaction mixture was noted down visually within 15 sec of incubation and also recorded spectrophotometrically at 585 and 635 nm against air. An approximate COD value of the sample was estimated from the color of standard reaction mixture prepared with glucose as reference material in the COD range of 200 to 100000 mg/L. OD 585 and OD 635 of the reaction mixture was read spectrophotometrically. OD 585 and 635 nm in the range of 0.130 to 0.175 and 0.120 to 0.160 respectively were considered and calculated with reference to glucose standard. (Glucose COD in the range of 320 to 535 mg/L corresponds to OD 585 and 635 nm in the range of 0.157 to 0.207 and 0.142 to 0.189, respectively).

COD of the sample was estimated to be 6210 mg/L. (By Merck method and read on photometer).

Based on the color of the reaction mixture, COD of the sample was estimated to be around 4000 mg/L. On the basis of the dilution method and the spectrophotometric reading at OD 585 and OD 635 of the sample, COD was estimated to be 4270 mg/l and 3370mg/l. A comparison with the COD value calculated through Merck method reveals the following relation i.e. Our COD value (at OD 585) $\times 1.45 = \text{COD value (Merck Method)}$. or our COD value (at OD 635) $\times 1.84 = \text{COD (Merck Method)}$.

Hence for a preliminary estimate of COD value color of the standard reaction mixture based on a range of glucose concentrations can be used as a reference. For samples, with COD concentration of 300 to 500 mg/L, OD 585 or OD 635 can be recorded for obtaining more precise values.

Table 37: Correlation between COD concentration, color of the reaction mixture and their OD at different wave lengths.

COD (mg/L)	Sample dilution (Times)	Color of the reaction mixture	OD	Estimated COD (mg/L)	Estimated final COD (mg/L) (Col. 2x5)
1	2	3	4	5	6
<u>Values at OD 585 nm</u>					
Glucose:					
320	NA	C300	0.165	NA	NA
426	NA	C400	0.193	NA	NA
535	NA	C550	0.213	NA	NA
Sample: Damaged wheat grain (as such)					
	Nil	C1100	0.491	≈1100	≈1100#
	10	C400	0.154	427	4270
<u>Values at OD 635 nm</u>					
Glucose:					
320	NA	C300	0.148	NA	NA
426	NA	C400	0.172	NA	NA
535	NA	C550	0.191	NA	NA
Sample: Damaged wheat grain (as such)					
	Nil	C1100	0.401	≈1100	≈1100#
	10	C400	0.130	337	3370

#: COD as deducible from Table 2 given in Example 2. D: Dark. T: Turbid.

*: OD out of range. NA: Not applicable. Samples are diluted to get an OD 585 in the range of 0.130 to 0.175 and OD 635 in the range of 0.120 to 0.160, to obtain a COD value in the range of 320 to 535 mg/L. (Based on the Tables 3 to 15 given in Examples 3 to 15).

EXAMPLE 38:

Dissolved 0.3 g glucose in 50 mL distilled water. Aliquots from the stock solution were taken and diluted with water to achieve different COD concentrations in the range of 320 to 535 mg/L. A 30 mL sample of Damaged wheat grain (ssp) was taken. Weighed 0.02 g HgSO_4 in separate 15 mL test tubes. 1.0 mL of glucose solution and sample were added to each test tube separately and mixed well. 0.5 mL

of $K_2Cr_2O_7$ was added to each sample. 1.5 mL H_2SO_4 silver sulphate reagent was added slowly with continuous shaking. The color of the reaction mixture was noted down visually within 15 sec of incubation and also recorded spectrophotometrically at 585 and 635 nm against air. An approximate COD value of the sample was estimated from the color of standard reaction mixture prepared with glucose as reference material in the COD range of 200 to 100000 mg/L. OD 585 and OD 635 of the reaction mixture was read spectrophotometrically. OD 585 and 635 nm in the range of 0.130 to 0.175 and 0.120 to 0.160 respectively were considered and calculated with reference to glucose standard. (Glucose COD in the range of 320 to 535 mg/L corresponds to OD 585 and 635 nm in the range of 0.157 to 0.207 and 0.142 to 0.189, respectively).

COD of the sample was estimated to be 4710.4 mg/L. (By Merck method and read on photometer).

Based on the color of the reaction mixture, COD of the sample was estimated to be around 4800 mg/L. On the basis of the dilution method and the spectrophotometric reading at OD 585 and OD 635 of the sample, COD was estimated to be 5220 and 4998 mg/l respectively. A comparison with the COD value calculated through Merck method reveals the following relation i.e. Our COD value (at OD 585) $\times 0.90 =$ COD value (Merck Method) or our COD value (at OD 635) $\times 0.94 =$ COD (Merck Method). Hence for a preliminary estimate of COD value color of the standard reaction mixture based on a range of glucose concentrations can be used as a reference. For samples, with COD concentration of 300 to 500 mg/L, OD 585 or OD 635 can be recorded for obtaining more precise values.

Table 38: Correlation between COD concentration, color of the reaction mixture and their OD at different wave lengths.

COD (mg/L)	Sample dilution (Times)	Color of the reaction mixture	OD	Estimated COD (mg/L)	Estimated final COD (mg/L) (Col. 2x5)
1	2	3	4	5	6
<u>Values at OD 585 nm</u>					
Glucose:					
320	NA	C300	0.163	NA	NA
426	NA	C400	0.176	NA	NA
535	NA	C550	0.226	NA	NA
Sample: Damaged wheat grain (ssp)					
	Nil	C1500	0.652	1500	≅1500
	10	C200	0.187	224	≅2240 #
	12	C400	0.147	435	5220
<u>Values at OD 635 nm</u>					
Glucose:					
320	NA	C300	0.146	NA	NA
426	NA	C400	0.159	NA	NA
535	NA	C550	0.205	NA	NA
Sample: Damaged wheat grain (ssp)					
	Nil	C1500	0.551	1500	≅1500
	10	C200	0.168	224	≅ 2240#
	12	C400	0.127	416	4998

#: COD as deducible from Table 2 given in Example 2. D: Dark. T: Turbid.

*: OD out of range. NA: Not applicable. Samples are diluted to get an OD 585 in the range of 0.130 to 0.175 and OD 635 in the range of 0.120 to 0.160, to obtain a COD value in the range of 320 to 535 mg/L. (Based on the Tables 3 to 15 given in Examples 3 to 15).

EXAMPLE 39:

Dissolved 0.15 g glucose in 25 mL distilled water. Aliquots from the stock solution were taken and diluted with water to achieve different COD concentrations in the range of 320 to 535 mg/L. A 30 mL sample of Palm oil mill effluent (as such) was

taken. Weighed 0.02 g HgSO_4 in separate 15 mL test tubes. 1.0 mL of glucose solution and sample were added to each test tube separately and mixed well. 0.5 mL of $\text{K}_2\text{Cr}_2\text{O}_7$ was added to each sample. 1.5 mL H_2SO_4 silver sulphate reagent was added slowly with continuous shaking. The color of the reaction mixture was noted down visually within 15 sec of incubation and also recorded spectrophotometrically at 585 and 635 nm against air. An approximate COD value of the sample was estimated from the color of standard reaction mixture prepared with glucose as reference material in the COD range of 200 to 100000 mg/L. OD 585 and OD 635 of the reaction mixture was read spectrophotometrically. OD 585 and 635 nm in the range of 0.130 to 0.175 and 0.120 to 0.160 respectively were considered and calculated with reference to glucose standard. (Glucose COD in the range of 320 to 535 mg/L corresponds to OD 585 and 635 nm in the range of 0.157 to 0.207 and 0.142 to 0.189, respectively).

COD of the sample was estimated to be 15,801 mg/L. (By Merck method and read on photometer).

Based on the color of the reaction mixture, COD of the sample was estimated to be around 12000 mg/L. On the basis of the dilution method and the spectrophotometric reading at OD 585 and OD 635 of the sample, COD was estimated to be 13389 and 13506 mg/l respectively. A comparison with the COD value calculated through Merck method reveals the following relation i.e. Our COD value at (OD 585) $\times 1.18 =$ COD value (Merck Method) or our COD value (at OD 635) $\times 1.17 =$ COD value (Merck Method).

Hence for a preliminary estimate of COD value color of the standard reaction mixture based on a range of glucose concentrations can be used as a reference. For samples, with COD concentration of 300 to 500 mg/L, OD 585 or OD 635 can be recorded for obtaining more precise values.

Table 39: Correlation between COD concentration, color of the reaction mixture and their OD at different wave lengths.

COD (mg/L)	Sample dilution (Times)	Color of the reaction mixture	OD	Estimated COD (mg/L)	Estimated final COD (mg/L) (Col. 2x5)
1	2	3	4	5	6
<u>Values at OD 585 nm</u>					
Glucose:					
320	NA	C300	0.165	NA	NA
426	NA	C400	0.193	NA	NA
535	NA	C550	0.213	NA	NA
Sample: Palm oil Mill effluent (as such)					
	Nil	C10000	-	>10,000	>10,000
	10	C1900	0.394	≈1900	≈19,000#
	30	C400	0.165	446.3	13389
<u>Values at OD 635 nm</u>					
Glucose:					
320	NA	C300	0.148	NA	NA
426	NA	C400	0.172	NA	NA
535	NA	C550	0.191	NA	NA
Sample: Palm oil Mill effluent (as such)					
	Nil	C10000	-	>10000	>10000
	10	C1900	0.360	≈ 1900	≈19000#
	30	C400	0.148	450.2	13506

#: COD as deducible from Table 2 given in Example 2. D: Dark. T: Turbid.

*: OD out of range. NA: Not applicable. Samples are diluted to get an OD 585 in the range of 0.130 to 0.175 and OD 635 in the range of 0.120 to 0.160, to obtain a COD value in the range of 320 to 535 mg/L. (Based on the Tables 3 to 15 given in Examples 3 to 15).

EXAMPLE 40:

Dissolved 0.15 g glucose in 25 mL distilled water. Aliquots from the stock solution were taken and diluted with water to achieve different COD concentrations in the range of 320 to 535 mg/L. A 30 mL sample of Palm oil mill effluent (ssp) was taken.

Weighed 0.02 g HgSO_4 in separate 15 mL test tubes. 1.0 mL of glucose solution and sample were added to each test tube separately and mixed well. 0.5 mL of $\text{K}_2\text{Cr}_2\text{O}_7$ was added to each sample. 1.5 mL H_2SO_4 silver sulphate reagent was added slowly with continuous shaking. The color of the reaction mixture was noted down visually within 15 sec of incubation and also recorded spectrophotometrically at 585 and 635 nm against air. An approximate COD value of the sample was estimated from the color of standard reaction mixture prepared with glucose as reference material in the COD range of 200 to 100000 mg/L. OD 585 and OD 635 of the reaction mixture was read spectrophotometrically. OD 585 and 635 nm in the range of 0.130 to 0.175 and 0.120 to 0.160 respectively were considered and calculated with reference to glucose standard. (Glucose COD in the range of 320 to 535 mg/L corresponds to OD 585 and 635 nm in the range of 0.157 to 0.207 and 0.142 to 0.189, respectively). COD of the sample was estimated to be 1129.6 mg/L. (By Merck method and read on photometer).

Based on the color of the reaction mixture, COD of the sample was estimated to be around 800 mg/L. On the basis of the dilution method and the spectrophotometric reading at OD 585 and OD 635 of the sample, COD was estimated to be 970 and 950 mg/l respectively. A comparison with the COD value calculated through Merck method reveals the following relation i.e. Our COD value (at OD 585) $\times 1.16 =$ COD value (Merck Method) or our COD value (at OD 635) $\times 1.19 =$ COD value (Merck Method).

Hence for a preliminary estimate of COD value color of the standard reaction mixture based on a range of glucose concentrations can be used as a reference. For samples, with COD concentration of 300 to 500 mg/L, OD 585 or OD 635 can be recorded for obtaining more precise values.

Table 40: Correlation between COD concentration, color of the reaction mixture and their OD at different wave lengths.

COD (mg/L)	Sample dilution (Times)	Color of the reaction mixture	OD	Estimated COD (mg/L)	Estimated final COD (mg/L) (Col. 2x5)
1	2	3	4	5	6
<u>Values at OD 585 nm</u>					
Glucose:					
320	NA	C300	0.163	NA	NA
426	NA	C400	0.176	NA	NA
535	NA	C550	0.226	NA	NA
Sample: Palm oil Mill effluent (ssp)					
	Nil	C1500	0.263	≈1500	≈1500#
	2	C400	0.171	485	970
<u>Values at OD 635 nm</u>					
Glucose:					
320	NA	C300	0.146	NA	NA
426	NA	C400	0.159	NA	NA
535	NA	C550	0.205	NA	NA
Sample: Palm oil Mill effluent (ssp)					
	Nil	C1500	0.242	≈1500	≈1500#
	2	C400	0.152	475	950

#: COD as deducible from Table 2 given in Example 2. D: Dark. T: Turbid.

*: OD out of range. NA: Not applicable. Samples are diluted to get an OD 585 in the range of 0.130 to 0.175 and OD 635 in the range of 0.120 to 0.160, to obtain a COD value in the range of 320 to 535 mg/L. (Based on the Tables 3 to 15 given in Examples 3 to 15).

EXAMPLE 41:

Dissolved 0.15 g glucose in 25 mL distilled water. Aliquots from the stock solution were taken and diluted with water to achieve different COD concentrations in the range of 320 to 535 mg/L. A 30 mL sample of Damaged wheat grain and POME (as such) was taken. Weighed 0.02 g HgSO_4 in separate 15 mL test tubes. 1.0 mL of glucose solution and sample were added to each test tube separately and mixed well. 0.5 mL of $\text{K}_2\text{Cr}_2\text{O}_7$ was added to each sample. 1.5 mL H_2SO_4 silver sulphate reagent was added slowly with continuous shaking. The color of the reaction mixture

was noted down visually within 15 sec of incubation and also recorded spectrophotometrically at 585 and 635 nm against air. An approximate COD value of the sample was estimated from the color of standard reaction mixture prepared with glucose as reference material in the COD range of 200 to 100000 mg/L. OD 585 and OD 635 of the reaction mixture was read spectrophotometrically. OD 585 and 635 nm in the range of 0.130 to 0.175 and 0.120 to 0.160 respectively were considered and calculated with reference to glucose standard. (Glucose COD in the range of 320 to 535 mg/L corresponds to OD 585 and 635 nm in the range of 0.157 to 0.207 and 0.142 to 0.189, respectively).

COD of the sample was estimated to be 9706 mg/L. (By Merck method and read on photometer).

Based on the color of the reaction mixture, COD of the sample was estimated to be around 9000 mg/L. On the basis of the dilution method and the spectrophotometric reading at OD 585 and OD 635 of the sample, COD was estimated to be 7660 and 7567 mg/l respectively. A comparison with the COD value calculated through Merck method reveals the following relation i.e. Our COD value (at OD 585) $\times 1.27$ = COD value (Merck Method) or our COD value (at OD 635) $\times 1.12$.

Hence for a preliminary estimate of COD value color of the standard reaction mixture based on a range of glucose concentrations can be used as a reference. For samples, with COD concentration of 300 to 500 mg/L, OD 585 or OD 635 can be recorded for obtaining more precise values.

Table 41: Correlation between COD concentration, color of the reaction mixture and their OD at different wave lengths.

COD (mg/L)	Sample dilution (Times)	Color of the reaction mixture	OD	Estimated COD (mg/L)	Estimated final COD (mg/L) (Col. 2x5)
1	2	3	4	5	6
<u>Values at OD 585 nm</u>					
Glucose:					
320	NA	C300	0.165	NA	NA
426	NA	C400	0.193	NA	NA
535	NA	C550	0.213	NA	NA
Sample: DWS (S3) + POME (as such)					
	Nil	C10000	1.231	>10,000	>10,000
	10	C1500	0.256	≅1500	≅15,000
	20	C450	0.186	≅450	≅9000 #
	25	C300	0.132	306.4	7660
<u>Values at OD 635 nm</u>					
Glucose:					
320	NA	C300	0.148	NA	NA
426	NA	C400	0.172	NA	NA
535	NA	C550	0.191	NA	NA
Sample: DWS (S3) + POME (as such)					
	Nil	C10000	0.994	10,000	>10,000
	10	C1500	0.221	≅1500	≅15,000
	20	C450	0.165	≅450	≅9000 #
	25	C300	0.117	302.6	7567

#: COD as deducible from Table 2 given in Example 2. D: Dark. T: Turbid.

*: OD out of range. NA: Not applicable. Samples are diluted to get an OD 585 in the range of 0.130 to 0.175 and OD 635 in the range of 0.120 to 0.160, to obtain a COD value in the range of 320 to 535 mg/L. (Based on the Tables 3 to 15 given in Examples 3 to 15).

EXAMPLE 42:

Dissolved 0.15 g glucose in 25 mL distilled water. Aliquots from the stock solution were taken and diluted with water to achieve different COD concentrations in the range of 320 to 535 mg/L. A 30 mL sample of Damaged wheat grain and POME (ssp) was taken. Weighed 0.02 g HgSO_4 in separate 15 mL test tubes. 1.0 mL of

glucose solution and sample were added to each test tube separately and mixed well. 0.5 mL of $K_2Cr_2O_7$ was added to each sample. 1.5 mL H_2SO_4 silver sulphate reagent was added slowly with continuous shaking. The color of the reaction mixture was noted down visually within 15 sec of incubation and also recorded spectrophotometrically at 585 and 635 nm against air. An approximate COD value of the sample was estimated from the color of standard reaction mixture prepared with glucose as reference material in the COD range of 200 to 100000 mg/L. OD 585 and OD 635 of the reaction mixture was read spectrophotometrically. OD 585 and 635 nm in the range of 0.130 to 0.175 and 0.120 to 0.160 respectively were considered and calculated with reference to glucose standard. (Glucose COD in the range of 320 to 535 mg/L corresponds to OD 585 and 635 nm in the range of 0.157 to 0.207 and 0.142 to 0.189, respectively).

COD of the sample was estimated to be 4427.5 mg/L. (By Merck method and read on photometer).

Based on the color of the reaction mixture, COD of the sample was estimated to be around 4000 mg/L. On the basis of the dilution method and the spectrophotometric reading at OD 585 and OD 635 of the sample, COD was estimated to be 4068 and 3936 mg/l respectively. A comparison with the COD value calculated through Merck method reveals the following relation i.e. Our COD value (at OD 585) $\times 1.09$ = COD value (Merck Method) or our COD value (at OD 635) $\times 1.12$.

Hence for a preliminary estimate of COD value color of the standard reaction mixture based on a range of glucose concentrations can be used as a reference. For samples, with COD concentration of 300 to 500 mg/L, OD 585 or OD 635 can be recorded for obtaining more precise values.

DWS (S3) + POME (ssp)

Table 42: Correlation between COD concentration, color of the reaction mixture and their OD at different wave lengths.

COD (mg/L)	Sample dilution (Times)	Color of the reaction mixture	OD	Estimated COD (mg/L)	Estimated final COD (mg/L) (Col. 2x5)
1	2	3	4	5	6
<u>Values at OD 585 nm</u>					
Glucose:					
320	NA	C300	0.165	NA	NA
426	NA	C400	0.193	NA	NA
535	NA	C550	0.213	NA	NA
Sample: DWS (S3) + POME (ssp)					
	Nil	C1500	0.520	>1500	>1500
	10	C400	0.185	≅400	≅4000#
	12	C300	0.146	339	4068
<u>Values at OD 635 nm</u>					
Glucose:					
320	NA	C300	0.148	NA	NA
426	NA	C400	0.172	NA	NA
535	NA	C550	0.191	NA	NA
Sample: DWS (S3) + POME (ssp)					
	Nil	C1500	0.427	>1500	>1500
	10	C400	0.165	≅400	≅4000#
	12	C300	0.127	328	3936

#: COD as deducible from Table 2 given in Example 2. D: Dark. T: Turbid.

*: OD out of range. NA: Not applicable. Samples are diluted to get an OD 585 in the range of 0.130 to 0.175 and OD 635 in the range of 0.120 to 0.160, to obtain a COD value in the range of 320 to 535 mg/L. (Based on the Tables 3 to 15 given in Examples 3 to 15).

EXAMPLE 43:

Dissolved 0.15 g glucose in 25 mL distilled water. Aliquots from the stock solution were taken and diluted with water to achieve different COD concentrations in the range of 320 to 535 mg/L. A 30 mL sample of Tea (as such) was taken. Weighed

0.02 g HgSO_4 in separate 15 mL test tubes. 1.0 mL of glucose solution and sample were added to each test tube separately and mixed well. 0.5 mL of $\text{K}_2\text{Cr}_2\text{O}_7$ was added to each sample. 1.5 mL H_2SO_4 silver sulphate reagent was added slowly with continuous shaking. The color of the reaction mixture was noted down visually within 15 sec of incubation and also recorded spectrophotometrically at 585 and 635 nm against air. An approximate COD value of the sample was estimated from the color of standard reaction mixture prepared with glucose as reference material in the COD range of 200 to 100000 mg/L. OD 585 and OD 635 of the reaction mixture was read spectrophotometrically. OD 585 and 635 nm in the range of 0.130 to 0.175 and 0.120 to 0.160 respectively were considered and calculated with reference to glucose standard. (Glucose COD in the range of 320 to 535 mg/L corresponds to OD 585 and 635 nm in the range of 0.157 to 0.207 and 0.142 to 0.189, respectively). COD of the sample was estimated to be 74,589 mg/L. (By Merck method and read on photometer).

Based on the color of the reaction mixture, COD of the sample was estimated to be around 90000 mg/L. On the basis of the dilution method and the spectrophotometric reading at OD 585 and OD 635 of the sample, COD was estimated to be 101750 and 102000 mg/l respectively. A comparison with the COD value calculated through Merck method reveals the following relation i.e. Our COD value at (OD 585 or OD 635) $\times 0.73 = \text{COD value (Merck Method)}$.

Hence, for a preliminary estimate of COD value color of the standard reaction mixture based on a range of glucose concentrations can be used as a reference. For samples, with COD concentration of 300 to 500 mg/L, OD 585 or OD 635 can be recorded for obtaining more precise values.

Tea (as such)

Table 43: Correlation between COD concentration, color of the reaction mixture and their OD at different wave lengths.

COD (mg/L)	Sample dilution (Times)	Color of the reaction mixture	OD	Estimated COD (mg/L)	Estimated final COD (mg/L) (Col. 2x5)
1	2	3	4	5	6
<u>Values at OD 585 nm</u>					
Glucose:					
320	NA	C300	0.159	NA	NA
426	NA	C400	0.186	NA	NA
535	NA	C550	0.223	NA	NA
Sample: TEA (as such)					
	Nil	C10000	-	>10,000	>10,000
	10	C10000	-	>10,000	>1,00,000
	100	C1900	0.267	≅1900	≅1,90,000
	111.1	C1300	0.245	≅1300	≅1,44,430
	166.7	C500	0.206	≅550	≅91,685
	200	C400	0.175	≅450	≅90,000#
	250	C300	0.145	407	101750
<u>Values at OD 635 nm</u>					
Glucose:					
320	NA	C300	0.141	NA	NA
426	NA	C400	0.164	NA	NA
535	NA	C550	0.198	NA	NA
Sample: TEA (as such)					
	Nil	C10000	-	>10,000	>10,000
	10	C10000	-	>10,000	>1,00,000
	100	C1900	0.237	≅1900	≅1,90,000
	111.1	C1300	0.214	≅1300	≅1,44,430
	166.7	C500	0.182	≅550	≅91,685
	200	C400	0.155	≅450	≅90,000#
	250	C300	0.128	408	102000

#: COD as deducible from Table 2 given in Example 2. D: Dark. T: Turbid.

*: OD out of range. NA: Not applicable. Samples are diluted to get an OD 585 in the range of 0.130 to 0.175 and OD 635 in the range of 0.120 to 0.160, to obtain

a COD value in the range of 320 to 535 mg/L. (Based on the Tables 3 to 15 given in Examples 3 to 15).

EXAMPLE 44:

Dissolved 0.15 g glucose in 25 mL distilled water. Aliquots from the stock solution were taken and diluted with water to achieve different COD concentrations in the range of 320 to 535 mg/L. A 30 mL sample of Tea (ssp) was taken. Weighed 0.02 g HgSO_4 in separate 15 mL test tubes. 1.0 mL of glucose solution and sample were added to each test tube separately and mixed well. 0.5 mL of $\text{K}_2\text{Cr}_2\text{O}_7$ was added to each sample. 1.5 mL H_2SO_4 silver sulphate reagent was added slowly with continuous shaking. The color of the reaction mixture was noted down visually within 15 sec of incubation and also recorded spectrophotometrically at 585 and 635 nm against air. An approximate COD value of the sample was estimated from the color of standard reaction mixture prepared with glucose as reference material in the COD range of 200 to 100000 mg/L. OD 585 and OD 635 of the reaction mixture was read spectrophotometrically. OD 585 and 635 nm in the range of 0.130 to 0.175 and 0.120 to 0.160 respectively were considered and calculated with reference to glucose standard. (Glucose COD in the range of 320 to 535 mg/L corresponds to OD 585 and 635 nm in the range of 0.157 to 0.207 and 0.142 to 0.189, respectively). COD of the sample was estimated to be 67,528 mg/L. (By Merck method and read on photometer).

Based on the color of the reaction mixture, COD of the sample was estimated to be around 75015 mg/L. On the basis of the dilution method and the spectrophotometric reading at OD 585 and OD 635 of the sample, COD was estimated to be 81380 and 18380 mg/l respectively. A comparison with the COD value calculated through Merck method reveals the following relation i.e. Our COD value (at OD 585 or OD 635) $\times 0.83 = \text{COD value (Merck Method)}$.

Hence for a preliminary estimate of COD value color of the standard reaction mixture based on a range of glucose concentrations can be used as a reference. For samples, with COD concentration of 300 to 500 mg/L, OD 585 or OD 635 can be recorded for obtaining more precise values.

Tea (ssp)

Table 44: Correlation between COD concentration, color of the reaction mixture and their OD at different wave lengths.

COD (mg/L)	Sample dilution (Times)	Color of the reaction mixture	OD	Estimated COD (mg/L)	Estimated final COD (mg/L) (Col. 2x5)
1	2	3	4	5	6
<u>Values at OD 585 nm</u>					
Glucose:					
320	NA	C300	0.159	NA	NA
426	NA	C400	0.186	NA	NA
535	NA	C550	0.223	NA	NA
Sample: TEA (ssp)					
	Nil	C10000	-	>10,000	≥>10,000
	10	C10000	-	>10,000	>1,00,000
	100	C1800	0.244	≈1800	≈1,80,000
	125	C550	0.222	≈550	≈68,750
	166.7	C450	0.173	≈450	≈75,015#
	200	C300	0.145	407	81380
<u>Values at OD 635 nm</u>					
Glucose:					
320	NA	C300	0.141	NA	NA
426	NA	C400	0.164	NA	NA
535	NA	C550	0.198	NA	NA
Sample: TEA(ssp)					
	Nil	C10000	-	>10,000	≥>10,000
	10	C10000	-	>10,000	>1,00,000
	100	C1800	0.216	≈1800	≈1,80,000
	125	C550	0.190	≈550	≈68,750
	166.7	C450	0.152	≈450	≈75,015#
	200	C300	0.128	408	81700

#: COD as deducible from Table 2 given in Example 2. D: Dark. T: Turbid.

*: OD out of range. NA: Not applicable. Samples are diluted to get an OD 585 in the range of 0.130 to 0.175 and OD 635 in the range of 0.120 to 0.160, to obtain a COD value in the range of 320 to 535 mg/L. (Based on the Tables 3 to 15 given in Examples 3 to 15).

EXAMPLE 45:

Dissolved 0.15 g glucose in 25 mL distilled water. Aliquots from the stock solution were taken and diluted with water to achieve different COD concentrations in the range of 320 to 535 mg/L. A 30 mL sample of Apple pomace (ssp) was taken. Weighed 0.02 g HgSO_4 in separate 15 mL test tubes. 1.0 mL of glucose solution and sample were added to each test tube separately and mixed well. 0.5 mL of $\text{K}_2\text{Cr}_2\text{O}_7$ was added to each sample. 1.5 mL H_2SO_4 silver sulphate reagent was added slowly with continuous shaking. The color of the reaction mixture was noted down visually within 15 sec of incubation and also recorded spectrophotometrically at 585 and 635 nm against air. An approximate COD value of the sample was estimated from the color of standard reaction mixture prepared with glucose as reference material in the COD range of 200 to 100000 mg/L. OD 585 and OD 635 of the reaction mixture was read spectrophotometrically. OD 585 and 635 nm in the range of 0.130 to 0.175 and 0.120 to 0.160 respectively were considered and calculated with reference to glucose standard. (Glucose COD in the range of 320 to 535 mg/L corresponds to OD 585 and 635 nm in the range of 0.157 to 0.207 and 0.142 to 0.189, respectively). COD of the sample was estimated to be 39,284 mg/L. (By Merck method and read on photometer).

Based on the color of the reaction mixture, COD of the sample was estimated to be around 30760 mg/L. On the basis of the dilution method and the spectrophotometric reading at OD 585 and OD 635 of the sample, COD was estimated to be 43000 and 43500 mg/l respectively. A comparison with the COD value calculated through Merck method reveals the following relation i.e. Our COD value (at OD 585) $\times 0.91 =$ COD value (Merck Method) or our COD value (at OD 635) $\times 435 =$ COD value (Merck Method).

Hence for a preliminary estimate of COD value color of the standard reaction mixture based on a range of glucose concentrations can be used as a reference. For samples, with COD concentration of 300 to 500 mg/L, OD 585 or OD 635 can be recorded for obtaining more precise values.

Apple Pomace (ssp)

Table 45: Correlation between COD concentration, color of the reaction mixture and their OD at different wave lengths.

COD (mg/L)	Sample dilution (Times)	Color of the reaction mixture	OD	Estimated COD (mg/L)	Estimated final COD (mg/L) (Col. 2x5)
1	2	3	4	5	6
<u>Values at OD 585 nm</u>					
Glucose:					
320	NA	C300	0.159	NA	NA
426	NA	C400	0.187	NA	NA
535	NA	C550	0.210	NA	NA
Sample: APPLE POMACE (ssp)					
	Nil	C10000	-	>10,000	>10,000
	10	C1600	0.392	≅1600	≅16,000
	16.7	C800	0.354	≅800	≅13,360
	50	C550	0.244	≅550	≅27,500
	76.9	C400	0.178	≅400	≅30,760#
	100	C300	0.154	430	43000
<u>Values at OD 635 nm</u>					
Glucose:					
320	NA	C300	0.141	NA	NA
426	NA	C400	0.167	NA	NA
535	NA	C550	0.191	NA	NA
Sample: APPLE POMACE (ssp)					
	Nil	C10000	-	>10,000	>10,000
	10	C1600	0.336	≅1600	≅16,000
	16.7	C800	0.317	≅800	≅13,360
	50	C550	0.217	≅550	≅27,500
	76.9	C400	0.159	≅400	≅30,760#
	100	C300	0.139	435	45000

#: COD as deducible from Table 2 given in Example 2. D: Dark. T: Turbid.

*: OD out of range. NA: Not applicable. Samples are diluted to get an OD 585 in the range of 0.130 to 0.175 and OD 635 in the range of 0.120 to 0.160, to obtain a COD value in the range of 320 to 535 mg/L. (Based on the Tables 3 to 15 given in Examples 3 to 15).

EXAMPLE 46:

Dissolved 0.15 g glucose in 25 mL distilled water. Aliquots from the stock solution were taken and diluted with water to achieve different COD concentrations in the range of 320 to 535 mg/L. A 30 mL sample of Tea and Apple pomace (as such) was taken. Weighed 0.02 g HgSO_4 in separate 15 mL test tubes. 1.0 mL of glucose solution and sample were added to each test tube separately and mixed well. 0.5 mL of $\text{K}_2\text{Cr}_2\text{O}_7$ was added to each sample. 1.5 mL H_2SO_4 silver sulphate reagent was added slowly with continuous shaking. The color of the reaction mixture was noted down visually within 15 sec of incubation and also recorded spectrophotometrically at 585 and 635 nm against air. An approximate COD value of the sample was estimated from the color of standard reaction mixture prepared with glucose as reference material in the COD range of 200 to 100000 mg/L. OD 585 and OD 635 of the reaction mixture was read spectrophotometrically. OD 585 and 635 nm in the range of 0.130 to 0.175 and 0.120 to 0.160 respectively were considered and calculated with reference to glucose standard. (Glucose COD in the range of 320 to 535 mg/L corresponds to OD 585 and 635 nm in the range of 0.157 to 0.207 and 0.142 to 0.189, respectively).

COD of the sample was estimated to be 56,948 mg/L. (By Merck method and read on photometer).

Based on the color of the reaction mixture, COD of the sample was estimated to be around 40000 mg/L. On the basis of the dilution method and the spectrophotometric reading at OD 585 and OD 635 of the sample, COD was estimated to be 68441 and 681334 mg/l respectively. A comparison with the COD value calculated through Merck method reveals the following relation i.e. Our COD value (at OD 585) $\times 0.83 =$ COD value (Merck Method) or our COD value (at OD 635) $\times 0.84 =$ COD value (Merck Method.)

Hence for a preliminary estimate of COD value color of the standard reaction mixture based on a range of glucose concentrations can be used as a reference. For samples, with COD concentration of 300 to 500 mg/L, OD 585 or OD 635 can be recorded for obtaining more precise values.



Tea + Apple Pomace (as such)

Table 46: Correlation between COD concentration, color of the reaction mixture and their OD at different wave lengths.

COD (mg/L)	Sample dilution (Times)	Color of the reaction mixture	OD	Estimated COD (mg/L)	Estimated final COD (mg/L) (Col. 2x5)
1	2	3	4	5	6

Values at OD 585 nm

Glucose:

320	NA	C300	0.157	NA	NA
426	NA	C400	0.176	NA	NA
535	NA	C550	0.202	NA	NA

Sample: TEA + APPLE POMACE (as such)

Nil	C10000	-	>10,000	>10,000
10	C1500	0.603	≅1500	≅15,000
50	C700	0.295	≅700	≅35,000
100	C400	0.187	≅400	≅40,000#
153.8	C450	0.150	445	68441

Values at OD 635 nm

Glucose:

320	NA	C300	0.142	NA	NA
426	NA	C400	0.155	NA	NA
535	NA	C550	0.182	NA	NA

Sample: TEA + APPLE POMACE (as such)

Nil	C10000	-	>10,000	>10,000
10	C1500	0.465	≅1500	≅15,000
50	C700	0.262	≅700	≅35,000
100	C400	0.168	≅400	≅40,000#
153.8	C400	0.122	443	68134

#: COD as deducible from Table 2 given in Example 2. D: Dark. T: Turbid.

*: OD out of range. NA: Not applicable. Samples are diluted to get an OD 585 in the range of 0.130 to 0.175 and OD 635 in the range of 0.120 to 0.160, to obtain a COD value in the range of 320 to 535 mg/L. (Based on the Tables 3 to 15 given in Examples 3 to 15).

EXAMPLE 47:

Dissolved 0.15 g glucose in 25 mL distilled water. Aliquots from the stock solution were taken and diluted with water to achieve different COD concentrations in the range of 320 to 535 mg/L. A 30 mL sample of Tea + Apple pomace (ssp) and POME was taken. Weighed 0.02 g HgSO_4 in separate 15 mL test tubes. 1.0 mL of glucose solution and sample were added to each test tube separately and mixed well. 0.5 mL of $\text{K}_2\text{Cr}_2\text{O}_7$ was added to each sample. 1.5 mL H_2SO_4 silver sulphate reagent was added slowly with continuous shaking. The color of the reaction mixture was noted down visually within 15 sec of incubation and also recorded spectrophotometrically at 585 and 635 nm against air. An approximate COD value of the sample was estimated from the color of standard reaction mixture prepared with glucose as reference material in the COD range of 200 to 100000 mg/L. OD 585 and OD 635 of the reaction mixture was read spectrophotometrically. OD 585 and 635 nm in the range of 0.130 to 0.175 and 0.120 to 0.160 respectively were considered and calculated with reference to glucose standard. (Glucose COD in the range of 320 to 535 mg/L corresponds to OD 585 and 635 nm in the range of 0.157 to 0.207 and 0.142 to 0.189, respectively).

COD of the sample was estimated to be 60,398 mg/L. (By Merck method and read on photometer).

Based on the color of the reaction mixture, COD of the sample was estimated to be around 45000 mg/L. On the basis of the dilution method and the spectrophotometric reading at OD 585 and OD 635 of the sample, COD was estimated to be 66,134 and 67,364 mg/l respectively. A comparison with the COD value calculated through Merck method reveals the following relation i.e. Our COD value (at OD 585) $\times 0.91 =$ COD value (Merck Method) or our COD value (at OD 635) $\times 0.90 =$ COD value (Merck Method.)

Hence for a preliminary estimate of COD value color of the standard reaction mixture based on a range of glucose concentrations can be used as a reference. For samples, with COD concentration of 300 to 500 mg/L, OD 585 or OD 635 can be recorded for obtaining more precise values.

Tea + Apple Pomace (ssp)

Table 47: Correlation between COD concentration, color of the reaction mixture and their OD at different wave lengths.

COD (mg/L)	Sample dilution (Times)	Color of the reaction mixture	OD	Estimated COD (mg/L)	Estimated final COD (mg/L) (Col. 2x5)
1	2	3	4	5	6

Values at OD 585 nm

Glucose:

320	NA	C300	0.157	NA	NA
426	NA	C400	0.176	NA	NA
535	NA	C550	0.202	NA	NA

Sample: TEA + APPLE POMACE (ssp)

Nil	C10000	-	>10,000	>10,000
10	C1400	0.621	≅1400	≅14,000
50	C800	0.284	≅800	≅40,000
100	C450	0.181	≅450	≅45,000
153.8	C400	0.145	430	66,134

Values at OD 635 nm

Glucose:

320	NA	C300	0.142	NA	NA
426	NA	C400	0.155	NA	NA
535	NA	C550	0.182	NA	NA

Sample: TEA + APPLE POMACE (ssp)

Nil	C10000	-	>10,000	>10,000
10	C1400	0.473	≅1400	≅14,000
50	C800	0.284	≅800	≅40,000
100	C450	0.159	≅450	≅45,000
153.8	C400	0.130	438	67,364

#: COD as deducible from Table 2 given in Example 2. D: Dark. T: Turbid.

*: OD out of range. NA: Not applicable. Samples are diluted to get an OD 585 in the range of 0.130 to 0.175 and OD 635 in the range of 0.120 to 0.160, to obtain a COD value in the range of 320 to 535 mg/L. (Based on the Tables 3 to 15 given in Examples 3 to 15).

	1980	1981	1982	1983	1984	1985	1986	1987	1988	1989	1990	1991	1992	1993	1994	1995	1996	1997	1998	1999	2000	2001	2002	2003	2004	2005	2006	2007	2008	2009	2010	2011	2012	2013	2014	2015	2016	2017	2018	2019	2020	2021	2022	2023	2024	2025	2026	2027	2028	2029	2030	2031	2032	2033	2034	2035	2036	2037	2038	2039	2040	2041	2042	2043	2044	2045	2046	2047	2048	2049	2050	2051	2052	2053	2054	2055	2056	2057	2058	2059	2060	2061	2062	2063	2064	2065	2066	2067	2068	2069	2070	2071	2072	2073	2074	2075	2076	2077	2078	2079	2080	2081	2082	2083	2084	2085	2086	2087	2088	2089	2090	2091	2092	2093	2094	2095	2096	2097	2098	2099	2100	2101	2102	2103	2104	2105	2106	2107	2108	2109	2110	2111	2112	2113	2114	2115	2116	2117	2118	2119	2120	2121	2122	2123	2124	2125	2126	2127	2128	2129	2130	2131	2132	2133	2134	2135	2136	2137	2138	2139	2140	2141	2142	2143	2144	2145	2146	2147	2148	2149	2150	2151	2152	2153	2154	2155	2156	2157	2158	2159	2160	2161	2162	2163	2164	2165	2166	2167	2168	2169	2170	2171	2172	2173	2174	2175	2176	2177	2178	2179	2180	2181	2182	2183	2184	2185	2186	2187	2188	2189	2190	2191	2192	2193	2194	2195	2196	2197	2198	2199	2200	2201	2202	2203	2204	2205	2206	2207	2208	2209	2210	2211	2212	2213	2214	2215	2216	2217	2218	2219	2220	2221	2222	2223	2224	2225	2226	2227	2228	2229	2230	2231	2232	2233	2234	2235	2236	2237	2238	2239	2240	2241	2242	2243	2244	2245	2246	2247	2248	2249	2250	2251	2252	2253	2254	2255	2256	2257	2258	2259	2260	2261	2262	2263	2264	2265	2266	2267	2268	2269	2270	2271	2272	2273	2274	2275	2276	2277	2278	2279	2280	2281	2282	2283	2284	2285	2286	2287	2288	2289	2290	2291	2292	2293	2294	2295	2296	2297	2298	2299	2300	2301	2302	2303	2304	2305	2306	2307	2308	2309	2310	2311	2312	2313	2314	2315	2316	2317	2318	2319	2320	2321	2322	2323	2324	2325	2326	2327	2328	2329	2330	2331	2332	2333	2334	2335	2336	2337	2338	2339	2340	2341	2342	2343	2344	2345	2346	2347	2348	2349	2350	2351	2352	2353	2354	2355	2356	2357	2358	2359	2360	2361	2362	2363	2364	2365	2366	2367	2368	2369	2370	2371	2372	2373	2374	2375	2376	2377	2378	2379	2380	2381	2382	2383	2384	2385	2386	2387	2388	2389	2390	2391	2392	2393	2394	2395	2396	2397	2398	2399	2400	2401	2402	2403	2404	2405	2406	2407	2408	2409	2410	2411	2412	2413	2414	2415	2416	2417	2418	2419	2420	2421	2422	2423	2424	2425	2426	2427	2428	2429	2430	2431	2432	2
--	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	------	---

Based on the color of the reaction mixture, COD of the sample was estimated to be around 1700 mg/L. On the basis of the dilution method and the spectrophotometric reading at OD 585 and OD 635 of the sample, COD was estimated to be 3120 and 3160 mg/l respectively. A comparison with the COD value calculated through Merck method reveals the following relation i.e. Our COD value (at OD 585) $\times 1.92 =$ COD value (Merck Method) or our COD value (at OD 635) $\times 1.90 =$ COD value (Merck Method).)

104



Digested POME (as such)

Table 48: Correlation between COD concentration, color of the reaction mixture and their OD at different wave lengths.

COD (mg/L)	Sample dilution (Times)	Color of the reaction mixture	OD	Estimated COD (mg/L)	Estimated final COD (mg/L) (Col. 2x5)
1	2	3	4	5	6
<u>Values at OD 585 nm</u>					
Glucose:					
320	NA	C300	0.162	NA	NA
426	NA	C400	0.189	NA	NA
535	NA	C550	0.207	NA	NA
Sample: Digested POME (as such)					
	Nil	C1700	1.021	≅1700	≅1700#
	10	C250	0.132	312	3120
<u>Values at OD 635 nm</u>					
Glucose:					
320	NA	C300	0.146	NA	NA
426	NA	C400	0.171	NA	NA
535	NA	C300	0.187	NA	NA
Sample: Digested POME (as such)					
	Nil	C1700	0.889	≅1700	≅1700#
	10	C250	0.120	316	3160

#: COD as deducible from Table 2 given in Example 2. D: Dark. T: Turbid.

*: OD out of range. NA: Not applicable. Samples are diluted to get an OD 585 in the range of 0.130 to 0.175 and OD 635 in the range of 0.120 to 0.160, to obtain a COD value in the range of 320 to 535 mg/L. (Based on the Tables 3 to 15 given in Examples 3 to 15).

EXAMPLE 49:

Dissolved 0.15 g glucose in 25 mL distilled water. Aliquots from the stock solution were taken and diluted with water to achieve different COD concentrations in the range of 320 to 535 mg/L. A 30 mL sample of Digested POME (ssp) was taken. Weighed 0.02 g HgSO_4 in separate 15 mL test tubes. 1.0 mL of glucose solution and sample were added to each test tube separately and mixed well. 0.5 mL of $\text{K}_2\text{Cr}_2\text{O}_7$

was added to each sample. 1.5 mL H₂SO₄ silver sulphate reagent was added slowly with continuous shaking. The color of the reaction mixture was noted down visually within 15 sec of incubation and also recorded spectrophotometrically at 585 and 635 nm against air. An approximate COD value of the sample was estimated from the color of standard reaction mixture prepared with glucose as reference material in the COD range of 200 to 100000 mg/L. OD 585 and OD 635 of the reaction mixture was read spectrophotometrically. OD 585 and 635 nm in the range of 0.130 to 0.175 and 0.120 to 0.160 respectively were considered and calculated with reference to glucose standard. (Glucose COD in the range of 320 to 535 mg/L corresponds to OD 585 and 635 nm in the range of 0.157 to 0.207 and 0.142 to 0.189, respectively). COD of the sample was estimated to be 588.4 mg/L. (By Merck method and read on photometer).

Based on the color of the reaction mixture, COD of the sample was estimated to be around 400 mg/L. On the basis of the dilution method and the spectrophotometric reading at OD 585 and OD 635 of the sample, COD was estimated to be 535 and 458 mg/l respectively. A comparison with the COD value calculated through Merck method reveals the following relation i.e. Our COD value (at OD 585) x 1.10 = COD value (Merck Method) or our COD value (at OD 635) x 1.30 = COD value (Merck Method).

Hence for a preliminary estimate of COD value color of the standard reaction mixture based on a range of glucose concentrations can be used as a reference. For samples, with COD concentration of 300 to 500 mg/L, OD 585 or OD 635 can be recorded for obtaining more precise values.

Digested POME (ssp)

Table 49: Correlation between COD concentration, color of the reaction mixture and their OD at different wave lengths.

COD (mg/L)	Sample dilution (Times)	Color of the reaction mixture	OD	Estimated COD (mg/L)	Estimated final COD (mg/L) (Col. 2x5)
1	2	3	4	5	6
<u>Values at OD 585 nm</u>					
Glucose:					
320	NA	C300	0.162	NA	NA
426	NA	C400	0.189	NA	NA
535	NA	C550	0.207	NA	NA
Sample: Digested POME (ssp)					
	Nil	C400	0.168	535	535
<u>Values at OD 635 nm</u>					
Glucose:					
320	NA	C300	0.146	NA	NA
426	NA	C400	0.171	NA	NA
535	NA	C550	0.187	NA	NA
Sample: Digested POME (ssp)					
	Nil	C400	0.148	458	458

#: COD as deducible from Table 2 given in Example 2. D: Dark. T: Turbid.

*: OD out of range. NA: Not applicable. Samples are diluted to get an OD 585 in the range of 0.130 to 0.175 and OD 635 in the range of 0.120 to 0.160, to obtain a COD value in the range of 320 to 535 mg/L. (Based on the Tables 3 to 15 given in Examples 3 to 15).

EXAMPLE 50:

Dissolved 0.15 g glucose in 25 mL distilled water. Aliquots from the stock solution were taken and diluted with water to achieve different COD concentrations in the range of 320 to 535 mg/L. A 30 mL sample of Digested POME and Tea (as such) was taken. Weighed 0.02 g HgSO_4 in separate 15 mL test tubes. 1.0 mL of glucose

solution and sample were added to each test tube separately and mixed well. 0.5 mL of $K_2Cr_2O_7$ was added to each sample. 1.5 mL H_2SO_4 silver sulphate reagent was added slowly with continuous shaking. The color of the reaction mixture was noted down visually within 15 sec of incubation and also recorded spectrophotometrically at 585 and 635 nm against air. An approximate COD value of the sample was estimated from the color of standard reaction mixture prepared with glucose as reference material in the COD range of 200 to 100000 mg/L. OD 585 and OD 635 of the reaction mixture was read spectrophotometrically. OD 585 and 635 nm in the range of 0.130 to 0.175 and 0.120 to 0.160 respectively were considered and calculated with reference to glucose standard. (Glucose COD in the range of 320 to 535 mg/L corresponds to OD 585 and 635 nm in the range of 0.157 to 0.207 and 0.142 to 0.189, respectively).

COD of the sample was estimated to be 59,593 mg/L. (By Merck method and read on photometer).

Based on the color of the reaction mixture, COD of the sample was estimated to be around 50000 mg/L. On the basis of the dilution method and the spectrophotometric reading at OD 585 and OD 635 of the sample, COD was estimated to be 41600 and 41900 mg/l respectively. A comparison with the COD value calculated through Merck method reveals the following relation i.e. Our COD value (at OD 585) $\times 1.43 =$ COD value (Merck Method) or our COD value (at OD 635) $\times 1.42 =$ COD value (Merck Method).

Hence for a preliminary estimate of COD value color of the standard reaction mixture based on a range of glucose concentrations can be used as a reference. For samples, with COD concentration of 300 to 500 mg/L, OD 585 or OD 635 can be recorded for obtaining more precise values.



Digested Pome + Tea (as such)

Table 50: Correlation between COD concentration, color of the reaction mixture and their OD at different wave lengths.

COD (mg/L)	Sample dilution (Times)	Color of the reaction mixture	OD	Estimated COD (mg/L)	Estimated final COD (mg/L) (Col. 2x5)
1	2	3	4	5	6
<u>Values at OD 585 nm</u>					
Glucose:					
320	NA	C300	0.162	NA	NA
426	NA	C400	0.187	NA	NA
535	NA	C550	0.208	NA	NA
Sample: Digested Pome + Tea (as such)					
	Nil	C10000	-	>10,000	>10,000
	10	C5000	0.416	≈5000	≈50,000#
	100	C400	0.149	416	41600
<u>Values at OD 635 nm</u>					
Glucose:					
320	NA	C300	0.141	NA	NA
426	NA	C400	0.162	NA	NA
535	NA	C500	0.192	NA	NA
Sample: Digested Pome + Tea (as such)					
	Nil	C10000	-	>10,000	>10,000
	10	C5000	0.353	≈5000	≈50,000#
	100	C400	0.130	419	41900

#: COD as deducible from Table 2 given in Example 2. D: Dark. T: Turbid.

*: OD out of range. NA: Not applicable. Samples are diluted to get an OD 585 in the range of 0.130 to 0.175 and OD 635 in the range of 0.120 to 0.160, to obtain a COD value in the range of 320 to 535 mg/L. (Based on the Tables 3 to 15 given in Examples 3 to 15).

EXAMPLE 51:

Dissolved 0.15 g glucose in 25 mL distilled water. Aliquots from the stock solution were taken and diluted with water to achieve different COD concentrations in the range of 320 to 535 mg/L. A 30 mL sample of Digested POME and Tea (ssp) was taken. Weighed 0.02 g HgSO_4 in separate 15 mL test tubes. 1.0 mL of glucose solution and sample were added to each test tube separately and mixed well. 0.5 mL

of $K_2Cr_2O_7$ was added to each sample. 1.5 mL H_2SO_4 silver sulphate reagent was added slowly with continuous shaking. The color of the reaction mixture was noted down visually within 15 sec of incubation and also recorded spectrophotometrically at 585 and 635 nm against air. An approximate COD value of the sample was estimated from the color of standard reaction mixture prepared with glucose as reference material in the COD range of 200 to 100000 mg/L. OD 585 and OD 635 of the reaction mixture was read spectrophotometrically. OD 585 and 635 nm in the range of 0.130 to 0.175 and 0.120 to 0.160 respectively were considered and calculated with reference to glucose standard. (Glucose COD in the range of 320 to 535 mg/L corresponds to OD 585 and 635 nm in the range of 0.157 to 0.207 and 0.142 to 0.189, respectively).

COD of the sample was estimated to be 39,928 mg/L. (By Merck method and read on photometer).

Based on the color of the reaction mixture, COD of the sample was estimated to be around 20000 mg/L. On the basis of the dilution method and the spectrophotometric reading at OD 585 and OD 635 of the sample, COD was estimated to be 38556 and 36985 mg/l respectively. A comparison with the COD value calculated through Merck method reveals the following relation i.e. Our COD value (at OD 585) $\times 1.43 =$ COD value (Merck Method) or our COD value (at OD 635) $\times 1.42 =$ COD value (Merck Method).

Hence for a preliminary estimate of COD value color of the standard reaction mixture based on a range of glucose concentrations can be used as a reference. For samples, with COD concentration of 300 to 500 mg/L, OD 585 or OD 635 can be recorded for obtaining more precise values.



Digested Pome + Tea (ssp)

Table 51: Correlation between COD concentration, color of the reaction mixture and their OD at different wave lengths.

COD (mg/L)	Sample dilution (Times)	Color of the reaction mixture	OD	Estimated COD (mg/L)	Estimated final COD (mg/L) (Col. 2x5)
1	2	3	4	5	6
<u>Values at OD 585 nm</u>					
Glucose:					
320	NA	C300	0.162	NA	NA
426	NA	C400	0.187	NA	NA
535	NA	C550	0.208	NA	NA
Sample: Digested Pome + Tea (ssp)					
	Nil	C10000	-	>10,000	≅10,000
	10	C2000	0.386	≅2000	≅20,000#
	100	C200	0.107	200	20,000
	71.4	C500	0.175	540	38556
<u>Values at OD 635 nm</u>					
Glucose:					
320	NA	C300	0.141	NA	NA
426	NA	C400	0.162	NA	NA
535	NA	C550	0.192	NA	NA
Sample: Digested Pome + Tea (ssp)					
	Nil	C10000	-	>10,000	≅10,000
	10	C2000	0.336	≅2000	≅20,000#
	100	C200	0.092	200	20,000
	71.4	C500	0.155	518	36985

#: COD as deducible from Table 2 given in Example 2. D: Dark. T: Turbid.

*: OD out of range. NA: Not applicable. Samples are diluted to get an OD 585 in the range of 0.130 to 0.175 and OD 635 in the range of 0.120 to 0.160, to obtain a COD value in the range of 320 to 535 mg/L. (Based on the Tables 3 to 15 given in Examples 3 to 15).

EXAMPLE 52:

Dissolved 0.15 g glucose in 25 mL distilled water. Aliquots from the stock solution were taken and diluted with water to achieve different COD concentrations in the

range of 320 to 535 mg/L. A 30 mL sample of Digested POME and Apple pomace (as such) was taken. Weighed 0.02 g HgSO_4 in separate 15 mL test tubes. 1.0 mL of glucose solution and sample were added to each test tube separately and mixed well. 0.5 mL of $\text{K}_2\text{Cr}_2\text{O}_7$ was added to each sample. 1.5 mL H_2SO_4 silver sulphate reagent was added slowly with continuous shaking. The color of the reaction mixture was noted down visually within 15 sec of incubation and also recorded spectrophotometrically at 585 and 635 nm against air. An approximate COD value of the sample was estimated from the color of standard reaction mixture prepared with glucose as reference material in the COD range of 200 to 100000 mg/L. OD 585 and OD 635 of the reaction mixture was read spectrophotometrically. OD 585 and 635 nm in the range of 0.130 to 0.175 and 0.120 to 0.160 respectively were considered and calculated with reference to glucose standard. (Glucose COD in the range of 320 to 535 mg/L corresponds to OD 585 and 635 nm in the range of 0.157 to 0.207 and 0.142 to 0.189, respectively).

COD of the sample was estimated to be 24,932 mg/L. (By Merck method and read on photometer).

Based on the color of the reaction mixture, COD of the sample was estimated to be around 15000 mg/L. On the basis of the dilution method and the spectrophotometric reading at OD 585 and OD 635 of the sample, COD was estimated to be 15800 and 16955 mg/l respectively. A comparison with the COD value calculated through Merck method reveals the following relation i.e. Our COD value (at OD 585) $\times 1.58 =$ COD value (Merck Method) or our COD value (at OD 635) $\times 1.47 =$ COD value (Merck Method).

Hence for a preliminary estimate of COD value color of the standard reaction mixture based on a range of glucose concentrations can be used as a reference. For samples, with COD concentration of 300 to 500 mg/L, OD 585 or OD 635 can be recorded for obtaining more precise values.



Digested Pome + Apple Pomace (as such)

Table 52: Correlation between COD concentration, color of the reaction mixture and their OD at different wave lengths.

COD (mg/L)	Sample dilution (Times)	Color of the reaction mixture	OD	Estimated COD (mg/L)	Estimated final COD (mg/L) (Col. 2x5)
1	2	3	4	5	6
<u>Values at OD 585 nm</u>					
Glucose:					
320	NA	C300	0.158	NA	NA
426	NA	C400	0.186	NA	NA
535	NA	C550	0.205	NA	NA
Sample: Digested Pome + Apple Pomace (as such)					
	Nil	C10000	-	>10,000	>10,000
	10	C1500	0.336	≅1500	≅15,000#
	50	C400	0.130	316	15800
<u>Values at OD 635 nm</u>					
Glucose:					
320	NA	C300	0.137	NA	NA
426	NA	C400	0.159	NA	NA
535	NA	C550	0.191	NA	NA
Sample: Digested Pome + Apple Pomace (as such)					
	Nil	C10000	-	>10,000	>10,000
	10	C1500	0.305	≅1500	≅15,000#
	50	C400	0.121	339	16955

#: COD as deducible from Table 2 given in Example 2. D: Dark. T: Turbid.

*: OD out of range. NA: Not applicable. Samples are diluted to get an OD 585 in the range of 0.130 to 0.175 and OD 635 in the range of 0.120 to 0.160, to obtain a COD value in the range of 320 to 535 mg/L. (Based on the Tables 3 to 15 given in Examples 3 to 15).

EXAMPLE 53:

Dissolved 0.15 g glucose in 25 mL distilled water. Aliquots from the stock solution were taken and diluted with water to achieve different COD concentrations in the range of 320 to 535 mg/L. A 30 mL sample of Digested Palm oil mill effluent (POME)

and Apple pomace (ssp) was taken. Weighed 0.02 g HgSO_4 in separate 15 mL test tubes. 1.0 mL of glucose solution and sample were added to each test tube separately and mixed well. 0.5 mL of $\text{K}_2\text{Cr}_2\text{O}_7$ was added to each sample. 1.5 mL H_2SO_4 silver sulphate reagent was added slowly with continuous shaking. The color of the reaction mixture was noted down visually within 15 sec of incubation and also recorded spectrophotometrically at 585 and 635 nm against air. An approximate COD value of the sample was estimated from the color of standard reaction mixture prepared with glucose as reference material in the COD range of 200 to 100000 mg/L. OD 585 and OD 635 of the reaction mixture was read spectrophotometrically. OD 585 and 635 nm in the range of 0.130 to 0.175 and 0.120 to 0.160 respectively were considered and calculated with reference to glucose standard. (Glucose COD in the range of 320 to 535 mg/L corresponds to OD 585 and 635 nm in the range of 0.157 to 0.207 and 0.142 to 0.189, respectively).

COD of the sample was estimated to be 21,666 mg/L. (By Merck method and read on photometer).

Based on the color of the reaction mixture, COD of the sample was estimated to be around 16000 mg/L. On the basis of the dilution method and the spectrophotometric reading at OD 585 and OD 635 of the sample, COD was estimated to be 32600 and 41310 mg/l respectively. A comparison with the COD value calculated through Merck method reveals the following relation i.e. Our COD value (at OD 585) $\times 0.66 =$ COD value (Merck Method) or our COD value (at OD 635) $\times 0.52 =$ COD value (Merck Method).

Hence for a preliminary estimate of COD value color of the standard reaction mixture based on a range of glucose concentrations can be used as a reference. For samples, with COD concentration of 300 to 500 mg/L, OD 585 or OD 635 can be recorded for obtaining more precise values.



Digested Pome + Apple Pomace (ssp)

Table 53: Correlation between COD concentration, color of the reaction mixture and their OD at different wave lengths.

COD (mg/L)	Sample dilution (Times)	Color of the reaction mixture	OD	Estimated COD (mg/L)	Estimated final COD (mg/L) (Col. 2x5)
1	2	3	4	5	6
<u>Values at OD 585 nm</u>					
Glucose:					
320	NA	C300	0.158	NA	NA
426	NA	C400	0.186	NA	NA
535	NA	C550	0.205	NA	NA
Sample: Digested Pome + Apple Pomace (ssp)					
	Nil	C10000	-	>10,000	>10,000
	10	C1600	0.357	≅1600	≅16,000#
	100	C300	0.134	326	32600
<u>Values at OD 635 nm</u>					
Glucose:					
320	NA	C300	0.137	NA	NA
426	NA	C400	0.159	NA	NA
535	NA	C550	0.191	NA	NA
Sample: Digested Pome + Apple Pomace (ssp)					
	Nil	C10000	-	>10,000	>10,000
	10	C1600	0.357	≅1600	≅16,000#
	100	C300	0.134	413	41300

#: COD as deducible from Table 2 given in Example 2. D: Dark. T: Turbid.

*: OD out of range. NA: Not applicable. Samples are diluted to get an OD 585 in the range of 0.130 to 0.175 and OD 635 in the range of 0.120 to 0.160, to obtain a COD value in the range of 320 to 535 mg/L. (Based on the Tables 3 to 15 given in Examples 3 to 15).

EXAMPLE 54:

Dissolved 0.15 g glucose in 25 mL distilled water. Aliquots from the stock solution were taken and diluted with water to achieve different COD concentrations in the range of 320 to 535 mg/L. A 30 mL sample of Digested POME, Apple pomace and

Tea (as such) was taken. Weighed 0.02 g HgSO_4 in separate 15 mL test tubes. 1.0 mL of glucose solution and sample were added to each test tube separately and mixed well. 0.5 mL of $\text{K}_2\text{Cr}_2\text{O}_7$ was added to each sample. 1.5 mL H_2SO_4 silver sulphate reagent was added slowly with continuous shaking. The color of the reaction mixture was noted down visually within 15 sec of incubation and also recorded spectrophotometrically at 585 and 635 nm against air. An approximate COD value of the sample was estimated from the color of standard reaction mixture prepared with glucose as reference material in the COD range of 200 to 100000 mg/L. OD 585 and OD 635 of the reaction mixture was read spectrophotometrically. OD 585 and 635 nm in the range of 0.130 to 0.175 and 0.120 to 0.160 respectively were considered and calculated with reference to glucose standard. (Glucose COD in the range of 320 to 535 mg/L corresponds to OD 585 and 635 nm in the range of 0.157 to 0.207 and 0.142 to 0.189, respectively).

COD of the sample was estimated to be 34,569 mg/L. (By Merck method and read on photometer).

Based on the color of the reaction mixture, COD of the sample was estimated to be around 45000 mg/L. On the basis of the dilution method and the spectrophotometric reading at OD 585 and OD 635 of the sample, COD was estimated to be 35774 and 36552 mg/l respectively. A comparison with the COD value calculated through Merck method reveals the following relation i.e. Our COD value (at OD 585) $\times 0.97 =$ COD value (Merck Method) or our COD value (at OD 635) $\times 0.94 =$ COD value (Merck Method).

Hence for a preliminary estimate of COD value color of the standard reaction mixture based on a range of glucose concentrations can be used as a reference. For samples, with COD concentration of 300 to 500 mg/L, OD 585 or OD 635 can be recorded for obtaining more precise values.



Digested Pome + Apple Pomace + Tea (as such)

Table 54: Correlation between COD concentration, color of the reaction mixture and their OD at different wave lengths.

COD (mg/L)	Sample dilution (Times)	Color of the reaction mixture	OD	Estimated COD (mg/L)	Estimated final COD (mg/L) (Col. 2x5)
1	2	3	4	5	6
<u>Values at OD 585 nm</u>					
Glucose:					
320	NA	C300	0.160	NA	NA
426	NA	C400	0.185	NA	NA
535	NA	C550	0.211	NA	NA
Sample: Digested Pome + Apple Pomace + Tea (as such)					
	Nil	C10000	-	>10,000	>10,000
	10	C4500	0.428	≈4500	≈45,000#
	100	C450	0.189	≈450	≈45,000#
	111.1	C300	0.134	322	35774
<u>Values at OD 635 nm</u>					
Glucose:					
320	NA	C300	0.141	NA	NA
426	NA	C400	0.162	NA	NA
535	NA	C550	0.197	NA	NA
Sample: Digested Pome + Apple Pomace + Tea (as such)					
	Nil	C10000	-	>10,000	>10,000
	10	C4500	0.375	≈4500	≈45,000#
	100	C450	0.173	≈450	≈45,000#
	111.1	C300	0.121	329	36552

#: COD as deducible from Table 2 given in Example 2. D: Dark. T: Turbid.

*: OD out of range. NA: Not applicable. Samples are diluted to get an OD 585 in the range of 0.130 to 0.175 and OD 635 in the range of 0.120 to 0.160, to obtain a COD value in the range of 320 to 535 mg/L. (Based on the Tables 3 to 15 given in Examples 3 to 15).

EXAMPLE 55:

Dissolved 0.15 g glucose in 25 mL distilled water. Aliquots from the stock solution were taken and diluted with water to achieve different COD concentrations in the range of 320 to 535 mg/L. A 30 mL sample of Digested POME, Apple pomace and

Tea (ssp) was taken. Weighed 0.02 g HgSO_4 in separate 15 mL test tubes. 1.0 mL of glucose solution and sample were added to each test tube separately and mixed well. 0.5 mL of $\text{K}_2\text{Cr}_2\text{O}_7$ was added to each sample. 1.5 mL H_2SO_4 silver sulphate reagent was added slowly with continuous shaking. The color of the reaction mixture was noted down visually within 15 sec of incubation and also recorded spectrophotometrically at 585 and 635 nm against air. An approximate COD value of the sample was estimated from the color of standard reaction mixture prepared with glucose as reference material in the COD range of 200 to 100000 mg/L. OD 585 and OD 635 of the reaction mixture was read spectrophotometrically. OD 585 and 635 nm in the range of 0.130 to 0.175 and 0.120 to 0.160 respectively were considered and calculated with reference to glucose standard. (Glucose COD in the range of 320 to 535 mg/L corresponds to OD 585 and 635 nm in the range of 0.157 to 0.207 and 0.142 to 0.189, respectively).

COD of the sample was estimated to be 34,040 mg/L. (By Merck method and read on photometer).

Based on the color of the reaction mixture, COD of the sample was estimated to be around 40000 mg/L. On the basis of the dilution method and the spectrophotometric reading at OD 585 and OD 635 of the sample, COD was estimated to be 46800 and 48100 mg/l. A comparison with the COD value calculated through Merck method reveals the following relation i.e. Our COD value (at OD 585) $\times 0.73 =$ COD value (Merck Method) or our COD value (at OD 635) $\times 0.71 =$ COD value (Merck Method). Hence for a preliminary estimate of COD value color of the standard reaction mixture based on a range of glucose concentrations can be used as a reference. For samples, with COD concentration of 300 to 500 mg/L, OD 585 or OD 635 can be recorded for obtaining more precise values.



Digested Pome + Apple Pomace + Tea (ssp)

Table 55: Correlation between COD concentration, color of the reaction mixture and their OD at different wave lengths.

COD (mg/L)	Sample dilution (Times)	Color of the reaction mixture	OD	Estimated COD (mg/L)	Estimated final COD (mg/L) (Col. 2x5)
1	2	3	4	5	6
<u>Values at OD 585 nm</u>					
Glucose:					
320	NA	C300	0.160	NA	NA
426	NA	C400	0.185	NA	NA
535	NA	C550	0.211	NA	NA
Sample: Digested Pome + Apple Pomace + Tea (ssp)					
	Nil	C10000	-	>10,000	>10,000
	10	C4000	0.373	≈4000	≈40,000#
	100	C400	0.166	468	46800
<u>Values at OD 635 nm</u>					
Glucose:					
320	NA	C300	0.141	NA	NA
426	NA	C400	0.162	NA	NA
535	NA	C550	0.197	NA	NA
Sample: Digested Pome + Apple Pomace + Tea (ssp)					
	Nil	C10000	-	>10,000	>10,000
	10	C4000	0.322	≈4000	≈40,000#
	100	C400	0.149	481	48100

#: COD as deducible from Table 2 given in Example 2. D: Dark. T: Turbid.

*: OD out of range. NA: Not applicable. Samples are diluted to get an OD 585 in the range of 0.130 to 0.175 and OD 635 in the range of 0.120 to 0.160, to obtain a COD value in the range of 320 to 535 mg/L. (Based on the Tables 3 to 15 given in Examples 3 to 15).

EXAMPLE 56:

Dissolved 0.15 g glucose in 25 mL distilled water. Aliquots from the stock solution were taken and diluted with water to achieve different COD concentrations in the

range of 320 to 535 mg/L. A 30 mL sample of Tea and Apple pomace (ssp) was taken. Weighed 0.02 g HgSO_4 in separate 15 mL test tubes. 1.0 mL of glucose solution and sample were added to each test tube separately and mixed well. 0.5 mL of $\text{K}_2\text{Cr}_2\text{O}_7$ was added to each sample. 1.5 mL H_2SO_4 silver sulphate reagent was added slowly with continuous shaking. The color of the reaction mixture was noted down visually within 15 sec of incubation and also recorded spectrophotometrically at 585 and 635 nm against air. An approximate COD value of the sample was estimated from the color of standard reaction mixture prepared with glucose as reference material in the COD range of 200 to 100000 mg/L. OD 585 and OD 635 of the reaction mixture was read spectrophotometrically. OD 585 and 635 nm in the range of 0.130 to 0.175 and 0.120 to 0.160 respectively were considered and calculated with reference to glucose standard. (Glucose COD in the range of 320 to 535 mg/L corresponds to OD 585 and 635 nm in the range of 0.157 to 0.207 and 0.142 to 0.189, respectively).

COD of the sample was estimated to be 75,716 mg/L. (By Merck method and read on photometer).

Based on the color of the reaction mixture, COD of the sample was estimated to be around 78540 mg/L. On the basis of the dilution method and the spectrophotometric reading at OD 585 and OD 635 of the sample, COD was estimated to be 86400 and 83600 mg/l respectively. A comparison with the COD value calculated through Merck method reveals the following relation i.e. Our COD value (at OD 585) $\times 0.88 =$ COD value (Merck Method) or our COD value (at OD 635) $\times 0.91 =$ COD value (Merck Method).

Hence for a preliminary estimate of COD value color of the standard reaction mixture based on a range of glucose concentrations can be used as a reference. For samples, with COD concentration of 300 to 500 mg/L, OD 585 or OD 635 can be recorded for obtaining more precise values.

Tea + Apple Pomace (ssp)

Table 56: Correlation between COD concentration, color of the reaction mixture and their OD at different wave lengths.

COD (mg/L)	Sample dilution (Times)	Color of the reaction mixture	OD	Estimated COD (mg/L)	Estimated final COD (mg/L) (Col. 2x5)
1	2	3	4	5	6
<u>Values at OD 585 nm</u>					
Glucose:					
320	NA	C300	0.159	NA	NA
426	NA	C400	0.186	NA	NA
535	NA	C550	0.205	NA	NA
Sample: Tea + Apple Pomace (ssp)					
	Nil	C10000	-	>10,000	≅>10,000
	10	C5500	0.372	≅5500	≅ 55,000
	100	C1600	0.225	≅1600	≅1,60,000
	142.8	C550	0.175	≅550	≅78,540#
	200	C400	0.157	432	86400
<u>Values at OD 635 nm</u>					
Glucose:					
320	NA	C300	0.143	NA	NA
426	NA	C400	0.168	NA	NA
535	NA	C550	0.186	NA	NA
Sample: Tea + Apple Pomace (ssp)					
	Nil	C10000	-	>10,000	≅>10,000
	10	C5500	0.295	≅5500	≅ 55,000
	100	C1600	0.201	≅1600	≅1,60,000
	142.8	C550	0.155	≅550	≅78,540#
	200	C400	0.136	418	83600

#: COD as deducible from Table 2 given in Example 2. D: Dark. T: Turbid.

*: OD out of range. NA: Not applicable. Samples are diluted to get an OD 585 in the range of 0.130 to 0.175 and OD 635 in the range of 0.120 to 0.160, to obtain a COD value in the range of 320 to 535 mg/L. (Based on the Tables 3 to 15 given in Examples 3 to 15).

EXAMPLE 57:

Dissolved 0.15 g glucose in 25 mL distilled water. Aliquots from the stock solution were taken and diluted with water to achieve different COD concentrations in the range of 320 to 535 mg/L. A 30 mL sample of Damaged wheat grain (ssp) was taken. Weighed 0.02 g HgSO_4 in separate 15 mL test tubes. 1.0 mL of glucose solution and sample were added to each test tube separately and mixed well. 0.5 mL of $\text{K}_2\text{Cr}_2\text{O}_7$ was added to each sample. 1.5 mL H_2SO_4 silver sulphate reagent was added slowly with continuous shaking. The color of the reaction mixture was noted down visually within 15 sec of incubation and also recorded spectrophotometrically at 585 and 635 nm against air. An approximate COD value of the sample was estimated from the color of standard reaction mixture prepared with glucose as reference material in the COD range of 200 to 100000 mg/L. OD 585 and OD 635 of the reaction mixture was read spectrophotometrically. OD 585 and 635 nm in the range of 0.130 to 0.175 and 0.120 to 0.160 respectively were considered and calculated with reference to glucose standard. (Glucose COD in the range of 320 to 535 mg/L corresponds to OD 585 and 635 nm in the range of 0.157 to 0.207 and 0.142 to 0.189, respectively).

COD of the sample was estimated to be 11,707 mg/L. (By Merck method and read on photometer).

Based on the color of the reaction mixture, COD of the sample was estimated to be around 1400 mg/L. On the basis of the dilution method and the spectrophotometric reading at OD 585 and OD 635 of the sample, COD was estimated to be 5160 and 5040 mg/l respectively. A comparison with the COD value calculated through Merck method reveals the following relation i.e. Our COD value (at OD 585) $\times 2.26 =$ COD value (Merck Method) or our COD value (at OD 635) $\times 2.31 =$ COD value (Merck Method).

Hence for a preliminary estimate of COD value color of the standard reaction mixture based on a range of glucose concentrations can be used as a reference. For samples, with COD concentration of 300 to 500 mg/L, OD 585 or OD 635 can be recorded for obtaining more precise values.

DWS (S3) (ssp)

Table 57: Correlation between COD concentration, color of the reaction mixture and their OD at different wave lengths.

COD (mg/L)	Sample dilution (Times)	Color of the reaction mixture	OD	Estimated COD (mg/L)	Estimated final COD (mg/L) (Col. 2x5)
1	2	3	4	5	6
<u>Values at OD 585 nm</u>					
Glucose:					
320	NA	C300	0.159	NA	NA
426	NA	C400	0.186	NA	NA
535	NA	C550	0.205	NA	NA
Sample: DWS (S3) (ssp)					
	Nil	C1400	0.636	≈1400	≈1400#
	10	C400	0.165	516	5160
<u>Values at OD 635 nm</u>					
Glucose:					
320	NA	C300	0.143	NA	NA
426	NA	C400	0.168	NA	NA
535	NA	C550	0.186	NA	NA
Sample: DWS (S3) (ssp)					
	Nil	C1400	0.494	≈1400	≈1400#
	10	C400	0.147	504	5040

#: COD as deducible from Table 2 given in Example 2. D: Dark. T: Turbid.

*: OD out of range. NA: Not applicable. Samples are diluted to get an OD 585 in the range of 0.130 to 0.175 and OD 635 in the range of 0.120 to 0.160, to obtain a COD value in the range of 320 to 535 mg/L. (Based on the Tables 3 to 15 given in Examples 3 to 15).

EXAMPLE 58:

Dissolved 0.15 g glucose in 25 mL distilled water. Aliquots from the stock solution were taken and diluted with water to achieve different COD concentrations in the range of 320 to 535 mg/L. A 30 mL sample of Tea, Apple pomace and Damaged

wheat grain (ssp) was taken. Weighed 0.02 g HgSO_4 in separate 15 mL test tubes. 1.0 mL of glucose solution and sample were added to each test tube separately and mixed well. 0.5 mL of $\text{K}_2\text{Cr}_2\text{O}_7$ was added to each sample. 1.5 mL H_2SO_4 silver sulphate reagent was added slowly with continuous shaking. The color of the reaction mixture was noted down visually within 15 sec of incubation and also recorded spectrophotometrically at 585 and 635 nm against air. An approximate COD value of the sample was estimated from the color of standard reaction mixture prepared with glucose as reference material in the COD range of 200 to 100000 mg/L. OD 585 and OD 635 of the reaction mixture was read spectrophotometrically. OD 585 and 635 nm in the range of 0.130 to 0.175 and 0.120 to 0.160 respectively were considered and calculated with reference to glucose standard. (Glucose COD in the range of 320 to 535 mg/L corresponds to OD 585 and 635 nm in the range of 0.157 to 0.207 and 0.142 to 0.189, respectively).

COD of the sample was estimated to be 22,689.5 mg/L. (By Merck method and read on photometer).

Based on the color of the reaction mixture, COD of the sample was estimated to be around 20000 mg/L. On the basis of the dilution method and the spectrophotometric reading at OD 585 and OD 635 of the sample, COD was estimated to be 22450 and 16950 mg/l respectively. A comparison with the COD value calculated through Merck method reveals the following relation i.e. Our COD value (at OD 585) $\times 2.26 =$ COD value (Merck Method) or our COD value (at OD 635) $\times 2.31 =$ COD value (Merck Method).

Hence for a preliminary estimate of COD value color of the standard reaction mixture based on a range of glucose concentrations can be used as a reference. For samples, with COD concentration of 300 to 500 mg/L, OD 585 or OD 635 can be recorded for obtaining more precise values.



Tea + Apple Pomace + DWS (S3) (ssp)

Table 58: Correlation between COD concentration, color of the reaction mixture and their OD at different wave lengths.

COD (mg/L)	Sample dilution (Times)	Color of the reaction mixture	OD	Estimated COD (mg/L)	Estimated final COD (mg/L) (Col. 2x5)
1	2	3	4	5	6
<u>Values at OD 585 nm</u>					
Glucose:					
320	NA	C300	0.161	NA	NA
426	NA	C400	0.185	NA	NA
535	NA	C550	0.209	NA	NA
Sample: Tea + Apple Pomace + DWS (S3) (ssp)					
	Nil	C10000	-	>10,000	>10,000
	10	C4000	0.356	≈4000	≈40,000
	66.7	C300	0.127	≈300	≈20,000#
	50	C400	0.162	449	22450
<u>Values at OD 635 nm</u>					
Glucose:					
320	NA	C300	0.142	NA	NA
426	NA	C400	0.163	NA	NA
535	NA	C550	0.198	NA	NA
Sample: Tea + Apple Pomace + DWS (S3) (ssp)					
	Nil	C10000	-	>10,000	>10,000
	10	C4000	0.319	≈4000	≈40,000
	66.7	C300	0.113	≈300	≈20,000#
	50	C400	0.142	339	16950

#: COD as deducible from Table 2 given in Example 2. D: Dark. T: Turbid.

*: OD out of range. NA: Not applicable. Samples are diluted to get an OD 585 in the range of 0.130 to 0.175 and OD 635 in the range of 0.120 to 0.160, to obtain a COD value in the range of 320 to 535 mg/L. (Based on the Tables 3 to 15 given in Examples 3 to 15).

EXAMPLE 59:

Dissolved 0.15 g glucose in 25 mL distilled water. Aliquots from the stock solution were taken and diluted with water to achieve different COD concentrations in the range of 320 to 535 mg/L. A 30 mL sample of Tea, Apple pomace, Damaged wheat grain and Peptone (ssp) was taken. Weighed 0.02 g HgSO_4 in separate 15 mL test tubes. 1.0 mL of glucose solution and sample were added to each test tube separately and mixed well. 0.5 mL of $\text{K}_2\text{Cr}_2\text{O}_7$ was added to each sample. 1.5 mL H_2SO_4 silver sulphate reagent was added slowly with continuous shaking. The color of the reaction mixture was noted down visually within 15 sec of incubation and also recorded spectrophotometrically at 585 and 635 nm against air. An approximate COD value of the sample was estimated from the color of standard reaction mixture prepared with glucose as reference material in the COD range of 200 to 100000 mg/L. OD 585 and OD 635 of the reaction mixture was read spectrophotometrically. OD 585 and 635 nm in the range of 0.130 to 0.175 and 0.120 to 0.160 respectively were considered and calculated with reference to glucose standard. (Glucose COD in the range of 320 to 535 mg/L corresponds to OD 585 and 635 nm in the range of 0.157 to 0.207 and 0.142 to 0.189, respectively).

COD of the sample was estimated to be 36,176 mg/L. (By Merck method and read on photometer).

Based on the color of the reaction mixture, COD of the sample was estimated to be around 20000 mg/L. On the basis of the dilution method and the spectrophotometric reading at OD 585 and OD 635 of the sample, COD was estimated to be 20900 and 21050 mg/l respectively. A comparison with the COD value calculated through Merck method reveals the following relation i.e. Our COD value (at OD 585) $\times 1.73 =$ COD value (Merck Method) or our COD value (at OD 635) $\times 1.72 =$ COD value (Merck Method).

Hence for a preliminary estimate of COD value color of the standard reaction mixture based on a range of glucose concentrations can be used as a reference. For samples, with COD concentration of 300 to 500 mg/L, OD 585 or OD 635 can be recorded for obtaining more precise values.



Tea + Apple Pomace + DWS (S3) + Peptone (ssp)

Table 59: Correlation between COD concentration, color of the reaction mixture and their OD at different wave lengths.

COD (mg/L)	Sample dilution (Times)	Color of the reaction mixture	OD	Estimated COD (mg/L)	Estimated final COD (mg/L) (Col. 2x5)
1	2	3	4	5	6

Values at OD 585 nm

Glucose:

320	NA	C300	0.161	NA	NA
426	NA	C400	0.185	NA	NA
535	NA	C550	0.209	NA	NA

Sample: Tea + Apple Pomace + DWS (S3) + Peptone (ssp)

Nil	C10000	-	>10,000	>10,000
10	C2000	0.348	≈2000	≈20,000#
50	C300	0.151	418	20900

Values at OD 635 nm

Glucose:

320	NA	C300	0.142	NA	NA
426	NA	C400	0.163	NA	NA
535	NA	C550	0.198	NA	NA

Sample: Tea + Apple Pomace + DWS (S3) + Peptone(ssp)

Nil	C10000	-	>10,000	>10,000
10	C20000	0.315	≈2000	≈20,000#
50	C300	0.134	421	21050

#: COD as deducible from Table 2 given in Example 2. D: Dark. T: Turbid.

*: OD out of range. NA: Not applicable. Samples are diluted to get an OD 585 in the range of 0.130 to 0.175 and OD 635 in the range of 0.120 to 0.160, to obtain a COD value in the range of 320 to 535 mg/L. (Based on the Tables 3 to 15 given in Examples 3 to 15).

EXAMPLE 60:

Dissolved 0.15 g glucose in 25 mL distilled water. Aliquots from the stock solution were taken and diluted with water to achieve different COD concentrations in the range of 320 to 535 mg/L. A 30 mL sample of Tea, Apple Pomace and Dabu (ssp) was taken. Weighed 0.02 g HgSO₄ in separate 15 mL test tubes. 1.0 mL of glucose

solution and sample were added to each test tube separately and mixed well. 0.5 mL of $K_2Cr_2O_7$ was added to each sample. 1.5 mL H_2SO_4 silver sulphate reagent was added slowly with continuous shaking. The color of the reaction mixture was noted down visually within 15 sec of incubation and also recorded spectrophotometrically at 585 and 635 nm against air. An approximate COD value of the sample was estimated from the color of standard reaction mixture prepared with glucose as reference material in the COD range of 200 to 100000 mg/L. OD 585 and OD 635 of the reaction mixture was read spectrophotometrically. OD 585 and 635 nm in the range of 0.130 to 0.175 and 0.120 to 0.160 respectively were considered and calculated with reference to glucose standard. (Glucose COD in the range of 320 to 535 mg/L corresponds to OD 585 and 635 nm in the range of 0.157 to 0.207 and 0.142 to 0.189, respectively).

COD of the sample was estimated to be 77,786 mg/L. (By Merck method and read on photometer).

Based on the color of the reaction mixture, COD of the sample was estimated to be around 55000 mg/L. On the basis of the dilution method and the spectrophotometric reading at OD 585 and OD 635 of the sample, COD was estimated to be 71714 and 74348 mg/l respectively. A comparison with the COD value calculated through Merck method reveals the following relation i.e. Our COD value (at OD 585) $\times 1.08$ = COD value (Merck Method) or our COD value (at OD 635) $\times 1.05$ = COD value (Merck Method).

Hence for a preliminary estimate of COD value color of the standard reaction mixture based on a range of glucose concentrations can be used as a reference. For samples, with COD concentration of 300 to 500 mg/L, OD 585 or OD 635 can be recorded for obtaining more precise values.



Tea + Apple Pomace + Dabu (ssp)

Table 60: Correlation between COD concentration, color of the reaction mixture and their OD at different wave lengths.

COD (mg/L)	Sample dilution (Times)	Color of the reaction mixture	OD	Estimated COD (mg/L)	Estimated final COD (mg/L) (Col. 2x5)
1	2	3	4	5	6

Values at OD 585 nm

Glucose:

320	NA	C300	0.159	NA	NA
426	NA	C400	0.187	NA	NA
535	NA	C550	0.206	NA	NA

Sample: Tea + Apple Pomace + Dabu (ssp)

Nil	C10000	-	>10,000	>10,000
10	C5000	-	≈5000	≈50,000
100	C550	0.216	≈550	≈55,000#
166.7	C400	0.157	430	71681

Values at OD 635 nm

Glucose:

320	NA	C300	0.141	NA	NA
426	NA	C400	0.162	NA	NA
535	NA	C550	0.195	NA	NA

Sample: Tea + Apple Pomace + Dabu (ssp)

Nil	C10000	-	>10,000	>10,000
10	C5000	-	≈5000	≈50,000
100	C550	0.195	≈550	≈55,000#
166.7	C400	0.141	446	74348

#: COD as deducible from Table 2 given in Example 2. D: Dark. T: Turbid.

*: OD out of range. NA: Not applicable. Samples are diluted to get an OD 585 in the range of 0.130 to 0.175 and OD 635 in the range of 0.120 to 0.160, to obtain a COD value in the range of 320 to 535 mg/L. (Based on the Tables 3 to 15 given in Examples 3 to 15).

EXAMPLE 61:

Dissolved 0.15 g glucose in 25 mL distilled water. Aliquots from the stock solution were taken and diluted with water to achieve different COD concentrations in the

range of 320 to 535 mg/L. A 30 mL sample of Damaged wheat grain and Dabu (ssp) was taken. Weighed 0.02 g HgSO_4 in separate 15 mL test tubes. 1.0 mL of glucose solution and sample were added to each test tube separately and mixed well. 0.5 mL of $\text{K}_2\text{Cr}_2\text{O}_7$ was added to each sample. 1.5 mL H_2SO_4 silver sulphate reagent was added slowly with continuous shaking. The color of the reaction mixture was noted down visually within 15 sec of incubation and also recorded spectrophotometrically at 585 and 635 nm against air. An approximate COD value of the sample was estimated from the color of standard reaction mixture prepared with glucose as reference material in the COD range of 200 to 100000 mg/L. OD 585 and OD 635 of the reaction mixture was read spectrophotometrically. OD 585 and 635 nm in the range of 0.130 to 0.175 and 0.120 to 0.160 respectively were considered and calculated with reference to glucose standard. (Glucose COD in the range of 320 to 535 mg/L corresponds to OD 585 and 635 nm in the range of 0.157 to 0.207 and 0.142 to 0.189, respectively).

COD of the sample was estimated to be 65,251 mg/L. (By Merck method and read on photometer).

Based on the color of the reaction mixture, COD of the sample was estimated to be around 45000 mg/L. On the basis of the dilution method and the spectrophotometric reading at OD 585 and OD 635 of the sample, COD was estimated to be 63625 and 58000 mg/l respectively. A comparison with the COD value calculated through Merck method reveals the following relation i.e. Our COD value (at OD 585) $\times 1.02 =$ COD value (Merck Method) or our COD value (at OD 635) $\times 1.12 =$ COD value (Merck Method).

Hence for a preliminary estimate of COD value color of the standard reaction mixture based on a range of glucose concentrations can be used as a reference. For samples, with COD concentration of 300 to 500 mg/L, OD 585 or OD 635 can be recorded for obtaining more precise values.

DWS(S3) + Dabu (ssp)

Table 61: Correlation between COD concentration, color of the reaction mixture and their OD at different wave lengths.

COD (mg/L)	Sample dilution (Times)	Color of the reaction mixture	OD	Estimated COD (mg/L)	Estimated final COD (mg/L) (Col. 2x5)
1	2	3	4	5	6

Values at OD 585 nm

Glucose:

320	NA	C300	0.159	NA	NA
426	NA	C400	0.187	NA	NA
535	NA	C550	0.206	NA	NA

Sample: DWS(S3) + Dabu (ssp)

Nil	C5500	-	5500	≈5500
10	C4000	-	4000	≈40,000
100	C450	0.168	450	≈45,000#
125	C300	0.163	509	63625

Values at OD 635 nm

Glucose:

320	NA	C300	0.141	NA	NA
426	NA	C400	0.162	NA	NA
535	NA	C550	0.195	NA	NA

Sample: DWS(S3) + Dabu (ssp)

Nil	C5500	-	5500	≈5500
10	C4000	-	4000	≈40,000
100	C450	0.150	450	≈45,000#
125	C300	0.145	464	58000

#: COD as deducible from Table 2 given in Example 2. D: Dark. T: Turbid.

*: OD out of range. NA: Not applicable. Samples are diluted to get an OD 585 in the range of 0.130 to 0.175 and OD 635 in the range of 0.120 to 0.160, to obtain a COD value in the range of 320 to 535 mg/L. (Based on the Tables 3 to 15 given in Examples 3 to 15).

EXAMPLE 62:

Dissolved 0.15 g glucose in 25 mL distilled water. Aliquots from the stock solution were taken and diluted with water to achieve different COD concentrations in the range of 320 to 535 mg/L. A 30 mL sample of Tea, Apple Pomace, Damaged wheat grain and NaCl (ssp) was taken. Weighed 0.02 g HgSO_4 in separate 15 mL test tubes. 1.0 mL of glucose solution and sample were added to each test tube separately and mixed well. 0.5 mL of $\text{K}_2\text{Cr}_2\text{O}_7$ was added to each sample. 1.5 mL H_2SO_4 silver sulphate reagent was added slowly with continuous shaking. The color of the reaction mixture was noted down visually within 15 sec of incubation and also recorded spectrophotometrically at 585 and 635 nm against air. An approximate COD value of the sample was estimated from the color of standard reaction mixture prepared with glucose as reference material in the COD range of 200 to 100000 mg/L. OD 585 and OD 635 of the reaction mixture was read spectrophotometrically. OD 585 and 635 nm in the range of 0.130 to 0.175 and 0.120 to 0.160 respectively were considered and calculated with reference to glucose standard. (Glucose COD in the range of 320 to 535 mg/L corresponds to OD 585 and 635 nm in the range of 0.157 to 0.207 and 0.142 to 0.189, respectively).

COD of the sample was estimated to be 27,163 mg/L. (By Merck method and read on photometer).

Based on the color of the reaction mixture, COD of the sample was estimated to be around 10000 mg/L. On the basis of the dilution method and the spectrophotometric reading at OD 585 and OD 635 of the sample, COD was estimated to be 16300 and 17450 mg/l respectively. A comparison with the COD value calculated through Merck method reveals the following relation i.e. Our COD value (at OD 585) $\times 1.67 =$ COD value (Merck Method) or our COD value (at OD 635) $\times 1.56 =$ COD value (Merck Method).

Hence for a preliminary estimate of COD value color of the standard reaction mixture based on a range of glucose concentrations can be used as a reference. For samples, with COD concentration of 300 to 500 mg/L, OD 585 or OD 635 can be recorded for obtaining more precise values.

Table 62: Correlation between COD concentration, color of the reaction mixture and their OD at different wave lengths.

COD (mg/L)	Sample dilution (Times)	Color of the reaction mixture	OD	Estimated COD (mg/L)	Estimated final COD (mg/L) (Col. 2x5)
1	2	3	4	5	6
<u>Values at OD 585 nm</u>					
Glucose:					
320	NA	C300	0.158	NA	NA
426	NA	C400	0.185	NA	NA
535	NA	C550	0.207	NA	NA
Sample: Tea + Apple Pomace + DWS (S3) + NaCl (ssp)					
	Nil	C10000	-	≈10,000	≈10,000#
	50	C300	0.134	326	16300
	66.7	C300	0.112	300	20,010
<u>Values at OD 635 nm</u>					
Glucose:					
320	NA	C300	0.132	NA	NA
426	NA	C400	0.160	NA	NA
535	NA	C550	0.195	NA	NA
Sample: Tea + Apple Pomace + DWS (S3) + NaCl (ssp)					
	Nil	C10000	-	≈10,000	≈10,000#
	50	C300	0.120	349	17450
	66.7	C300	0.099	300	20,010

#: COD as deducible from Table 2 given in Example 2. D: Dark. T: Turbid.

*: OD out of range. NA: Not applicable. Samples are diluted to get an OD 585 in the range of 0.130 to 0.175 and OD 635 in the range of 0.120 to 0.160, to obtain a COD value in the range of 320 to 535 mg/L. (Based on the Tables 3 to 15 given in Examples 3 to 15).

EXAMPLE 63:

Dissolved 0.15 g glucose in 25 mL distilled water. Aliquots from the stock solution were taken and diluted with water to achieve different COD concentrations in the range of 320 to 535 mg/L. A 30 mL sample of Tea, Apple Pomace, Damaged wheat

grain and Glucose (ssp) was taken. Weighed 0.02 g HgSO_4 in separate 15 mL test tubes. 1.0 mL of glucose solution and sample were added to each test tube separately and mixed well. 0.5 mL of $\text{K}_2\text{Cr}_2\text{O}_7$ was added to each sample. 1.5 mL H_2SO_4 silver sulphate reagent was added slowly with continuous shaking. The color of the reaction mixture was noted down visually within 15 sec of incubation and also recorded spectrophotometrically at 585 and 635 nm against air. An approximate COD value of the sample was estimated from the color of standard reaction mixture prepared with glucose as reference material in the COD range of 200 to 100000 mg/L. OD 585 and OD 635 of the reaction mixture was read spectrophotometrically. OD 585 and 635 nm in the range of 0.130 to 0.175 and 0.120 to 0.160 respectively were considered and calculated with reference to glucose standard. (Glucose COD in the range of 320 to 535 mg/L corresponds to OD 585 and 635 nm in the range of 0.157 to 0.207 and 0.142 to 0.189, respectively).

COD of the sample was estimated to be 53,866 mg/L. (By Merck method and read on photometer).

Based on the color of the reaction mixture, COD of the sample was estimated to be around 50000 mg/L. On the basis of the dilution method and the spectrophotometric reading at OD 585 and OD 635 of the sample, COD was estimated to be 517000 and 41600 mg/l respectively. A comparison with the COD value calculated through Merck method reveals the following relation i.e. Our COD value (at OD 585) $\times 1.04 =$ COD value (Merck Method) or our COD value (at OD 635) $\times 1.30 =$ COD value (Merck Method).

Hence for a preliminary estimate of COD value color of the standard reaction mixture based on a range of glucose concentrations can be used as a reference. For samples, with COD concentration of 300 to 500 mg/L, OD 585 or OD 635 can be recorded for obtaining more precise values.

Tea + Apple Pomace + DWS (S3) + Glucose (ssp)

Table 63: Correlation between COD concentration, color of the reaction mixture and their OD at different wave lengths.

COD (mg/L)	Sample dilution (Times)	Color of the reaction mixture	OD	Estimated COD (mg/L)	Estimated final COD (mg/L) (Col. 2x5)
1	2	3	4	5	6

Values at OD 585 nm

Glucose:

320	NA	C300	0.158	NA	NA
426	NA	C400	0.185	NA	NA
535	NA	C550	0.207	NA	NA

Sample: Tea + Apple Pomace + DWS (S3) + Glucose (ssp)

Nil	C10000	-	>10,000	>10000
10	C5000	-	5000	50,000#
100	C500	0.167	517	51700

Values at OD 635 nm

Glucose:

320	NA	C300	0.132	NA	NA
426	NA	C400	0.160	NA	NA
535	NA	C550	0.195	NA	NA

Sample: Tea + Apple Pomace + DWS (S3) + Glucose (ssp)

Nil	C10000	-	>10,000	>10000
10	C5000	-	5000	50,000#
100	C400	0.150	416	41600

#: COD as deducible from Table 2 given in Example 2. D: Dark. T: Turbid.

*: OD out of range. NA: Not applicable. Samples are diluted to get an OD 585 in the range of 0.130 to 0.175 and OD 635 in the range of 0.120 to 0.160, to obtain a COD value in the range of 320 to 535 mg/L. (Based on the Tables 3 to 15 given in Examples 3 to 15).

EXAMPLE 64:

Dissolved 0.15 g glucose in 25 mL distilled water. Aliquots from the stock solution were taken and diluted with water to achieve different COD concentrations in the

range of 320 to 535 mg/L. A 30 mL sample of Damaged wheat grain (ssp) was taken. Weighed 0.02 g HgSO_4 in separate 15 mL test tubes. 1.0 mL of glucose solution and sample were added to each test tube separately and mixed well. 0.5 mL of $\text{K}_2\text{Cr}_2\text{O}_7$ was added to each sample. 1.5 mL H_2SO_4 silver sulphate reagent was added slowly with continuous shaking. The color of the reaction mixture was noted down visually within 15 sec of incubation and also recorded spectrophotometrically at 585 and 635 nm against air. An approximate COD value of the sample was estimated from the color of standard reaction mixture prepared with glucose as reference material in the COD range of 200 to 100000 mg/L. OD 585 and OD 635 of the reaction mixture was read spectrophotometrically. OD 585 and 635 nm in the range of 0.130 to 0.175 and 0.120 to 0.160 respectively were considered and calculated with reference to glucose standard. (Glucose COD in the range of 320 to 535 mg/L corresponds to OD 585 and 635 nm in the range of 0.157 to 0.207 and 0.142 to 0.189, respectively).

COD of the sample was estimated to be 6649.3 mg/L. (By Merck method and read on photometer).

Based on the color of the reaction mixture, COD of the sample was estimated to be around 500 mg/L. On the basis of the dilution method and the spectrophotometric reading at OD 585 and OD 635 of the sample, COD was estimated to be 5260 and 5190 mg/l respectively. A comparison with the COD value calculated through Merck method reveals the following relation i.e. Our COD value (at OD 585) $\times 1.26 =$ COD value (Merck Method) or our COD value (at OD 635) $\times 1.30 =$ COD value (Merck Method).

Hence for a preliminary estimate of COD value color of the standard reaction mixture based on a range of glucose concentrations can be used as a reference. For samples, with COD concentration of 300 to 500 mg/L, OD 585 or OD 635 can be recorded for obtaining more precise values.

DWS (S3) (ssp)

Table 64: Correlation between COD concentration, color of the reaction mixture and their OD at different wave lengths.

COD (mg/L)	Sample dilution (Times)	Color of the reaction mixture	OD	Estimated COD (mg/L)	Estimated final COD (mg/L) (Col. 2x5)
1	2	3	4	5	6
<u>Values at OD 585 nm</u>					
Glucose:					
320	NA	C300	0.159	NA	NA
426	NA	C400	0.187	NA	NA
535	NA	C550	0.210	NA	NA
Sample: DWS (S3) (ssp)					
	Nil	C550	-	550	550#
	10	C400	0.172	526	5260
<u>Values at OD 635 nm</u>					
Glucose:					
320	NA	C300	0.142	NA	NA
426	NA	C400	0.167	NA	NA
535	NA	C550	0.189	NA	NA
Sample: Digested POME (ssp)					
	Nil	C550	-	550	550#
	10	C400	0.151	531	5310

#: COD as deducible from Table 2 given in Example 2. D: Dark. T: Turbid.

*: OD out of range. NA: Not applicable. Samples are diluted to get an OD 585 in the range of 0.130 to 0.175 and OD 635 in the range of 0.120 to 0.160, to obtain a COD value in the range of 320 to 535 mg/L. (Based on the Tables 3 to 15 given in Examples 3 to 15).

EXAMPLE 65:

Dissolved 0.15 g glucose in 25 mL distilled water. Aliquots from the stock solution were taken and diluted with water to achieve different COD concentrations in the range of 320 to 535 mg/L. A 30 mL sample of Tea, Apple Pomace, Damaged wheat grain and Glucose (ssp) was taken. Weighed 0.02 g HgSO_4 in separate 15 mL test tubes. 1.0 mL of glucose solution and sample were added to each test tube separately and mixed well. 0.5 mL of $\text{K}_2\text{Cr}_2\text{O}_7$ was added to each sample. 1.5 mL

H₂SO₄ silver sulphate reagent was added slowly with continuous shaking. The color of the reaction mixture was noted down visually within 15 sec of incubation and also recorded spectrophotometrically at 585 and 635 nm against air. An approximate COD value of the sample was estimated from the color of standard reaction mixture prepared with glucose as reference material in the COD range of 200 to 100000 mg/L. OD 585 and OD 635 of the reaction mixture was read spectrophotometrically. OD 585 and 635 nm in the range of 0.130 to 0.175 and 0.120 to 0.160 respectively were considered and calculated with reference to glucose standard. (Glucose COD in the range of 320 to 535 mg/L corresponds to OD 585 and 635 nm in the range of 0.157 to 0.207 and 0.142 to 0.189, respectively).

COD of the sample was estimated to be 30,866 mg/L. (By Merck method and read on photometer).

Based on the color of the reaction mixture, COD of the sample was estimated to be around 20000 mg/L. On the basis of the dilution method and the spectrophotometric reading at OD 585 and OD 635 of the sample, COD was estimated to be 20,850 and 20,450 mg/l respectively. A comparison with the COD value calculated through Merck method reveals the following relation i.e. Our COD value (at OD 585) x 1.48 = COD value (Merck Method) or our COD value (at OD 635) x 1.51 = COD value (Merck Method).

Hence for a preliminary estimate of COD value color of the standard reaction mixture based on a range of glucose concentrations can be used as a reference. For samples, with COD concentration of 300 to 500 mg/L, OD 585 or OD 635 can be recorded for obtaining more precise values.



Tea + Apple Pomace + DWS (S3) + Glucose (ssp)

Table 65: Correlation between COD concentration, color of the reaction mixture and their OD at different wave lengths.


COD (mg/L)	Sample dilution (Times)	Color of the reaction mixture	OD	Estimated COD (mg/L)	Estimated final COD (mg/L) (Col. 2x5)
1	2	3	4	5	6
<u>Values at OD 585 nm</u>					
Glucose:					
320	NA	C300	0.159	NA	NA
426	NA	C400	0.187	NA	NA
535	NA	C550	0.210	NA	NA
Sample: Tea + Apple Pomace + DWS (S3) + Glucose (ssp)					
	Nil	C10000	-	-	>10,000
	10	C2000	0.359	2000	20,000#
	50	C400	0.152	417	20,850
<u>Values at OD 635 nm</u>					
Glucose:					
320	NA	C300	0.142	NA	NA
426	NA	C400	0.167	NA	NA
535	NA	C550	0.189	NA	NA
Sample: Tea + Apple Pomace + DWS (S3) + Glucose (ssp)					
	Nil	C10000	-	-	>10,000
	10	C2000	0.359	2000	20,000#
	50	C400	0.152	409	20,450

#: COD as deducible from Table 2 given in Example 2. D: Dark. T: Turbid.

*: OD out of range. NA: Not applicable. Samples are diluted to get an OD 585 in the range of 0.130 to 0.175 and OD 635 in the range of 0.120 to 0.160, to obtain a COD value in the range of 320 to 535 mg/L. (Based on the Tables 3 to 15 given in Examples 3 to 15).

The main advantages of the present invention are

1. The test can be carried out rapidly.
2. The test does not require any major instrument.
3. The test is very cheap.

- 
4. The test requires very small amount of reagents.
 5. The test is easy to perform.
 6. The test can be done without the need for specially trained technician.
 7. The test can be performed with the help of a kit also.

106260" 8E302360